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Material description

Cenit L is a light-curing PMMA-resin in dosage syringes for using in the laboratory in the laboratory.

Field of application:

Cenit L is excellently suitable as space maintainer and for the application of reservoirs during the production of medication supports resp. bleaching splints due to the especially adjusted thixotropy. An exact determination of the thickness of the layer is easily possible thanks to the increased colour intensity. Further fields of application are the filling of air bubbles in the working model as well as the blocking-out of undercut areas and the building-up of broken teeth.

Processing:

- Apply the required quantity of Cenit L onto the model area that is to be repaired or covered and mould it with a brush if necessary.
- Place the model into a light-curing unit (e. g. EyeVolution®) for polymerisation which disposes of a spectral area of 300–500 nm. According to the light intensity the polymerisation time should vary between 3 and 5 minutes.
- After the polymerisation of Cenit L the inhibition layer may be removed by means of a soft cloth and an alcohol containing cleaning agent.
- Empty the dosage tips after having used the syringe by means of pulling the piston back carefully in order to avoid a polymerisation of the material in the tip.

Contains

tetrahydrofurfuryl-2-methacrylate, Diphenyl(2,4,6-trimethylbenzoyl)phosphinoxid

Hazard statements

H317 May cause an allergic skin reaction. H360 May damage fertility or the unborn child.

Keep out of the reach of children!

Only for dental use!

Storage conditions:

Store at room temperature, protected from light.

Shelf life: see date of expiry

Lot number / Expiry date

The lot number and the expiry date are indicated on each Cenit L packaging. In case of claims, please always indicate the lot number of the product. Do not use the product.

Delivery form:

REF **D3815** 4 syringes of 2.5 g each
+ 8 dosage tips

The product has been developed for the dental field and has to be processed according to the patient information. For damages which can result from ulterior use or faulty processing the manufacturer does not assume any liability. Furthermore, the user is bound to check the material to aptitude and purpose before use, especially when the purpose is not mentioned in the patient information.

Precautionary statements

P261 Avoid breathing dust / fume / gas / mist / vapours / spray. P280 Wear protective gloves / protective clothing / eye protection / face protection. P308+P313 If exposed or concerned: Get medical advice / attention. P333+P313 If skin irritation or rash occurs: Get medical advice / attention. P405 Store locked up. P501 Dispose of contents / container in accordance with local / regional / national / international regulations.

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