**Bedienungsanleitung · Working instructions**  
**Mode d'emploi · Instrucciones de trabajo**  
**Istruzioni d'uso · Gebruiksaanwijzing**

# Drufosmart

**Basis Druckformgerät 2,5 bar**

**Basic pressure thermoforming unit 36 psi**

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## Contents (EN)

1.	Description of unit	22
2.	Technical data	22
3.	Declaration of conformity	23
4.	Safety advice	24
5.	Initiation / scope of delivery	25
6.	Functional components	26
7.	Operation	27
7.1	Adjustment heating time	27
7.2	Start of programme	28
7.3	Thermoforming	29
7.4	Opening of the pressure chamber	30
8.	Deviation to the rule	31
8.1	Modification of given times	31
8.2	Plasticizing of blanks with unknown heating times	31
8.3	Abort / termination of programme	31
9.	Additional optional functions	32
9.1	Fabrication of cast form parts	32
10.	Exchange of fuses	33
11.	General advice	34
12.	Care and maintenance	34
13.	Spare part list	35
14.	Trouble shooting	36
15.	Explosion sketch	102
16.	Operating diagram	103
17.	Pneumatic plan	104

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# Drufosmart

## 1. Description of unit

The Drufosmart is a basic pressure thermoforming unit as an alternative to common vacuum thermoforming units for daily use in the dental laboratory and dental surgery.

Simplest handling as well as continuous digital and acoustic control of all

working steps ensure problem-free working, even with poor precognition.

A forming pressure of 2.5 bars in connection with a vertically forming blank guarantee for very precise results.

## 2. Technical data

Dimensions (H x W x D):	420 x 320 x 290 mm
Weight:	15 kg
Supply voltage:	230 V / 50–60 Hz, 115 V / 50–60 Hz
Max. power consumption:	315 W
Fuse of unit:	T 3.15 A
Operating pressure:	2.5 bars
Min. operating pressure:	2 bars
Max. pressure in pipe:	10 bars

### 3. Declaration of Conformity

#### According to Machinery Directive (2006 / 42 / EC) attachment II 1.A

Herewith we declare that the machine described below is in its conception and design and in the shape delivered by us in accordance with the fundamental requirement for safety and health as prescribed in the applicable EC-guidelines. In the case of any change or modification of the machine not authorised by us this declaration becomes invalid.

**Designation of machine:** Drufosmart  
D3200 / D3200A / D32001 / D32001A

**Type of machine:** Pressure thermoforming unit

**Applicable EC directives:**  
2006 / 42 / EC Machinery Directive  
2014 / 35 / EU Low Voltage Directive  
2014 / 30 / EU EMV Directive

#### Applied harmonised standards:

##### 2006 / 42 / EC

EN ISO 12100:2010  
+ Corrigendum to  
DIN EN ISO 12100:  
2011-03

##### According to Machinery Directive

Safety of machinery – General principles for design – Risk assessment and risk reduction (ISO 12100:2010); German version EN ISO 12100:2010 + Corrigendum to DIN EN ISO 12100:2011-03.

EN ISO 13732-1:2008

Ergonomics of the thermal environment – Methods for the assessment of human responses to contact with surfaces – Part 1: Hot surfaces (ISO 13732-1:2006); German version EN ISO 13732-1:2008-12.

DIN EN 60204-1:  
2006 / A1:2009  
+ Corrigendum to  
DIN EN 60204-1:2006

Safety of machinery – Electrical equipment of machines – Part 1: General requirements (IEC 60204-1:2005/A1:2008); German version EN 60204-1:2006/A1:2009 + Corrigendum to EN 60204-1:2006; German version CENELEC-Cor.: 2010 to EN 60204-1:2006.

##### 2014 / 35 / EU

##### Low Voltage Directive

EN 61010-1:2010

Safety requirements for electrical for measurement, control and laboratory use – Part 1: General requirements (IEC 61010-1:2010 + Cor.: 2011); German version EN 61010-1:2010.

##### 2014 / 30 / EU

##### EMV Directive

EN 61000-3-2:2014

Electromagnetic compatibility (EMC) – Part 3-2: Limits – Limits for harmonic current emissions (equipment input current  $\leq$  16 A per phase) (IEC 61000-3-2 :2014); German version EN 61000-3-2:2014.

EN 61000-3-3:2013

Electromagnetic compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current  $\leq$  16 A per phase and not subject to conditional connection (IEC 61000-3-3:2013; German version EN 61000-3-3:2013.

EN 61326-1:2013

Electrical equipment for measurement, control and laboratory use – EMC requirements; Part 1: General requirements (IEC 61326-1: 2012); German version EN 61326-1:2013. Emission according to living area, business and industrial undertakings as well as small enterprises; interference stability according to industrial area.

## 4. Safety advice

**Attention!** Read the following notes before connecting and putting the unit into service. The operating service and the function of the unit can only be guaranteed for, if the general safety and accident prevention instructions of the legislation as well as the safety instructions of the working instructions are observed.

1. This unit must only be used according to the present working instructions. We cannot be held liable for damages which have been caused by improper treatment resp. faulty manners of operation.
2. Put the unit on a stable and even basis which is suitable for the operation weight of approx. 16 kg of the unit.
3. In order to avoid the intrusion of water into the unit (e. g. splash water) the unit is to be placed into a dry surrounding.
4. Do not store any flammable substances in the close surrounding of the unit.
5. The indicated voltage on the type plate has to meet the voltage of the power source.
6. The unit must only be connected with a socket with protective conductor. Never touch the plug with wet hands.
7. Authorised operators: The operating company of the machine has to hold the working instructions at the operators' disposal and has to make sure that authorised operators have read and understood the manual. Only in that case the operators may operate the unit.
8. The unit has to be examined before operation in view of proper condition and operating safety. If the unit is not in proper condition, it must not be used and has to be marked accordingly.
9. The pressure of the compressed air line must not exceed 10 bars.
10. Do not intrude any objects into the unit.
11. Do not lock-up the operation elements of the unit.
12. Do not block the downwards-moving piston.
13. Do not touch the downwards-moving piston.
14. The swivelling housing of the radiator reaches high temperatures and must therefore not be touched – not even for short moments.
15. If the unit is not used, please disconnect the mains.
16. Plates and labels have to be maintained in a well legible condition and must not be removed.
17. Before cleaning and maintenance of the unit or the exchange of parts the plug has to be disconnected.
18. The opening of the unit and repairs are only to be carried out by accordingly trained experts.

19. Only those appliances and spare parts may be used which are released by the manufacturer. We are not responsible for damages caused by the use of foreign parts.

20. Arbitrary rebuilding or modifications are not permitted due to safety reasons.

## 5. Initiation

Before initiation make sure that the voltage indicated on the type plate meets the voltage of the power source.

Connect the unit to an air compressor. Push the pressure hose (included in the delivery) into the compressed air connection (9) until it is locked in place. The hose does not need any further fastening. To release the pressure hose push the blue ring upwards and draw the hose off by pulling downwards. Assemble the fast coupling at the other end of the tube and connect it with the pressure line.

**The pressure of the compressed air line must not exceed 10 bars.**

Put the power cord into the female connector (33) and create a connection with the Euro/USA plug.

21. The required operation and maintenance conditions mentioned in these working instructions are to be complied with imperatively. When working with the Drufosmart the general accident prevention regulations have to be observed.

**Important:** Before the first initiation please read the working instructions thoroughly. In case of comprehension problems contact your responsible dealer or Dreve Dentamid GmbH, if necessary.


### Scope of delivery

- 1 Drufosmart
- 1 drawer for granules
- 1 insert bowl and nub plate
- 1 tension ring with exhaust valve
- 1 base range of blanks
- 1 power cord
- 1 air connection hose
  - incl. hose clamps and fast coupling
- 1 working instructions in 6 languages
- 1 information folder thermoforming foils

## 6. Functional components

No.	Description
5	Heating handle
6	Nub plate and insert bowl
7	Drawer for granules
9	Compressed air supply
10	Handle plate reception
11	Plate reception
12	Ventilation valve
13	Tension ring
20	Heating
27	Display
33	Mains supply / fuse carrier
34	Activation lever
43	Main switch
44	Turnable pressure controller

## 7. Operation

Turn the unit on by means of the main switch (43) at the back of the unit. By this the display (27) at the front of the

unit will be activated and shows 0:00. An acoustic signal shows that the unit is ready for operation.

### 7.1 Adjustment heating time


By means of the turnable pressure controller (44) every time between 0:05 min. and 9:55 min. can be adjusted by turn-ing the controller clockwise. Adjustment is effected in 5 sec. steps.

Please refer to the enclosed information folder in order to find out the optimum times for Dreve thermoforming foils.

## 7.2 Start of programme

Before starting the adjusted heating time put a corresponding blank into the plate reception (11) and fix it by means of the tension ring (13).

The model can either be put onto the nub plate (6) or fixed by taking out the nub plate and embedding it into granules directly in the insert bowl.

### Info

**For an optimum adaptation of the blank towards the model always place the models in a way that the anterior teeth are directed towards the centre of the blank. When using hard blanks it is recommended to embed the models into granules. Models that are placed onto the nub plate should not exceed a max. height of 23 mm.**

Swivel the plate reception inwards at the handle (10) just up to the mechanical interlock in the centre of the unit (magnetic interlocking). Also swivel in the heating element (20) at the corresponding handle (5) just up to the

mechanical interlock in the centre of the unit. Close the ventilation valve (12) of the tension ring (13). In order to start the heating process press the turnable pressure controller (44). The display shows a countdown towards 0:00 of the adjusted time.

### Info

**For safety reasons the heating is equipped with a final switch and only works in swivelled in position after having started a programme. If the heating carrier (20) is swivelled outwards at the handle (5) after having started a programme, the heating process will be interrupted automatically. The indication in the display will stop and show the remaining time. By swivelling the heating carrier (20) back in the programme will be re-activated.**

**Attention! The swivelling housing of the heating element reaches high temperatures. Therefore it must not be touched – even not for short moments.**

**After the adjusted time has elapsed, there will be an acoustic signal, the indicator in the display will blink for 4 seconds and the cooling time will be counted up again automatically.**

## 7.3 Thermoforming

After the heating time has elapsed there will be an acoustic signal; the heating time in the display will be 0:00. The heating will be switched off automatically. Pull the heating carrier (20) with the right hand at the handle (5) completely out of the centre of the unit. Then push the activation lever (34) with the left hand downwards as far as it will go.

By pulling the heating handle (5) again for approx. 3 seconds as far as it will go pressurising of the piston will be carried out automatically. This is confirmed by an acoustic signal.

### Info


**For safety reasons the pressure chamber can only be pressurised, if it has been moved downwards completely manually. Always make sure that no objects are in the contact area between piston and insert bowl while the piston is closing! In order to protect the operating personnel securely against bruises of the fingers, pressurising can only be effected, if the unit is operated with two hands at the same time (two-hand control).**

The forming pressure will be verified continuously and regulated, if necessary. The display shows the duration since turning off the heating (cooling time). The maximum time that can be shown is 9:55 min. Please refer to the enclosed information folder in order to find out the optimum cooling times for the used Dreve thermoforming foil.

## 7.4 Opening of the pressure chamber

In order to attain a complete cool down of the thermoforming blank the valve at the tension ring (12) is to be opened slightly after 2 / 3 of the adjusted cooling time. Through this the warm air inside the piston will be released and replaced by cool air due to the automatic pressure regulation. After the cooling time has elapsed pull the heating handle (5) as far as it will go and maintain it in this position (approx. 3 seconds) until an acoustic signal can be heard. Release the heating handle and the pressure will be relieved automatically.

Then swivel the pressure chamber into its original position by leading back the activation lever (34).

Swivel out the plate reception (11) at the corresponding handle (10) and remove tension ring and thermoformed result.

## 8. Deviation to the rule

### 8.1 Modification of given times

The specified times may be prolonged or shortened individually, if required.

#### Info

**The performance of the heating might decrease slightly after 1000 operating hours. By prolonging the heating times by 5–10 seconds**

**optimum thermoforming results can be attained again. The given times are average times and designed for general use. In case of extreme undercut areas or very high models a prolongation of the heating time can optimise the results.**

### 8.2 Plasticizing of blanks with unknown heating times

Create a heating programme corresponding to the foil thickness in question with a sufficient prolongation of the heating time (see point 7.1). Start the programme (see point 7.2). You may test the plasticization degree of the blank by means of an edgeless metal instrument (e. g. wax knife).

#### Info

**The complete deformability of the blank can only be tested at the edges, as that is where warmth will leak through the metal con-**

**tacts and the necessary flexibility will be attained at last!**

After the blank has attained the necessary flexibility initiate the thermoforming process (see point 7.3). The cooling time can be terminated at any time by pushing the heating handle (5). However, a sufficient cooling down time is to be observed, as warm foils may deform during removing them from the model and therefore lead to bad resp. useless thermoforming results.

### 8.3 Abort / termination of programme


It is possible to abort / terminate all actions of the Drufosmart by pushing

the turnable pressure controller (44) for at least 3 seconds.

## 9. Optional additional functions

### 9.1 Fabrication of cast form parts

With the K7 insert it is possible to thermoform up to 7 resin bases for the crown and bridge technique at the same time. Fix the single dies by means of a flexible blocking out material (e. g. Fillin) in the holes of the K7 insert.


All further working steps as described from point 7.1 on.

## 10. Exchange of fuses

Disconnect the unit from mains supply. The fuse carrier is situated at the rear of the unit, integrated in the mains supply.



In order to control or exchange the fuses push the push catch of the carrier, keep it pressed and remove the carrier from the unit.

Pull the defective fuse out of the clamp carrier applying slight traction and replace it by a new, equivalent one.

**Important: Only use the fuse strength indicated on the rear.**

After exchanging the fuse insert the fuse carrier into the housing applying slight pressure until the fixation latches noticeably.

## 11. General advice

Repairs on the Drufosmart units may only be executed by experts. Disconnect the unit from mains supply at any rate. Dreve company is only responsible for the safety, reliability and performance of the unit, if:

- the unit is exclusively used for the corresponding field of application.
- operation is only carried out by accordingly trained personnel.

- the unit is operated in accordance with these working instructions.
- extensions, new adjustments, changes or repairs are executed by persons authorised by the company.

As we are currently developing our products, we reserve the right for technical changes.

## 12. Care and maintenance

Generally disconnect the unit from mains supply before executing care and maintenance.

Ideally cleaning is to be effected dryly with a soft piece of cloth, if necessary with a slightly humid sponge and a mild cleaning liquid. Water and cleaning liquids must not get into the unit.

The O-rings should be greased regularly with silicone grease.

### **Removal and cleaning of the filter sieve**

Disconnect the power cable from the socket and shut down the compressed air supply. Push the blue ring of the compressed air connection upwards and draw the hose off the unit by pulling it downwards. Now unscrew the connection using a suitable tool. Then carefully remove the filter sieve with a slotted screwdriver by turning anti-clockwise. Clean with stream of compressed air and reassemble in reverse order.

### 13. Spare part list

No.	Description	REF
5	Heating handle	D6505519
6	Insert bowl with nub plate	50591/50592
7	Drawer for granules	5565
9	Compressed air connection, complete	55637
10	Handle plate reception	D6505519
11	Plate reception	D6505512
12	Ventilation valve	50761
13	Tension ring	D50038
17	Housing red	D65555R
19	Micro-switch	55790
20	Heating housing	D6505513
21	Quartz-infrared-radiator	51515
22	Temperature monitor	5542
23	Pillar heating element	55773
25	Magnetic valve	52021230
26	Cylinder	52019
27	Display foil, adhesive	55820
28	Limit plate reception	D6555344
29	Support plate reception	D6505514
30	Pillar plate reception	55772
31	Bottom plate	D6505511
33	Mains supply / fuse carrier	51208
34	Activation lever	D65885
35	Connector lever / axis	D65017
36	Bushings	55708
37	Final switch activation lever	E5521
38	Final switch heating	55790
39	Recuperating spring	55861/55862
40	Lever axis	D6555847
43	Main switch	50841
<b>without drawing no.</b>		
	Power cord Europe	51021
	Power cord England	51284
	Power cord America	51283
	Electronics incl. display, key pad, connection cable	55201
	Specify supply voltage or serial number.	
	Fuse F3.15 AT	51321
	O-rings 112 x 3 mm	51077
	Cable harness complete	55521
	Limit heating element	D65016

## 14. Trouble shooting

<b>Problem</b>	<b>Cause</b>	<b>Aid</b>
Unit does not work	Power line is not correct	Control power supply
	Main breaker is off	Switch on main breaker
	Miniature fuse is defect	Check miniature fuse and exchange it, if necessary
Main breaker on, electronics show no function	Electronics are out of order	Electronics have to be exchanged by an expert
Piston cannot move downwards	Heating is not swivelled outwards	Swivel heating outwards
Piston has moved downwards, unit does not take over self-maintenance	Minimum pressure of 2 bars has not been reached	Hold escapement for longer, see working instructions
	Electronics are out of order	Electronics have to be checked by an expert
	Unit is not connected to compressed air	Connect unit to compressed air
Piston does not move upwards	Wrong operation of push button	Pull push button for 3 sec.
	Contact pressure has not been relieved	Relieve contact pressure, see working instructions
	Final switch is defect	Final switch has to be exchanged by an expert
	Magnet valve or electronics are defect	Unit has to be controlled by an expert

<b>Problem</b>	<b>Cause</b>	<b>Aid</b>
Electronics can be set, but radiator does not heat	Heating element is not swivelled-in completely	Swivel the heating element in correctly
	Micro-switch defect	Micro-switch has to be exchanged by an expert
Heating time elapses, but radiator does not heat	Quartz radiator defect	Exchange quartz radiator
	Temperature monitor was activated	Exchange temperature monitor
Heating time elapsed, blank is not plasticized sufficiently	Lifetime of quartz radiator exceeded	Exchange quartz radiator
	Heating time was too short	Prolong heating time by 5–10 sec.
Unit does not reach the pre-adjusted pressure	Line pressure too low	Increase line pressure
Unit makes fizzy sounds during thermoforming	Ventilation valve not closed	Close ventilation valve
	Tension ring has been placed wrongly into plate reception	Put tension ring into correct position
	Connection hose leaky between valve and cylinder	Unit has to be controlled by an expert, connection hose has to be exchanged, if necessary

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