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Stand der Information: 10.2021

**Intended use**

Flexible gingiva mask material for the C&B technique

**Material description**

**GumQuick scan** is an addition-vulcanising vinylploysi-loxane for the production of flexible gingiva masks on the master model.
**GumQuick scan** reaches a final hardness of approx. 63 Shore A and adapts the colour naturally. The material is characterised by its high tear-resistance and is easy to treat with rotating instruments.

**Field of application**

Gingiva masks produced in the laboratory for the c+b technique and for high-precision implant works.

**Technical data according to ISO 4823**

<b>GumQuick scan</b>	
Consistency:	Type 3
Processing time:	1: 30 min.
Curing time:	≥ 5: 00 min.
Final hardness:	63 ± 3 Shore A
Recovery from deformation:	99.6 <span> </span> %
Linear dimensional change:	< -0,1 <span> </span> %

The mixing and processing times refer to a room temperature of 23°C ± 1°C and a relative air humidity of 50 %. Whereas lower temperatures prolong these times, higher temperatures shorten them.

**Processing**

**Indirect method**

A matrix of hard kneadable silicone material (e. g. **Laborsil**) is to be applied onto the dies of the unsawn die model which should be provided with the gingiva mask. After curing the kneadable silicone should be cut so that the top of the dies will appear (approx. 1 mm).

The model can be sawn and gingiva parts ground as

usual. Afterwards the correct position of the matrix has to be checked. In order to guarantee a faultless surface of the gingiva mask **DuoSep** is to be applied onto the contact surface of the kneadable silicone / **GumQuick scan**.

**GumQuick scan** is to be applied directly from the cartridge onto the model and into the matrix. Now the matrix has to be positioned at the model so that the surplus of material can flow out through the holes at the top of the die. Ideally, the gingiva mask should be cured under pressure.

After a period of 5–6 min. the matrix can be removed and the mask can be treated with scissors, scalpel or with rotating grinding instruments.

**Direct method**

At first **DuoSep** has to be applied onto the contact surface impression material / **GumQuick scan** in order to guarantee a faultless surface of the gingiva mask.

**Attention**

In order to avoid a damage of the gingiva mask surface when applying **GumQuick scan** in conjunction with polyether impressions **DuoSep** has to be used! Before the production of the gingiva mask the polyether impression should ventilate for at least 3 hours.

Afterwards **GumQuick scan** is to be applied directly from the cartridge into the impression. After a period of 5–6 min. **GumQuick scan** will be cured and the impression can be filled with plaster. The gingiva mask can easily be treated with scissors, scalpels or with rotating grinding instruments.

**Important advices**

- Due to technical filling reasons an unequal transportation of the components can occur when using a new cartridge. For this reason the first transported quantity – approx. size of a pea – should not be used.
- Never try to squeeze the material through a cured mixing canula. This may lead to leakage in the cartridge system and therefore to transportation of un-mixed and uncured single components.
- Before each application the cartridge has to be checked for consistent transportation of material.
- The vulcanisation can be disturbed by the use of latex gloves, hand creams, cleaning agents etc.
- The product is exclusively designed for the above mentioned field of application and should only be used by qualified and instructed persons.

**Storage conditions**

At room temperature, i. e. 64.4 °F–82.4 °F (18 °C–28 °C)

Shelf life: 2 years in the closed original container

<b>Delivery form</b>	
<span><span><span></span></span></span> <b>[REF]</b> <b>D504890</b>	2 x 50 ml double cartridges, 12 mixing canulas, 10 ml <b>DuoSep</b>
<span><span><span></span></span></span> <b>[REF]</b> <b>D4970X3</b>	40 mixing canulas, blue
<span><span><span></span></span></span> <b>[REF]</b> <b>D49701X3</b>	100 mixing canulas, blue
<span><span><span></span></span></span> <b>[REF]</b> <b>D4306</b>	10 ml <b>Duo Sep</b>

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