Uso previsto

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Intended use

Dental splinting material long-term

Extended purpose

Soft/hard thermoforming foil made of polyethylene terephthalate and ethyl vinyl acetate for the production of functional splints such as occlusion or bruxism splints.

Indications for use

Suited for interim splints in case of myoarthropathy (muscular and mandibular discomfort) and accompanying dysfunctions resulting from it, like pressing, bruxism, etc. The splint provides for a very good bite lock and does not cause any tension. Furthermore the production is very easy and there is no need for anchorages, the splint fits very well and secure due to the snap-in effect.

Contraindication

In case of allergies to any of the ingredients do not use Kombiplast. Not desired side effects of this medical device are not to be expected in case of proper processing and application. However, immune reactions (e. g. allergies) or irritations principally cannot be excluded. Should you get to know non-desired side effects – also in doubtful cases – we appreciate in any case to be informed with a detailed description of the attendant circumstances and symptoms. We will follow each hint.

Target patient group

Patients who need dental care.

Intended users

Trained staff in dental field

Processing

All Kombiplast materials are delivered with an insulating foil. This extremely thin foil (0,07 mm) is thermoformed together with the Kombiplast on the model side and guarantees a decisively improved transparency of the splint. Thus an additional insulation of the model surface is no longer necessary. The heating-up of the Kombiplast foils is done from the hard side (the foil is only suitable for pressure forming units). The adjustment of the heating and cooling times can either be taken from the survey of thermoforming foils, or can be

programmed directly into the unit with the help of the bar-code on the Kombiplast packaging. If a Drufomat scan or a Drufosmart scan is used, please take note of the special information included in this package.

Cleaning and desinfection

The adhesive surfaces must not come into contact with PMMA or other alcoholic cleaning agents. Alcohol-based disinfection may also lead to a delamination of the two foils. For cleaning, the use of MaxiFresh™ cleaning tablets is recommended. Please be aware of the MaxiFresh™ instruction for use.

Note

While working out on the material, it is recommended to use an exhaust system and protective glass or goggles as well as a face mask.

Serious incidents

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/ or patient is settled.

Storage conditions

At room temperature, i. e. 18–28 °C. Protect against humidity, always close foil bag when packs are in use.

Lot number / Expiry date

The lot number and the expiry date are indicated on each packaging. In case of claims, please always indicate the lot number of the product. Do not use the product after the expiry date.

Disposal

Disposal of contents / container in accordance with local / regional / national/international regulations

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This product was developed for dental use and has to be used according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence

before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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