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

The mixing and processing times refer to a room temperature of 23 °C and a relative air humidity of 50%. Lower temperatures prolong, higher temperatures shorten these times.

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Description of material

Laborsil pro 70/90 are 1:1 addition-vulcanising kneadable silicones for the universal use in the dental laboratory with a hardness of 70 and 90 Shore A.

Indication

For the production of copings, of repair models with under-growing areas, to block out undercuts and embedding silicone for the denture pressure procedure.

Processing mixing

Take two equal-sized amounts of components A and B and knead them manually until a homogenous colour develops.

Important advice

Polymerisation may be disturbed by use of e. g. latex gloves, creams, cleaning agents, resins etc.

Keep away from children!**Only for use in the dental laboratory!****Storage condition**

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This product was developed for dental use and has to be used according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.