**Bedienungsanleitung
Instructions for use
Mode d'emploi**

deutsch

english

français

PCU LED N₂ Pioneering Technology

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**Light polymerization device for the
CAD / CAM 3D printing process**

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Contents (EN)

1. Device description	24
2. Technical data	26
3. Declaration of conformity	26
4. Safety instructions	27
5. Commissioning	29
6. Information on the use of the protective gas	30
7. Functional elements	30
8. Operation	31
9. Maintenance	37
10. Monitoring	38
11. Troubleshooting plan	38
12. Exploded view	40
13. Spare parts list	41
14. General information	42
15. Decommissioning of the product	42

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PCU LED N₂

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1. Device description

The PCU LED N₂ (Post Curing Unit LED) is a light polymerization device for the curing of components which have been created by means of 3D printing. The following device properties are used in the curing process:

Intended use:

The PCU LED N₂ may only be used for curing of resins and lacquers which are crosslinking in a wavelength range between 365 and 425 nm.

The following device properties are used in the curing process:

- LEDs must not be covered.
- Use only in dry environments.
- Operate only with closed lid.
- The needle valve may only be adjusted by trained maintenance personnel (manufacturer).
- The unit may only be used with nitrogen (N₂) as the inert gas.
- A suitable coupling socket must be used on the N₂ connection.
- Only UBS sticks with a maximum memory size of 8 GB may be used.
- The maximum permissible length of the connection cable is 3 m.
- Position the post-curing insert only in the intended orientation in the exposure chamber.
- Do not lift the unit by the handle of the lid.

Device properties:

In the curing process, the following device properties are applied:

a) LED technology

The advantages of the LED technology are that a permanent light of one wavelength is emitted. The permanently emitted light radiation can be captured by a standard UV sensor and thus allows for the measurement and the observation of the light output in the exposure chamber. Furthermore, LED lamps do not produce any heat which could cause the material properties of the plastic component to be cured to be negatively influenced. Another advantage of the LED technology is that the stated operating life of the LED lamps is approx. 20,000 hours and is accompanied by low operating and maintenance costs.

b) Electronic control (micro-controller)

The use of a micro-controller allows for the control and reading of the current, pressure and temperature sensors.

c) Logging and Monitoring

The information which is recorded by the current, pressure and temperature sensors represents the used process parameters (e.g. light on/off; vacuum on/off; nitrogen on/off). These process parameters are recorded, saved and transferred to a USB storage medium. This allows for the transparent traceability of the production chain of the medical device.

d) Simple control concept

The device is operated by the removable rotary button. The corresponding option is selected by turning the button. The selected option is operated by pressing. Longer pressing opens the parent and sub-folder menu. Both elements can be cleaned with the appropriate cleaning agents after removal.

e) Process safety

To be able to guarantee a high process liability, the control software is subdivided into two levels. Level 1 is accessible to the user and process owner. Level 2 is only accessible to the process owner. In level 1 only a prede-

defined process can be selected. In level 2, for example, a process can be defined and an option can be switched on and/or off.

f) Curing environment

The components to be cured can be cured under different environments. The environment to be set depends on the application and can be found in table 1.

g) Open system

In the PCU LED N₂, 10 freely programmable storage spaces can be used.

Table 1: Curing environment for different uses

Environment	Use	PCU LED N ₂ nitrogen
Vacuum	Curing of laboratory products	✓
Nitrogen	Curing of medical devices	✓

PCU LED N₂

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2. Technical data

Weight of the overall package	11.1 kg
Weight of the device	9.3 kg
Dimensions of the overall package	H 220 x W 455 x D 370 mm
Dimensions of the device	H 110 x W 389 x D 276 mm
Polymerization chamber	H 65 x W 150 x D 150 mm
Voltage	100–240 V
Nominal frequency	50–60 Hz
Max. Power consumption	110 W
Max. altitude	2000 m
Ambient temperature	5–40 °C max. relative humidity 80 % at temperatures up to 3126 °C, decreasing linearly up to 50 % relative humidity at 40 °C

Fluctuations in the mains supply voltage up to ± 10 % of the nominal voltage transient overvoltages up to the values of overvoltage category II temporary overvoltages.

Permissible degree of soiling 2

3. Declaration of conformity

According to Machinery Directive (2006 / 42 / EG), Low Voltage Directive (2014 / 35 / EU) and EMV Directive (2014 / 30 / EU)

Herewith we declare that the machine described below in its conception and design and in the shape delivered by us is in accordance with the fundamental requirements for safety and health as prescribed in the applicable EC directives. In the case of any change or modification of the machine, not authorised by us, this declaration becomes invalid.

Denomination of the machine: PCU LED N₂
REF 4317

Type of machine: Light polymerization unit

Applicable EC directives: 2006 / 42 / EC According to Machinery Directive
2014 / 35 / EC Low Voltage Directive
2014 / 30 / EC EMV Directive

Applied harmonized standards:

2006 / 42 / EC According to Machinery Directive

EN ISO 12100:2010 Safety of machinery – General principles for design – Risk assessment and risk reduction
+ Corrigendum to (ISO 12100:2010); German Version EN ISO 12100:2010 + Corrigendum to DIN EN ISO 12100:2011-03

2014 / 35 / EC Low Voltage Directive

EN 61010-1:2010 + A1:2019

+ A1:2019/AC:2019

Safety requirements for electrical equipment for measurement, control and laboratory use control and laboratory equipment – Part 1: General requirements. Requirements [IEC 61010-1:2010 + Cor.:2011]

2014 / 30 / EC EMV Directive

EN 61000-3-2:2014

Electromagnetic compatibility (EMC) – Part 3-2: Limits – Limits for harmonic current emissions (equipment Input current ≤ 16 A per phase) (IEC 61000-3-2:2014); German version EN 61000-3-2:2014

EN 61000-3-3:2013

Electromagnetic compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current ≤ 16 A per phase and not subject to conditional connection (IEC 61000-3-3:2013); German version EN 61000-3-3:2013

EN 61326-1:2013

Electromagnetic compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current ≤ 16 A per phase and not subject to conditional connection (IEC 61000-3-3:2013); German version EN 61000-3-3:2013

Dreve Dentamid GmbH • Max-Planck-Straße 31 • 59423 Unna/Germany • www.dreve.com

4. Safety instructions

Attention! Please read these notes carefully before the connection and commissioning of the device. The operational safety and functioning of the device can only be ensured if both the general safety and accident prevention regulations as specified by the legislator and the safety instructions in the operating instructions are observed.

**Disconnect the mains plug before cleaning and maintenance work.
Risk of electric shock.**

In case of long and intensive use or leakage of the inert gas, the oxygen might be displaced. The unit may therefore only be used in well-ventilated environments.

- The post-curing insert must be inserted into the appliance with the handle facing forwards, any other orientation may result in damage to the appliance.
- Warning notices and other labels on the appliance must not be removed.

1. The device must only be used in accordance with the present operating instructions. We do not assume any liability for damages which result from improper use and/or faulty operation.
2. Only operate the device at a socket with a protective conductor. Never touch the plug with wet hands.
3. The device must be placed on a stable and level surface which is sufficiently stable and loadable for the operating weight of approximately 9.5 kg.
4. Do not insert any objects into the device apart from the polymerization chamber.
5. Do not lock the operating and safety elements of the device.
6. Please check the device to ensure its proper condition and operating safety before each operation. If its condition is not flawless, the device must not be used and must be labelled accordingly.
7. Before the cleaning and maintenance of the device or before replacing parts, the power plug must be pulled out.
8. Signs and stickers must always be kept in an easy-to-read condition and must not be removed.

9. The opening of the device and repairs may only be performed by authorized professionals.

10. Only use nitrogen as protective gas. Do not use flammable or toxic gases under any circumstances. Please observe the general safety and accident prevention regulations when handling nitrogen.

11. Only use approved gas pressure bottles.

12. When connecting the gas pressure bottle to the device, be sure to check the tightness of the feed tube. Please make sure that the gas pressure bottle cannot tip over.

13. Never connect the gas bottle without using a pressure reducer; the device may be operated with a maximum pressure of 5 bar (70 psi).

14. Please close the gas supply, when the device is not used. Do not operate the device unattended.

15. Be sure to provide for sufficient ventilation in the work area.

16. Do not keep any highly inflammable substances within the immediate vicinity of the device.

17. Authorized operators: The operator of the device must make the operating instructions accessible to the operator and must also make sure that s/he has read and understood them. Only then may the operator put the device into operation.

18. To prevent the ingress of water into the device (e. g. splash water), the device should be installed in a dry environment.

19. Please disconnect the device from the mains when it is not used.

20. To avoid the risk of getting caught, only touch the lid at the handle.

21. Do not operate the device if the housing of the device is not closed and screwed on.

22. Unauthorized modifications and changes are inadmissible for safety reasons.

23. The prescribed operating and maintenance conditions specified in these operating instructions must be complied with. Please observe the general accident prevention regulations when working with the PCU LED N₂.

5. Commissioning

Please ensure before the commissioning that no damage has occurred during transport. Check the housing for external damages and control the glass plates on the floor and at the lid of the polymerization chamber. The lid is secured for transport reasons. The transport locks must be removed before the commissioning of the device. Please ensure that the indicated voltage on the type plate corresponds to the voltage of the power source. Insert the power cable into the appliance socket and create a network connection with the connection cable. Attach the protective gas supply to the connection of the pressure reducer of the gas pressure bottle with the use

of the hose clamp. Establish a connection to the device with the quick-lock coupling. Additional warning about the max. pressure of the shielding gas.

Scope of delivery:

1x PCU LED N₂

1x Accessory package consisting of:

3x Power cable

Only connect the cable with the plug that matches the local socket system.

1x Operating instructions

1x Gas connection line kit including quick coupling and 2 hose clamps

1x Post curing insert


6. Information on the use of the protective gas

The use of the protective gas allows for the completion of the polymerization in an atmosphere without oxygen. The protective gas replaces the air in the polymerization chamber and thus reduces the oxygen content in the polymerization space of the PCU LED N₂, which leads to the hardening of components without an inhibition layer. The protective gas supply is automatically controlled when the polymerization process starts. The flow rate is preset; the adjustment can be changed at the control needle valve at the back of the device (please observe the safety instructions).


The outlet pressure of the bottle can be mounted to the pressure reducer of the gas pressure bottle so that the flow is changed without changing the control needle valve. The gas flow is factory-set. The use of a flow meter is recommended for monitoring the gas flow during operation. (Please refer to the instructions for use of the relevant material for exact information on type and flow rate.) The connection pressure must be set to 5 bars.

7. Functional elements

Overview

Rear view

8.1. Operation: Starting the device

The software of the PCU LED N₂ is subdivided into two levels: The first level is the user level, where the operator can select and confirm different options. The second level is the process owner level, where the process owner can perform different device configurations.


(1) Switching on the power plug at the back of the device brings you directly to the user level.

(2) Switching on the power switch at the back of the device and the simultaneous pressing of the operating button brings you to the process level.


8.2. Functions of the rotary button

Selection

Confirmation

Selection of the sub-menu

PCU LED N₂

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8.3. Operation: User level

Start screen after switching on the device


Welcome screen

Firmware version

Selection of the material parameters


If a parameter set is selected, it will be highlighted on the display. Additional confirmation brings you to the selected parameter-set in the status window.

Status window


Additional confirmation starts the exposure process.

- 1 – Programme name
- 2 – Return button
- 3 – Exposure time
- 4 – Pressure (mbar)
- 5 – Exposure icon (ON/OFF)
- 6 – Exposure atmosphere (air/vacuum/nitrogen)

Ventilation of the exposure room


Wait until the vacuum is released so that e. g. the sealing for the exposure chamber is not destroyed!

Final process screen

After the exposure process, the result is shown together with the used parameter sets.

8.4. Operation: Process level


Access to the configuration level

Switching on the power switch at the rear of the device and simultaneously pressing the control button takes you to the process level. (See chapter 8.1)

Entry of the access code 217 (Status 03 / 2017)

Please confirm the code after entry. The access code cannot be freely selected and can also change in the case of a firm-ware update.

Configuration level in overview

1 – Orientation line:
A vertically oriented narrow bar at the left and the right margin of the selected window is provided for guidance purposes


Selection of the process parameters

10 programme blocks available for programming

Green – activated for the user level


White – not activated for the user level

Setting the process parameters – press the rotary button > 5 seconds

Setting the process parameters


- 1 – Return button
- 2 – Programme name
(freely definable by turning the rotary button)
- 3 – Intensity setting for LED sources
(freely definable by turning the rotary button)
- 4 – Exposure medium
 - 4.1 **Green** – nitrogen active
 - 4.2 **Green** – vacuum active
 - 4.1 and 4.2 **not marked** – air active
- 5 – Exposure time

Setting the LED intensities


The intensity can be shown in percent per panel. There are 3 LED spots on each panel.

Attention: If the predetermined values are modified for the appropriate material, the final material properties indicated by the manufacturer can be changed.


Checking the vacuum function

The vacuum function is checked on the basis of a special vacuum value being reached after a defined time. The result is shown on this screen after the test.

Time setting

Setting the date and time for the time stamp of the process recording which is shown, inter alia, in the "final screen" in the user level.

Factory setting

When reverting back to the factory settings, the individual parameter settings are deleted and the device is reset to the delivery state.

8.5. Placing the components

With the nitrogen version, the supplied post-curing insert is placed in the exposure chamber and the components

are positioned on the glass plate of the post-curing insert.

9. Maintenance

Maintenance work to be performed by the user

Maintenance unit	Maintenance frequency
Clean the glass panes in the bottom and the lid of the polymerization chamber with isopropyl alcohol	daily*
Before commissioning check whether all connections of the device are positioned correctly: <ul style="list-style-type: none"> • gas connection • mains connection 	daily*
Cleaning of the device with a mild detergent	weekly*
Check safety functions at regular intervals: <ul style="list-style-type: none"> • Opening the lid while the programme is running: UV LEDs must go out immediately and the programme must stop 	weekly*

Maintenance work to be performed by the service commissioner

Maintenance unit	Maintenance frequency
Check whether the LED panels must be or changed	weekly* with program LED-Test*
Check whether the LED panels must be calibrated	monthly*
Check whether the gas valve works reasonably	weekly*
Check the lid closing function	weekly*
Check whether all seals are still functional and/or fit correctly, whether the lid closes correctly	monthly*
Check whether all seals are still functional and fit correctly and check for leaks with the vacuum test programme	monthly*

* Recommendation, actual need depends on use intensity and other requirements.

PCU LED N₂

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10. Monitoring

The different process parameters are recorded by means of the sensors used in the PCU LED N₂. These are shown in the final screen in the user level.

Explanation of the parametric representation in the final screen in the user level

Presentation	Meaning
216	Sequence number (number of exposure jobs)
01/01/2000	Date
01:17:35	Time
040	LED – intensity of the upper panel in %
040	LED – intensity of the lower panel in %
0800	Exposure time [place 1 and 2: minutes; place 3 and 4: seconds]
1	Vacuum ON [1 – ON; 0 – OFF]
0	Nitrogen OFF [1 – ON; 0 – OFF]

Breakdown of the error codes in the final screen in the user level

Code	Status	Background display
00	Accurate work flow	Green
01	Pump runs during the exposure process	Red
02	Pump cannot draw any vacuum	Red
03	Lid opened	Red
04	No nitrogen flow	Red
05	Run time of program = 0	Red
06	LED error in program LEDTest	Red

11. Troubleshooting plan


Code	Defect	Cause	Troubleshooting
00	Accurate work flow	Not applicable	Not applicable
01	Pump is running during the exposure process	Cover gasket of exposure chamber leaky	Replacement of the lid sealing
		Hoses of the vacuum system leaky	Replacement of the hoses

Code	Defect	Cause	Troubleshooting
02	Pump cannot draw any vacuum Cutout after 150 seconds	Cover gasket for exposure chamber leaky Hoses of the vacuum system leaky Vacuum pump defective	Replacement of the lid sealing Replacement of the hoses Exchange of the vacuum pump
03	Lid is opened LED do not work	Lid not closed Lid control (micro switch) defective LED defective	Close lid Exchange lid control (micro switch) Exchange of the LED
04	No nitrogen flow	Nitrogen valve at the gas cylinder not opened. Warning: The device is only to be operated with a maximum pressure of 5 bar (70 psi)	Open nitrogen valve gas cylinder not opened
05	Run time of program=0	No run time of program was set	Set a run time > 0
06	LED error in program LEDTest	LED defective	LED has to be changed by trained personnel

Further

Device does not show any function	Network connection not correct Main switch off Main switch defective	Control network connection Switch on main switch Exchange main switch
-----------------------------------	--	---

12. Exploded view

13. Spare parts list

REF

1	Button cell	51090CR2032
2	PCU LED N ₂ control	521022
3	Power supply	521023
4	Glass plate	501024
5	Mirror plate, rear side lower part	521025
6	Mirror plate, rear side cover	521026
7	Mirror plate, sideways right	521027
8	Mirror plate, sideways left	521028
9	Rotary button	521029
10	LED panel	521030
11	Mirror plate, front side	521031
12	Adaptor for rotary button	521044
13	Vacuum seal	521050
14	Voilframe - display cap	521057
15	Cap for measuring feed-through	521063
16	Vacuum pump	521064KM32
17	UK power supply angled	51284
18	Plate	521065
19	Cable harness PCU LED N ₂ with USB bushing	521066
20	Angle piece	521075
21	Compression spring	521078
22	Cold device power cable	521512
23	Locking screw	521083
24	Magnetic valve	521560
25	Gas compression spring	521565
26	Encoder	521566
27	Power input + filter	521568
28	Fabric tube	2004
29	Microfuse	50908
30	Choke valve	50937
31	Sealing ring	51012
32	Mecanyl pipe	51030
33	Threaded connection	51050
34	Enclosure stand	51079
35	USA power supply straight	51283
36	Glass plate for PCU LED N ₂ insert	521092
37	Frame of PCU LED N ₂ insert	521091

PCU LED N₂

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14. General information

Repair and maintenance work at the PCU LED N₂ may only be performed by specialists. The power supply to the device must be disconnected in any case. The company Dreve will be responsible for the safety, reliability and performance of the device, if

- the device is only used for the appropriate application purpose.
- the operation is only carried out by professionally trained personnel.

- the device is used in accordance with these operating instructions.
- extensions, new settings, changes or repairs are only performed by authorised persons.

Since we at Dreve are constantly working on the further development of our products, we cannot exclude the possibility of technical changes.

14. Decommissioning of the product

The DIRECTIVE 2012/19/EU OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 4 July 2012 on waste electrical and electronic equipment (WEEE) lays down measures to protect the environment and human health. Waste electrical and electronic equipment marked as such may not be disposed of as unsorted municipal waste within the European Union. For detailed information regarding the professional disposal of WEEE please contact Dreve or your distributor.

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