Bedienungsanleitung · Working instructions
Mode d'emploi · Instrucciones de trabajo
Istruzioni d'uso · Gebruiksaanwijzing

Polymax 1, 1 / 120°C, 3, 5

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Pressure polymerisation unit

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Contents (EN)

1.	Description of unit	30
2.	Technical data	30
3.	Declaration of conformity	
3.1	Polymax 1	32
3.2	Polymax 1 / 120 °C	33
3.3	Polymax 3	34
3.4	Polymax 5	35
4.	Safety advice	36
5.	Initiation / scope of delivery	38
6.	Keyboard unit	39
7.	Operation Polymax 1, 3 and 5	40
7.1	Adjustment of water temperature	40
7.2	Adjustment of polymerisation time	40
7.3	Polymerisation	40
7.4	The query of the present temperature	41
7.5	Pressure adjustment	42
8.	Operation Polymax 1 / 120 °C	42
9.	Abort / termination of programme	43
10.	Change between °C and °F	43
11.	Reactivation of temperature monitor	44
12.	Additional fields of application	44
13.	Care and maintenance	45
14.	Exchange of fuses	46
15.	General notes	47
16.	Guarantee	47
17.	Trouble shooting	48
18.	Spare part list	50
19.	Explosion sketch	
19.1	Polymax 1	134
19.2	Polymax 1 / 120 °C	135
19.3	Polymax 3	136
19.4	Polymax 5	137
20.	Operation diagram	
20.1	Polymax 1, 3	138
20.2	Polymax 1 / 120 °C	139
20.3	Polymax 5	140
21.	Pneumatic plan	
21.1	Polymax 1, 3	141
21.2	Polymax 1 / 120 °C	142
21.3	Polymax 5	143

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Polymax 1, 1 / 120 °C, 3 and 5

1. Description of unit

Polymax units are machines for pressure-polymerisation of powder / liquid resins which polymerise in a tempera-

ture area from 30 °C to 95 °C (120 °C in case of Polymax 1 / 120 °C) under pressure within 90 minutes.

2. Technical data

Polymax 1

H x W x D	230 x 290 x 330 mm
Empty weight	10 kg
Power consumption	450 W
Min. Operating pressure	3 bar
Max. Operating pressure	5.2 bar
Max. Operating temperature	95 °C
Pressure pot	Ø 160 x 110 mm

Polymax 1 / 120 °C

H x W x D	230 x 290 x 330 mm
Empty weight	10 kg
Power consumption	450 W
Min. Operating pressure	3 bar
Max. Operating pressure	4.2 bar
Max. Operating temperature	120 °C
Pressure pot	Ø 160 x 110 mm

Polymax 3

H x W x D	450 x 290 x 330 mm
Empty weight	13.5 kg
Power consumption	900 W
Min. Operating pressure	3 bar
Max. Operating pressure	5.2 bar
Max. Operating temperature	95 °C
Pressure pot	Ø 160 x 300 mm

Polymax 5

H x W x D	310 x 340 x 400 mm
Empty weight	15.6 kg
Power consumption	900 W
Min. Operating pressure	3 bar
Max. Operating pressure	6 bar
Max. Operating temperature	95 °C
Pressure Pot	Ø 220 x 174 mm

Applicable for all Polymax units:

Supply voltage	230 V / 50 Hz 220 V / 60 Hz 115 V / 60 Hz (Voltage / frequency see type plate)
----------------	---

Min. compressor pressure	4 bar
Max. compressor pressure	10 bar

deutsch

english

français

español

italiano

nederlands

3. Declaration of conformity

3.1 Polymax 1

Low Voltage Directive (2014/35/EU).

We herewith declare that the machine described below is in its conception and design and in the shape delivered by us in accordance with the fundamental requirement for safety and health as prescribed in the applicable EEC-guidelines.

In case of any change or modification of the machine not authorised by us this declaration becomes invalid.

Denomination of the machine: Polymax 1/95°C
D3429/D3429A/D34291/D34291A

Type of machine: Pressure polymerisation unit

Applicable EC directives:
2014/35/EU Low Voltage Directive
2014/30/EU EMC Directive

Applied harmonised standards:

2014/35/EU	Low Voltage Directive
EN 61010-1:2010	Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1; General requirements (IEC 61010-1:2010 + Cor.: 2011); German version EN 61010-1:2010.
EN 61010-2-010:2014	Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 2-010: Particular requirements for laboratory equipment for the heating of materials (IEC 61010-2-010:2014); German version EN 61010-2-010:2014.
2014/30/EU	EMC Directive
EN 61000-3-2:2014	Electromagnetic compatibility (EMC) – Part 3-2: Limits – Limits for harmonic current emissions (equipment input current ≤ 16 A per phase) (IEC 61000-3-2:2014); German version EN 61000-3-2:2014.
EN 61000-3-3:2013	Electromagnetic compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems. for equipment with rated current ≤ 16 A per phase and not subject to conditional connection (IEC 61000-3-3:2013); German version EN 61000-3-3:2013.
EN 61326-1:2013	Electrical equipment for measurement, control and laboratory use – EMC requirements Part 1. General requirement (IEC 61326-1:2012); German version EN 61326-1:2013 Emission according to living area, business and industrial undertakings as well as small enterprises interference stability according to industrial area.

3.2 Polymax 1/120°C

Low Voltage Directive (2014/35/EU).

We herewith declare that the machine described below is in its conception and design and in the shape delivered by us in accordance with the fundamental requirement for safety and health as prescribed in the applicable EEC-guidelines.

In case of any change or modification of the machine not authorised by us this declaration becomes invalid.

Denomination of the machine: Polymax 1/120°C
D3428/D3428A/D34281/D34281A

Type of machine: Pressure polymerisation unit

Applicable EC directives:
2014/35/EU Low Voltage Directive
2014/30/EU EMC Directive

Applied harmonised standards:

2014/35/EU	Low Voltage Directive
EN 61010-1:2010	Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1; General requirements (IEC 61010-1:2010 + Cor.: 2011); German version EN 61010-1:2010.
EN 61010-2-010:2014	Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 2-010: Particular requirements for laboratory equipment for the heating of materials (IEC 61010-2-010:2014); German version EN 61010-2-010:2014.
2014/30/EU	EMC Directive
EN 61000-3-2:2014	Electromagnetic compatibility (EMC) – Part 3-2: Limits – Limits for harmonic current emissions (equipment input current ≤ 16 A per phase) (IEC 61000-3-2:2014); German version EN 61000-3-2:2014.
EN 61000-3-3:2013	Electromagnetic compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems. for equipment with rated current ≤ 16 A per phase and not subject to conditional connection (IEC 61000-3-3:2013); German version EN 61000-3-3:2013.
EN 61326-1:2013	Electrical equipment for measurement, control and laboratory use – EMC requirements Part 1. General requirement (IEC 61326-1:2012); German version EN 61326-1:2013 Emission according to living area, business and industrial undertakings as well as small enterprises interference stability according to industrial area.

3.3 Polymax 3

Low Voltage Directive (2014/35/EU).

We herewith declare that the machine described below is in its conception and design and in the shape delivered by us in accordance with the fundamental requirement for safety and health as prescribed in the applicable EEC-guidelines.

In case of any change or modification of the machine not authorised by us this declaration becomes invalid.

Denomination of the machine: Polymax 3
D3422/D3422A/D34221/D34221A

Type of machine: Pressure polymerisation unit

Applicable EC directives:
2014/35/EU Low Voltage Directive
2014/30/EU EMC Directive

Applied harmonised standards:
2014/35/EU
EN 61010-1:2010 **Low Voltage Directive**
Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1; General requirements (IEC 61010-1:2010 + Cor.: 2011); German version EN 61010-1:2010.
EN 61010-2-010:2014 Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 2-010: Particular requirements for laboratory equipment for the heating of materials (IEC 61010-2-010:2014); German version EN 61010-2-010:2014.
2014/30/EU
EN 61000-3-2:2014 **EMC Directive**
Electromagnetic compatibility (EMC) – Part 3-2: Limits – Limits for harmonic current emissions (equipment input current ≤ 16 A per phase) (IEC 61000-3-2:2014); German version EN 61000-3-2:2014.
EN 61000-3-3:2013 Electromagnetic compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems. for equipment with rated current ≤ 16 A per phase and not subject to conditional connection (IEC 61000-3-3:2013); German version EN 61000-3-3:2013.
EN 61326-1:2013 Electrical equipment for measurement, control and laboratory use – EMC requirements Part 1. General requirement (IEC 61326-1:2012); German version EN 61326-1:2013 Emission according to living area, business and industrial undertakings as well as small enterprises interference stability according to industrial area.

3.4 Polymax 5

Low Voltage Directive (2014/35/EU).

We herewith declare that the machine described below is in its conception and design and in the shape delivered by us in accordance with the fundamental requirement for safety and health as prescribed in the applicable EEC-guidelines.

In case of any change or modification of the machine not authorised by us this declaration becomes invalid.

Denomination of the machine: Polymax 5
D3423/D3423A/D34231/D34231A

Type of machine: Pressure polymerisation unit

Applicable EC directives:
2014/35/EU Low Voltage Directive
2014/30/EU EMC Directive

Applied harmonised standards:
2014/35/EU
EN 61010-1:2010 **Low Voltage Directive**
Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1; General requirements (IEC 61010-1:2010 + Cor.: 2011); German version EN 61010-1:2010.
EN 61010-2-010:2014 Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 2-010: Particular requirements for laboratory equipment for the heating of materials (IEC 61010-2-010:2014); German version EN 61010-2-010:2014.
2014/30/EU
EN 61000-3-2:2014 **EMC Directive**
Electromagnetic compatibility (EMC) – Part 3-2: Limits – Limits for harmonic current emissions (equipment input current ≤ 16 A per phase) (IEC 61000-3-2:2014); German version EN 61000-3-2:2014.
EN 61000-3-3:2013 Electromagnetic compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems. for equipment with rated current ≤ 16 A per phase and not subject to conditional connection (IEC 61000-3-3:2013); German version EN 61000-3-3:2013.
EN 61326-1:2013 Electrical equipment for measurement, control and laboratory use – EMC requirements Part 1. General requirement (IEC 61326-1:2012); German version EN 61326-1:2013 Emission according to living area, business and industrial undertakings as well as small enterprises interference stability according to industrial area.

4. Safety advice

Attention! Read the following notes before connecting and putting the unit into service. The operating service and the function of the unit can only be guaranteed for, if the general safety and accident prevention instructions of the legislation as well as the safety instructions of the working instruction are observed.

1. This unit must only be used according to the present working instructions. We cannot be held liable for damages which have been caused by improper treatment resp. faulty manners of operation.

2. The unit has to be placed onto a stable and even surface, which is suitable for the operation weight of the unit.

3. The indicated voltage on the type plate has to meet the voltage of the power source.

4. Connect the unit only to a socket with protective conductor. Never touch the plug with wet hands.

5. The pressure of the compressed air line must not exceed 10 bar.

6. The heating of the unit may only be switched on, if a sufficient quantity of water has been poured.

7. Do not lock the operation elements of the unit.

8. Plates and labels have to be maintained in a well legible condition and must not be removed.

9. Authorised operators: The operating company of the machine has to hold

the working instructions at the operators' disposal and has to make sure that authorised operators have read and understood the manual. Only in that case the operators may operate the unit.

10. The unit has to be examined before operation in view of proper condition. If the unit is not in proper condition, it must not be used and has to be marked accordingly.

11. In case of longer works at maximum operation temperature the lid of the unit may heat up to 70 °C (100 °C in case of the Polymax 1 ° / 120 °C). Therefore it must not be touched – not even for a short moment.

12. At all processing steps which are executed under pressure please wait with the opening of the lid until the compressed air has completely escaped from the pressure pot.

13. Please do not bend over the unit while opening the lid. Due to the escape of hot steam there is the danger of burns.

14. Before cleaning and maintenance of the unit or the changing of parts the plug has to be disconnected.

15. Do not intrude any objects into the unit.

16. The opening of the unit and repairs are only to be carried out by accordingly trained experts.

17. Arbitrary rebuilding or modifications are not permitted due to safety reasons.

18. Only those appliances and spare parts may be used which are released by the manufacturer. We are not responsible for damages caused by the use of outside parts.

19. If the unit is not used for a longer time, please disconnect the mains.

20. According to the pressure tank regulation the Polymax is a pressure tank of the test group I. These pressure tanks do not need to be taken a test again. But we recommend to check the pressure pot every 3 years by the manufacturer.

21. The required operation and maintenance conditions mentioned in these working instructions are to be complied with imperatively. When working with the Polymax the general accident prevention regulations have to be observed.

5. Initiation

Before initiation make sure that the voltage indicated on the type plate meets the voltage of the power source.

When selecting a suitable place for the machine, please make sure that the storage space is even, resistant to splash water as well as sufficiently stable for a minimum weight of 21 kg.

Connect the unit to an air compressor. Push the pressure hose (included in the delivery) into the compressed air connection (25) until it is locked in place. The hose does not need any further fastening. To release the pressure hose push the blue ring upwards and draw the hose off by pulling downwards.

The pressure of the compressed air line must not exceed 10 bars.


Put the power cord into the female connector (15) and create a connection with the Euro / USA plug. Open the lid (2) by means of the handle (1) and fill the unit with water. (Maximum water level: 30 mm below the upper edge).

Important: Before first initiation please read the working instructions thoroughly. In case of comprehension problems contact your responsible dealer or Drevo, if necessary.

Scope of delivery

- 1 Polymax
- 1 power cord
- 1 air connection hose incl. hose clamps and fast coupling
- 1 water drain hose incl. coupling
- 1 working instructions in 6 languages

6. Keyboard unit

No. Description

- a Manometer
- b Adjustment knob
- c LED – heating
- d Temperature indication display
- e Temperature switch
- f Time indication display
- g LED – polymerisation
- h only Polymax 1 / 120°C – lid opening control

7. Operation Polymax 1, 3 and 5

7.1 Adjustment of water temperature

Switch on the unit at the main switch (14). At the first startup the display for temperature (d) shows 45°C, which has been adjusted by the factory. The display for time (f) shows 0. To get to the adjustment area for polymerisation temperature press temperature switch (e) for more than 3 seconds. The temperature indication display (d) as well as the LED for heating (c) are blinking. By turning the knob (b) to the right or left every temperature between 5 and 95°C (41 and 203°F) can be adjusted in intervals of 5°C (9°F). The adjusted value will be confirmed by pushing the temperature switch (e) again.

On the right besides the display (d) an LED indication (c) shows the current heating process.

The adjusted value will be indicated blinking in the display (d) until the pre-adjusted temperature will be attained in the pressure pot. If the temperature of the water bath meets the adjusted temperature, the LED indication (c) disappears and the blinking in the display indication (d) stops.

The water temperature in the polymerisation pot is measured constantly electronically by means of a temperature sensor. In case of a deviation of 2°C the heating switches on again automatically until the set temperature is reached.

Tip: For a quick and energy efficient heating of the water bath the lid of the polymerisation pot (2) is to be closed during this phase.

7.2 Adjustment of polymerisation time

Push the adjustment knob (b) one time and adjust the desired polymerisation time by turning the knob to the

right/left. Adjustment can be realised in 5-minute intervals up to a maximum value of 90 minutes.

7.3 Polymerisation

Open the lid (2) and put the objects to be polymerised into the warm water bath. Handle sharp-edged objects (e. g. metal flasks or articulators) very carefully, as scratches in the surface coating of the polymerisation pot may lead

to calcification and pitting corrosion in the long run.

Control the correct position of the O-Ring (7) in the groove and correct it, if necessary by pressing onto it with the

fingers. Close the lid (2) until the block and push adjustment knob (b) one time in order to pressurize. The lid will seal automatically and cannot be opened during active pressurizing.

Important: Polymerisation pot and lid can become very hot in case of longer use at high temperatures. Therefore contact to metal parts is to be avoided. The lid should only be opened and closed by means of the isolated handle (1)!

On the right besides the time indication display (f) an LED indication (g) shows the active pressure polymerisation. The manometer (a) shows the pressure within the polymerisation pot and the pre-adjusted time runs backwards in one-minute steps towards 0.

After the polymerisation time has run out the overpressure will be released automatically. Only open the lid after complete ventilation at the handle (1) and remove the objects from the water bath. The temperature adjustment is saved for future polymerisation processes.

Tip: Polymerisation can also be started, if the final temperature has not been reached.

By pushing the adjustment knob (b) for approx. 3 seconds a running polymerisation programme can be cancelled (see point 9).

7.4 The query of the present temperature

To ensure the right condition for polymerisation the active temperature can be checked at all times. For this press the temperatures switch (e) and keep pressing. During this action, but at the most for 3 seconds the actual value is shown in the display (d). The proceeding programs (heating, recooling or polymerisation) will not be interrupted due to this activity.

If pressing the temperature switch (e) less than 3 seconds, the electronic turns back to the initial position and show the adjusted set temperature.

If pressing the temperature switch longer than 3 seconds, the Polymax unit will vary to the adjustment area of temperature automatically (see point 7.1).

7.5 Pressure adjustment

The pressure inside the polymerisation pot can be adjusted infinitely variable between 3 bar and the maximum value of the unit in question. For this purpose pull out the reduction valve (24) at the rear of the unit by means of the knob and turn it anti-clockwise until there is no resistance anymore (valve closed).

Start pressure-polymerisation (see 7.3) and turn the pressure reduction valve clockwise. The manometer (a) shows the adjusted pressure. When attaining the desired pressure, let knob latch again.

Important: Do not adjust the pressure below 2.5 bar! Lower pressure may lead to leakage of the pneumatic sealing between polymerisation pot and lid.

8. Operation Polymax 1 / 120 °C

Besides the general temperature area of 0–95 °C the Polymax 1 / 120 °C also permits to increase the water temperature to 120 °C. As the boiling point of water is 100 °C, higher temperatures may only be reached under pressure. In order to avoid that hot water vapour of 120 °C escapes from the pressure pot and in order to reduce the risk of scalding when opening the lid, the Polymax 1 / 120 °C is equipped with a safety lock which only allows the lid to open when the operating pressure in the unit has fallen and the temperature is lower than 95 °C (203 °F). This is also applicable in case of power outage.

This means that in the application process the pressure pot has to remain closed when the unit is turned off and pressure is exerted.

When the unit is turned on at the main switch (14) the pressure pot exhausts automatically and the lid (2) can be opened.

Adjustment of temperature and polymerisation time is to be executed as described from point 7 et seqq. If a temperature > 95 °C is required, 95 °C should be selected as pre-warming temperature; the final temperature is to be adjusted directly before the polymerisation process.

According to the temperature of the water bath the objects have to be inserted and taken out with suitable instruments!

The time until the final temperature is reached has to be added to the general polymerisation time. This can differ between 10 and 20 minutes, depending on the water filling level. As long as the value in the temperature indication display (d) blinks, the final temperature has not been reached. The present temperature of the water bath can be checked at all times by pressing the temperature switch (e), (see point 7.4).

When the time has passed the unit ventilates automatically. The lid of the Polymax 1 / 120 °C can, however, only be opened, if the temperature of the water bath has fallen below 95 °C.

If the pressure pot has got a temperature of more than 95 °C after the time has run out, it is cooled down by cold air until the temperature has fallen

below 95 °C. During this time the cooling air is released by means of the ventilation valve. When operating at 120 °C the cool down phase can take 5–8 minutes from time expiration to the opening of the lid.

The possibility of opening the lid is indicated by the green control lamp (h). Open the lid and remove the objects.

9. Abort / termination of programme

The polymerisation process can be terminated prematurely. For this purpose push the adjustment knob (b) for 3 seconds. The Polymax ventilates automatically, the display for time shows 0, the pre-adjusted temperature is maintained.

For following polymerisation processes with identical temperature continue as described from point 7.2 on.

10. Change between °C and °F

At delivery the Polymax units are adjusted to °C. In order to change to °F push the temperature switch (e) and the main switch (14) at the same time. The temperature indication display (d) shows L (which stands for °C), the time indication display (f) shows the maximum value of 95. By turning the adjustment knob (b) it is possible to change between L (°C) and H (°F). When the machine is adjusted to H the time indication display shows the value 203.

Exception Polymax (1 / 120 °C): The time indication display (f) shows the maximum values 120 °C or 248 °F. Push the adjustment knob (b) one time in order to confirm the selected temperature denomination. The time and temperature indication displays (d + f) indicate 0. Then go on as described from point 7 on.

11. Reactivation of temperature monitor

In order to avoid an overheating and consequently the destruction of important Polymax components, all units have been equipped with a temperature control (11). This part is situated from outside at the bottom of the polymerisation pot and interrupts the circuit of the heating in case of too much temperature development. In order to reactivate the temperature monitor disconnect the unit from the mains and let it cool down. Drain the water completely through the drain tap (23). Lay the Poly-

max on its side and use a soft tissue as underlay, if necessary.

Unscrew the perforated metal plate at the bottom of the unit and push the red reset knob of the temperature monitor until it latches completely. Insert perforated metal plate again and tighten it. Put the Polymax into normal position, close drain tap (23), fill the unit with water and connect to the mains by means of the power cord. Then go on as described from point 7 on.

12. Additional fields of application

All Polymax units are suitable for pressure compression of embedding masses, plasters and silicones. For this purpose drain water through the drain tap (23), close the tap again and adjust the temperature to 0 °C (32 °F) (see 7.1).

Important: Heating up the polymerisation pot without water may lead to damage of the heating cartridge.

Open the lid (2) at the handle (1) and put the objects into the dry polymerisation pot.

Tip: In order to avoid impurities and/or unevenness at the bottom of the pot caused by plaster or embedding masses we recommend an underlay of paper or the like.

Close the lid of the pot up to the block and start pressure-polymerisation (see 7.3).

13. Care and maintenance

Disconnect the unit from mains supply before care and maintenance! The water inside the pressure pot has to be changed regularly (at least once a week), because the lime particulates harm the surface of the pressure pot.

Polymax 1, 3 and 5

For the cleaning of the pressure pot open the drain tap (23) at the rear of the unit and let the water flow out. Remove dirt thoroughly.


Attention: Do not use scouring agents, steel-wool or sponge scourers.

Clean the interior pot with fresh water, a neutral washing-up liquid and a soft scarf or sponge. Afterwards let the water run out. Close the drain tap (23) and fill with fresh water. The O-ring should be greased regularly with silicone grease (35 g tube, [REF](#) D3350).

14. Exchange of fuses

Disconnect the unit from mains supply. The fuse carrier is situated at the rear of the unit, integrated in the mains supply.



In order to control or exchange the fuses push the push catch of the carrier, keep it pressed and remove the carrier from the unit.

Pull the defective fuse out of the clamp carrier applying slight traction and replace it by a new, equivalent one.

Important: Only use the fuse strength indicated on the rear.

After exchanging the fuse insert the fuse carrier into the housing applying slight pressure until the fixation latches noticeably.

15. General notes

Repairs on the Polymax units may only be executed by experts. Disconnect the unit from mains supply at any rate. Dreve company is only responsible for the safety, reliability and performance of the unit, if:

- the unit is exclusively used for the corresponding field of application.
- operation is only carried out by accordingly trained personnel.

16. Guarantee

The guarantee time amounts to two years after the date of delivery. Guarantee is only accepted in case of material or fabrication defects. The guarantee obligation does not refer to natural wear and tear or to damage occurring after the transfer of risk due to improper or unsuitable handling, excessive stressing, unsuitable operating material, and chemical, electrochemical or electrical influences of a nature not provided for according to the contract. On improper alterations, repairs or maintenance or breach of seals carried out by the purchaser or third parties our warranty obligations and liability for consequences caused thereby are rendered void.

- the unit is operated in accordance with these working instructions.
- extensions, new adjustments, changes or repairs are executed by persons authorised by the company.

As we are currently developing our products, we reserve the right for technical changes.

Liability

The right of the purchaser to bring actions arising out of defects shall in all cases be in lapse one year after the date of delivery. The warranty obligation does not refer to natural wear and tear or to damage occurring after the transfer of risk due to improper or unsuitable handling, excessive stressing, unsuitable operating material, and chemical, electrochemical or electrical influences of a nature not provided for according to the contract. On improper alterations, repairs or maintenance or breach of seals carried out by the purchaser or third parties our warranty obligations and liability for consequences caused thereby are rendered void.

17. Trouble shooting

Problem	Cause	Aid
Unit does not show any function	Mains supply is not correct	Control the mains supply
	Main switch is off	Turn on the main switch
	Fuses in fuse carrier are defective	Exchange fuses; in case of unsuccessful exchange look for other causes
Unit does not heat	Temperature regulator shows 0 °C	Adjust the temperature regulator
	Heating is defective	Heating has to be exchanged by an expert
	Temperature monitor has switched off the heating	Push reset button at temperature monitor, temperature monitor is situated underneath the pressure pot (perforated metal plate) or Temperature monitor is defective and has to be exchanged
Unit does not create pressure	Unit does not have compressed air	Control compressed air connection
	Timer is not turned on	Turn on the timer
	Magnetic valve is defective	Magnetic valve has to be exchanged by an expert

Problem	Cause	Aid
Unit loses compressed air under the lid	Lid O-ring is polluted	Clean and grease the lid O-ring
	Lid O-ring is defective	Renew the lid O-ring
Unit blows air out of the security valve	Input pressure is too high	Reduce the input pressure

Addition only for Polymax 1 / 120 °C

Problem	Cause	Aid
Lid cannot be opened	Temperature in the pressure tank is too high	Unit has to be cooled down
	Temperature monitor type: 1822 L 20-4 98 °C (temperature monitor) is defective	Exchange the temperature monitor
	Magnetic valve is defective	Magnetic valve has to be exchanged by an expert
Unit does not cool down automatically	Magnetic valve is defective	Magnetic valve has to be exchanged by an expert
Unit does not keep the holding temperature	Temperature monitor type: 18922 L 20-4 98 °C (temperature monitor) is defective	Exchange the temperature monitor

18. Spare part list

No.	Description	REF
Polymax 1 230V/50 Hz		
1	Handle	D50275
2	Lid	D50479KPL
3	Electronics / digital displays	55200
4	Counter bearing	D50417KPL
5	Seal housing 0,6 mm	D6550617
6	Housing red	D6552571R
	Housing silver	D6552571S
7	O-ring 164 x 4 mm	51031
8	Pressure screw connection manometer	D6550183
9	Pressure screw	D6550183
10	Polymerisation pot PM 1	D55551PTFE
11	Temperature monitor	55743
13	Bottom plate PM 1-3	510975
14	Main switch	50848
15	Mains supply/fuse carrier	51208
16	Pressure control valve	50792
17	Manometer	50702
18	Turning knob	51322
19	Lid block	50305
20	Spacer lid bolt mm: 3/51307; 3,5/50076; 4/51331; 4,2/50072	50230
21	Rear panel of unit PM 1, 1/120°C	510977
22	Magnetic valve complete / conversion kit	D65741PM1
23	Drain tap	D6550744
24	Pressure reduction valve	D6550726
25	Compressed air supply	50340
26	Compressed air outlet	50171
27	Metal sheet cover	510964
Without drawing no.		
	Power cord EU	51021
	Power cord UK	51284
	Power cord USA / JP	51283
	Heating cartridge 200 W	510542
	Water drain hose	D20041
	Adhesive display foil PM 1	510209
	Fine fuse 4 AT	51172
	Valve plug with rectifier from SN 21700133 (always indicate the serial number of the device, otherwise conversion kit complete)	507451

No.	Description	REF
from SN 21700133		
	Angular screw-in connection 5-1/8 at the Pressure control valve, Rear panel and pressure outlet	50673 QM 1/8-1/8 50269
	Rotary encoder	55200B
	Ahesive pad	51079
	„T“ Screw-connection manometer 5-1/8	50674 QM 1/8-1/8 50269
	Angular screw-in connection 5-1/8	
	Pressure reduction valve and Magnetic valve	50673
	Angular swivelling screw connection, Air release Magnetic valve	50170
	Temperature sensor	55744
Older devices need in addition:		
Only Polymax 1 115 V / 60 Hz		
	Heating cartridge	513022
	Fine fuse 8 AT	51061
	Magnetic valve with rectifier (Maybe screws are needed!)	D6550741110
	Valve plug with rectifier from SN 11700107 (always indicate the serial number of the device)	507451
Only Polymax 1 / 120°C 230 V / 50 Hz		
1a	Handle	D50319
12/22	Magnetic valve complete / Conversion kit	D65120PM1
7a	O-ring 164 x 4 / 120 °C	51168
11a	Temperature monitor lid	50701
11b	Temperature monitor maintaining temperature	50935
27	Signal lamp green	50853
Only Polymax 1 / 120°C 115 V / 60 Hz		
	Heating cartridge	513022
12/22	Magnetic valve complete / Conversion kit	D6555742110
Only Polymax 1 / 120°C		
6	Housing red	D6552571R120
	Housing silver	D6552571S120
Without drawing no.		
	„T“-Screw connection manometer 5-1/8	50674 QM 1/8-1/8 50269
	Angular screw-in connection- 5-1/8	
	Pressure reduction valve and Magnetic valve at the air outlet and Pressure control valve additionally with	510164 QM 1/8-1/8 50269
Only Polymax 3 230 V / 50 Hz		
6a	Housing red	D6552573R
	Housing silver	D6552573S
10a	Polymerisation pot	D55553PTFE
21a	Rear panel of unit	510954
	Heating cartridge 400 W	50737
Only Polymax 3 115 V / 60 Hz		
	Heating cartridge 400 W	50714
	Fine fuse 8 AT	51061
	Magnetic valve with rectifier (always indicate the serial number of the device. Maybe screws are needed!)	D6550741110
	Valve plug with rectifier from SN 11700107 (always indicate the serial number of the device)	507451

Polymax 1, 1 / 120°C, 3 and 5

No.	Description	REF
Polymax 5 230V/50 Hz		
1a	Handle	D50234
2a	Lid	D50466KPL
3	Electronics / digital displays	55200
4a	Counter bearing	D50467KPL
5a	Seal housing 0.75 m	D6550617
6b	Housing red	D6552575R
	Housing silver	D6552575S
7a	O-ring 226 x 4 mm	50663
8	Pressure screw connection manometer	D6550183
9	Pressure screw connection PM 5	D6550183
10b	Polymerisation pot PM 5	D55555PTFE
11	Temperature monitor	55743
13a	Bottom plate	510714
14	Main switch	50848
15	Mains supply / fuse carrier	51208
16	Pressure control valve	50792
17	Manometer	50702
18	Tuning knob	51322
19	Lid block	D6550305
20	Spacer lid bolt	50230
21b	Rear panel of unit PM 5	510715
22	Magnetic valve complete / conversion kit	D65745PM5
23	Drain tap	D6550744
24	Pressure reduction valve	D6550726
25	Compressed air supply	50340
26	Compressed air outlet	50171
27	Metal sheet cover	510716
Without drawing no.		
	Power cord EU	51021
	Power cord UK	51284
	Power cord USA / JP	51283
	Heating cartridge	510713
	Water drain hose	D20041
	Adhesive display foil PM 5	510210
	Fine fuse 4 AT	51172
	Valve plug with rectifier from SN 21700001 (always indicate the serial number of the device, otherwise conversion kit complete)	507451

No.	Description	REF
from SN 21700001		
	Angular screw-in connection 5–1/8 at pressure control valve, back panel	50673 QM 1/8–1/8 50269
	Rotary encoder	55200B
	Adhesive pad	51079
	T-screw connection manometer 5–1/8	50674 QM 1/8–1/8 50269
	Angular screw-in connection 6–1/8 at pressure reduction valve- magnetic valve and pressure outlet	50672 QM 1/8–1/8 50269
	Angular swivelling screw connection, air outlet	
	Air release magnetic valve	50267
	Temperature sensor	55744
Only Polymax 5 115 V/60 Hz		
	Heating cartridge	510722
	Fine fuse 8 AT	51061
	Magnetic valve with rectifie (always indicate the serial number of the device. Maybe screws are needed !)	D6550741110
	Valve plug with rectifier from SN 11700107 (Always indicate the serial number of the device)	507451

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Dreve Dentamid GmbH