Bedienungsanleitung · Working instructions
Mode d'emploi · Instrucciones de trabajo
Istruzioni d'uso

Dosper evo

Dosier- und Mischgerät für Silikonmassen
Dosing and mixing unit for silicone masses
Appareil de mélange et de dosage pour silicones
Dispositivo dosificador y mezclador para masas de silicona
Dosatore e miscelatore per masse di silicone

deutsch

english

français

español

italiano

Inhalt (DE)

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Contents (GB)

1.	Technical description	22
2.	Technical data	22
3.	Declaration of conformity	23
4.	Safety advice	24
5.	Initiation – Scope of delivery	25
6.	Scheme of functional components	26
7.	Handling	27
7.1	Connection of the silicone containers	27
7.1.1	Connection of 2 x 1 kg (850 ml) bottles	27
7.1.2	Connection of 2 x 6 kg (5.1 l) and 2 x 25 kg (21.2 l)	28
7.1.3	Connection of the material reservoirs	28
7.2	Connection of the Dosper evo to the electric circuit	28
7.3	Preliminaries	29
7.4	Applying the mixing tips	29
7.5	Duplication	30
7.6	Bottle-/canister change	30
8.	Deviations from the regular process	31
8.1	Mixing Head only delivers one component	31
8.2	Proportion of the components clearly deviates from a 1:1 mixing	32
8.3	The pump is in operation, the LED's are flashing, no silicone is being delivered	32
8.4	The machine shows no function	32
9.	Optional accessories	33
10.	Care and maintenance	33
11.	Replacement of the fuse	34
12.	General advice	35
13.	Guarantee	35
14.	Spare parts list	36
15.	Exploded drawing	37
16.	Wiring scheme	86-87

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
Dosper evo

1. Technical description

The Dosper evo is a mixing and dosing unit for all thin-flowing 1:1 duplicating silicones (addition-vulcanising or VPS) that are daily used in dental or orthodontic laboratories.

The machine delivers, doses and mixes the silicones fully automatic and makes

the exactly required material quantity available at the push of a button.

Quick couplings allow the connection of all original packaging units (1 kg, 6 kg, 25 kg) as well as the optionally available clear-transparent material reservoirs.

2. Technical data

Dimensions (H x W x D):

320 x 225 x 405 mm

Weight:

approx. 10 kg (without material)

Mains voltage:

220–230 V / 50–60

Hz 115 V / 60 Hz

(Voltage / Frequency see type plate)

Power consumption:

approx. 280 W

Fuse:

1.6 A

Output:

150 ml / min.

Gear motor with overheat protection

Manual adjustability of delivered quantities

Illuminated work area

Optional:

Auto-stop before complete emptying
(only when used with material reservoirs)

3. EC-Declaration of conformity

according to Low-Voltage Directive (2006/95/EG)

Herewith we declare that the machine described below in its conception and design and in the shape delivered by us is in accordance with the fundamental requirements for safety and health as prescribed in the applicable EC directives.

In the case of any change or modification of the machine, not authorized by us, this declaration becomes invalid.

Denomination of the machine:

Dosper evo
D-5050 / D-5050-A

Type of machine: Dosing and mixing unit for silicone masses

Applicable harmonised standards:

2006/95/EG Low voltage guideline
2004/108/EG EMV-guideline

Applicable harmonised standards:

2006/95/EG Low Voltage Directive

EN 61010-1:2001 Safety requirements for electrical equipment for measurement, control and laboratory use – part 1: General requirements [IEC 61010-1:2001]

2004/108/EG

EN 61326-1:2006 **EMV-guideline**
Electrical process measuring and control devices, laboratory units – EMV-requirements – part 1: General requirements (IEC 61326-1:2005) German Version 61326-1:2006
Emission according to business sector and small enterprises

EN 55014-1S

Measuring of the interference voltage according to EN 55014-1:2006 frequency range 150 kHz–30 MHz
Conducted emissions special evaluation of the interference

Dreve Dentamid GmbH
Max-Planck-Straße 31 · 59423 Unna/Germany
Tel.: +49 2303 8807-40 · Fax: +49 2303 8807-55
E-Mail: info@dreve.de · www.dreve.com

Dosper evo

4. Safety advice

Attention! Please read these instructions carefully before connecting and putting into service. The reliability and the functionality of the machine can only be ensured, if the general governmental accident prevention regulations as well as the safety information in the manual are respected.

1. The machine may only be used according to the operating instructions on hand. We accept no liability for damages that are caused by improper handling or incorrect operation.
2. The Dosper must stand on an even undersurface which has to be sufficiently stable and resilient for the operating weight of 12–60 kg (depending on the used container sizes).
3. To avoid penetration of water (e.g. splashing water), the machine should be set up in a dry place.
4. The voltage mentioned on the type plate has to correspond with the voltage of the power source.
5. The machine shall only be connected to a socket with protective conductor. Never connect or disconnect the main plug with wet or moist hands.
6. Approved operators: the machine owner must make the user instructions accessible to the operators. He also has to make sure, that any approved operator has read and understood them. Only after that the operator is allowed to put the machine into operation.
7. The Dosper must be checked for its impeccable condition and work safety before use. If it is not in an impeccable condition, the machine shall not be used and has to be marked accordingly.
8. Do not put any objects into the machine.
9. Do not lock the pushbutton of the machine.
10. When not in use, the power supply has to be disconnected.
11. The signs and stickers have to be kept readable and shall not be removed.
12. Before cleaning and maintenance or replacement of spare parts disconnect the main plug.
13. The machine may only be opened and maintained by skilled personnel.
14. Only use accessories and spare parts which have been authorized by the manufacturer. We cannot be held liable for any damage occurring from the use of unsuitable parts.

15. For safety reasons any unauthorized modification or change of the machine is prohibited.

16. The operating and maintenance instructions of this manual must abso-

lutely be followed. While working with the Dosper evo the general accident prevention regulations have to be kept to.

5. Initiation – Scope of delivery

Before taking the machine into service, make sure that the mains voltage indicated on the type plate complies with the voltage of the power source.



Put the power cord into the appliance plug socket (3) and connect the Euro/UK/USA plug to the power source.

Scope of delivery:

- 1 x Dosper evo
- 1 x power cord
- 1 x operating instructions in 5 languages
- 5 x mixing tips
- 1 x set bottle adapter

Important: Read the operating instructions carefully before initial startup. In case of uncertainties consult the responsible distributor or Dreve Dentamid GmbH directly.

6. Scheme of functional components

No. Description

- | | | | |
|----|--|----|-------------------|
| 1 | Main switch | 7 | Activation button |
| 2 | Fuse holder | 8 | LED lighting |
| 3 | Mains connection | 9 | Setting screw |
| 4a | Connector socket for autostop
(only if material reservoirs
are used) | 10 | Setting screw |
| 4b | Connector socket for autostop
(only if material reservoirs
are used) | | |
| 5a | Adapter for material container | | |
| 5b | Adapter for material container | | |
| 6 | Mixing head | | |

7. Handling of the Dosper ev0

Use high quality duplication silicones only. Unequal viscosities and distribution of filler materials in the individual components can lead to an abrasion in the pump housing and thus to a

deviation of the ideal 1:1 mixture (see point. 8.2)

We recommend the use of Dreve DUBLISIL®-Silicones.

7.1 Connection of the silicone containers

With the help of adapters all commonly used material containers as well as the optionally available material reservoirs can be connected to the Dosper ev0.

within the machine, causing major repairs which will have to be executed in our plant in Unna.

Important notice:

Please make sure that the component connections are never mixed up. On replacement of bottles or the refill of the reservoirs always ensure that the catalyst (standard colour white) as well as the base (coloured) are always connected to the same side of the machine.

Important notice:

The filling level of the containers must be controlled constantly! They shall never be emptied completely. Otherwise air might get into the system and the machine must be de-aerated in a complex procedure (see. point 8.1). A minimum level of about 2–3 cm always has to remain in the containers.

Already through minimal blending of the components the silicone hardens

7.1.1 Connection of 2 x 1 kg (850 ml) bottles

Replace the original cap of the bottle by the adapter, which is included in the delivery, and fasten it hand-tight. Note that each adapter must absolutely be used for the same silicone component again when changing bottles in the future.

attach them to the respective counterparts of the adapters (5) on the machine. Press down until they noticeably snap into place.

The bottles can now be connected to the machine with the help of the adapters. Turn the bottles upside down and

To allow a problem-free flow of the silicone mass in the Dosper ev0, and in order to avoid the development of a vacuum in the bottles, pierce tiny holes into their bottoms.

Dosper eva

7.1.2 Connection of 2 x 6 kg (5.1 l) and 2 x 25 kg (21.2 l)

Replace the original caps of the canisters by the optionally available connection hoses (see point 9) and fasten them hand-tight.

Note that each connection hose must absolutely be used for the same silicone component again when changing canisters in the future.

The other side of the connection hose can now be connected to the machine with the help of the adapters. Turn the

canisters upside down and position them above the machine (a compatible support stand is also available (see point 9)). Press the connection hoses onto the adapters (5) until they noticeably snap into place.

To allow a problem-free flow of the silicone mass in the Dosper eva, and in order to avoid the development of a vacuum in the canisters, pierce tiny holes into their bottoms.

7.1.3 Connection of the material reservoirs

The optionally available material reservoirs (see point 9) are clear-transparent containers, marked with two different colours for distinction. They are connected to the machine with the help of the pre-assembled quick couplings. An electronically controlled float lever avoids the unwanted complete emptying of the reservoirs. To activate the system connect the plug-in con-

tacts mounted to the reservoirs with the respective sockets (4a/4b) of the machine.

Take off the lids of the reservoirs and fill in the respective silicone components. Each of the two reservoirs has a capacity of up to 1.5 l (corresponds to approx. 1.8 kg). Then close with the lids.

7.2 Connection of the Dosper eva to the electric circuit

Connect the provided power cord to the respective socket of the machine (3),

then connect to the electric circuit.

7.3 Preliminaries

From the factory the Dosper evo is delivered with an inactive test liquid inside the system. This test liquid must completely be replaced by silicone before taking the machine into operation.

Turn the machine on at the main switch (1). The illuminated activation button (7) will light up.

Place an appropriate container underneath the mixing head and press the activation button (7).

The internal pump will press the test mass out of the machine and will replace it automatically by the connected silicone.

Keep the button (7) pressed until both silicone components are delivered evenly out of the mixing head (6).

The Dosper evo is now ready for use.


7.4 Applying the mixing tips

The Dosper evo is equipped for the application of mixing tips with pink closure cap (D-4970-1).

There is a groove at the front of the mixing head (6). This groove has to

match with the triangular bulge of the mixing tip.

Snap the mixing tip into the guiding and apply to the mixing head (6) by a clockwise rotation of 90°.

Dosper evo

7.5 Duplication

Place your duplicating flask under the mixing tip and press the activation button (7). The two LEDs (8) behind the mixing head will light up. On the one hand this is to control the operation of the pump, on the other hand the view in the work area is improved. As long as the silicone in the mixing tip is not cured, an arbitrary number of flasks can be filled one after the other. Leave the tip on the mixing head after the filling process, it will serve as a closure of the

mixing head. Next time the machine is going to be used just remove the old mixing tip by an anti-clockwise turn of 90° and replace it by a new one.

Important: Never activate the machine when a mixing tip with cured material is still locked to the mixing head!
In case the Dosper evo is not used for a longer period of time we recommend to turn the machine off at the main switch (1).

7.6 Bottle-/canister change

Release the quick couplings from the machine by pressing the locking button and pulling upwards.
Screw the quick couplings to the new bottle/canister as described above.

Important: Never swap the connection sides of the bottles/canisters (see chapter 8.1).


8. Deviations from the regular processes

8.1 Mixing head only delivers one component

Cause: As a bottle or a canister was completely emptied air has been sucked in by the pump. This might also happen when the minimum filling height falls below approx. 2 cm.

Correction: Replace the emptied silicone container and disassemble the

mixing head (see sequence of pictures) from the machine. Then put an appropriate container underneath the conveying hoses. Press the activation button until silicone comes out of both hoses equally. Reassemble the mixing head.

Dosper evo

8.2 Proportion of the components clearly deviates from a 1:1 mixing.

Cause: The volume of one of the pump chambers has considerably increased due to abrasion – caused by quartz particles in the silicone.

Correction: Bring the filling height of the reservoirs (bottles/canisters) onto the same level and extract about 250 ml (corresponds to about 2–3 duplications). If the extraction volume after this try noticeably deviates in one of the components, the output can be adjusted to the required 1:1 dosage with the help of the setting screws (9/10).

Important notice: Exclusively change the output of the component, that shows a higher volume output (lower filling level).

Release the cover of the setting screw (9/10) that needs adjustment and carefully insert a slotted screwdriver into the opening until it noticeably snaps into the slot of the setting screw.

The volume output of this component can be reduced through anti-clockwise rotation.

Change the volume output only gradually. A ½ turn will already effect a volume reduction of about 4 %. Repeat this process, if necessary, after you have checked the mixing ratio.

After adjustment re-lock the setting screw with its cover.

8.3 The pump is in operation, the LEDs are flashing, no silicone is being delivered.

Cause: The silicone in the mixing tip is cured.

Correction: Replace the cured mixing tip by a new one.

8.4 The machine shows no function

Causes:

- a) The machine is not connected to the power supply. Check the connection.
- b) The main switch (1) is turned off. Turn the machine on.
- c) The sealing plugs respectively the connector plugs for the reservoirs (4a/4b) do not have contact. Pull the connectors from the machine, clean them, if necessary, and firmly insert them into their sockets.

d) The fuse is blown. Take the fuse out (see Chapter 11) and exchange, if necessary.

In case the Dosper evo still does not show any function, disconnect from the power supply and contact the responsible service technician or get in touch directly with the supplier Drewe.

9. Optional accessories

Material reservoirs: D-5022

Connection hose

5.1 l canister: D-49366

Connection hose

21.2 l canister: D-49367

Storage rack for

5.1 l canister: D-49364

Storage rack for

21.2 l canister: D-49365

Dublisil® flask:

D-3805

Mixing tips:

40 pieces

D-4970-1

100 pieces

D-49701-1

10. Care and maintenance

Generally disconnect the unit from mains supply before executing care and maintenance.

Ideally cleaning is to be effected dryly with a soft piece of cloth, if necessary

with a slightly humid sponge and a mild cleaning liquid. Water and cleaning liquids must not get into the unit.

11. Replacement of the fuse

Disconnect the power cord from the current supply.


The fuse holder is integrated in the mains connection on the back of the machine.

For the control or the change of the fuse press the upper and the lower latch of the fuse holder simultaneously and take the holder out.

Pull the defective fuse out of the clamp holder by pulling gently. Replace with a new, identical fuse.

Important notice: Only use fuses with protection power (T 1, 6A) as stated on the back of the machine.

After the change insert the fuse holder again by pressing gently until the fixation locks noticeably.

12. General advice

Repairs and maintenance shall only be executed by authorized and qualified personnel. In any case the electricity supply of the machine must be disconnected. The company Dreve is only responsible for the safety, reliability and performance of the device, if

- the machine is only used according to its purpose
- the machine is only operated by trained personnel

- the machine is used in compliance with this instructions manual
- enhancements, re-adjustments, changes or repairs are executed by persons authorized by the manufacturer Dreve

As we are constantly working on the advancement of our products, we reserve the right for technical modifications.

13. Guarantee

The right of the purchaser to bring actions arising out of defects shall in all cases be in lapse 2 years after the date of delivery. The guarantee obligation is only valid for material or production defects. It does not refer to natural wear and tear nor to damage occurring after the transfer of risk due to improper or unsuitable handling, excessive stressing, unsuitable operating material, and chemical, electrochemical or electrical influences of a nature not provided for according to the contract. On improper alterations or repairs or maintenance carried out by the purchaser or third parties, the guarantee obligations of the supplier and his liability for consequences caused thereby are rendered void.

Liability

The right of the purchaser to bring actions arising out of defects shall in all cases be in lapse one year after the date of delivery. The warranty obligation does not refer to natural wear and tear or to damage occurring after the transfer of risk due to improper or unsuitable handling, excessive stressing, unsuitable operating material, and chemical, electrochemical or electrical influences of a nature not provided for according to the contract. On improper alterations, repairs or maintenance or breach of seals, carried out by the purchaser or third parties, our warranty obligations and liability for consequences caused thereby are rendered void.


Dosper evo

14. Spare parts list

No.	Denomination	Art.No.		
2	Fuse	50048		
3	Combi element / mains connection	58514 58515		
4	Connector socket for auto stop (only if material reservoirs are used)	58516		
5a	Adapter for material container	58508		
5b	Adapter for material container	58509		
6	Mixing head	D-6505045		
7	Activation button	58513		
8	Set of LEDs	D-65058LED		
10	Setting screw / restrictor valve	58501		
11	Housing red	6558550R		
11a	Housing silver	D-6558550S		
11b	Housing special colour + printing (indicate Ral-colour)	D-65058+RAL		
12	Bottom plate	58555		
13	Appliance base foot (6 pcs.)	D-6551079		
14	Motor & pump (indicate voltage)	50755		
15	Rubber buffer (4 pcs)	D-6558560		
16	Gasket ring for housing	D-6558561		
17	Reducing nipple G3/8 a - G1/4i	58505		
18	Elbow connector 3/8 for hose 12/10	50264		
19	Polyamide hose 12x1	D-65028		
20	Straight screw joint G3/8	58502		
21	T-piece RP 3/8 I/I/I	58500		
22	Double nipple 1/4 a - 3/8 a	50197		
23	Angular swivel joint F-WSI-1/4-MSv	58504		
24	Hollow bolt single 1/4	58506		
25	Straight screw joint-V.120614 6/4-1/4	50154		
26	Tecalan tube TR 6/4 mixing head	D-6551068M		
27	Tecalan tube TR 6/4 bypass	D-6551068B		
28	Angular screw joint G 3/8	58503		
29	Condensator 5µf 220-230V 50-60Hz EC	60808		
29a	Condensator 20µf 115V 60Hz US	50796		
29b	Condensator 25µf 100V 50/60Hz Japan	50797		
30	Electronics (indicate voltage)	58510		
31	Connecting terminal plate	51067		
32	Auto cut-off plug	D-50502		
Without drawing No.				
	Denomination		Art.No.	
	Power cord Europe		51021	
	Power cord UK		51284	
	Power cord USA/Japan/Korea		51283	
	PVC gasket 1/4		50746	
	FAG balls for mixing head		50271	
	O-ring 18x2 for bottle adapter		50769	
	Pump chamber complete	D-65201		
	Set of gear wheels including lock ring	D-65202		
	Lid for pump flange complete	D-65203		
	Set of springs for mixing head	D-65207		
	Anti-interference capacitor	50799		

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