**Bedienungsanleitung · Working instructions
Mode d'emploi · Instrucciones de trabajo
Istruzioni d'uso**

Drufomat scan

**Universelles Druck-Tiefziehgerät
Universal pressure thermo-forming unit
Appareil de thermoformage universel
Termoadaptador a presión universal
Apparecchio universale per termostampaggio**

deutsch

english

français

español

italiano

Dreve

Inhalt (DE)

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Contents (EN)

1.	Description of unit	30
2.	Technical data	30
3.	Declaration of conformity	31
4.	Safety advice	32
5.	Initiation – scope of delivery	33
6.	Functional components	34
7.	Display	35
8.	Operation	36
8.1	Scan modus	36
8.2	Pressure adjustment	37
8.3	Start of programme	37
8.4	Thermoforming	38
8.5	Opening of the pressure chamber	39
9.	Deviation to the rule	40
9.1	Modification of given times	40
9.2	Creation and storage of manual time programmes	41
9.3	Calling of individual time programmes	42
9.4	Plasticizing of blanks with unknown heating times	42
9.5	Abort/termination of programme	42
9.6	Exceeding of entered heating/cooling times	43
9.7	Programming of a wrong Barcode	43
10.	Optional additional functions	44
10.1	Pressure polymerisation	44
10.2	Flask press	45
10.3	Fabrication of cast form parts	46
11.	Care and maintenance	46
12.	Exchange of fuses	47
13.	Accessories	48
14.	General advice	49
15.	Guarantee	49
16.	Trouble shooting	50
17.	Spare part list	52
18.	Explosion sketch	53
19.	Operating diagram	132
20.	Pneumatic plan	133
21.	Barcodes	134

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1. Description of unit

The Drufomat scan is a universally usable pressure thermoforming unit for daily use in the dental laboratory. An integrated scanner allows a direct taking over of product names, blank thicknesses as well as individual heating and cooling times by means of a barcode that can be found on each blank package.

A variably adjustable forming pressure of 29–87 psi (2–6 bar) in connection with a vertically forming blank guarantees for very precise results.

In combination with optionally available inserts the Drufomat scan can also be used as flask press (P insert) and as pressure pot (D insert).

2. Technical data

Dimensions (H x W x D):	460 x 290 x 290 mm
Weight:	16 kg
Supply Voltage:	100–240 V / 50–60 Hz
Max. power consumption:	115V–270 W, 230 V–315 W
Fuse of unit:	T 3, 15 A
Operating pressure:	29–87 psi (2–6 bar)
Min. operating pressure:	29 psi (2 bar)
Max. pressure in pipe:	145 psi (10 bar)

3. Declaration of Conformity

According to Machinery Directive (2006/42/EC)

Herewith we declare that the machine described below is in its conception and design and in the shape delivered by us in accordance with the fundamental requirement for safety and health as prescribed in the applicable EC-guidelines. In the case of any change or modification of the machine not authorised by us this declaration becomes invalid.

Designation of machine:	Drufomat scan D3300 / D3300X1 / D3300X2
Type of machine:	Pressure thermoforming unit
Applicable EC-guidelines:	
2006 / 42 / EC	Machinery Directive
2014 / 35 / EU	Low Voltage Directive
2014 / 30 / EU	EMC Directive

Applied harmonised standards:

2006 / 42 / EC	Machinery Directive
EN ISO 12100:2010 + Corrigendum to DIN EN ISO 12100:2011-03	Safety of machinery – General principles for design – Risk assessment and risk reduction (ISO 12100:2010); German version EN ISO 12100:2010 + Corrigendum to DIN EN ISO 12100:2011-03.
EN ISO 13732-1:2008	Ergonomics of the thermal environment – Methods for the assessment of human responses to contact with surfaces – Part 1: Hot surfaces (ISO 13732-1:2006); German version EN ISO 13732-1:2008-12.
EN 60204-1:2006/A1:2009 + Corrigendum to EN 60204-1:2006	Safety of machinery – Electrical equipment of machines – Part 1: General requirements (IEC 60204-1:2005/A1 2008); German version EN 60204-1:2006/A1:2009 + Corrigendum to EN 60204-1:2006; German version CENELEC-Cor.: 2010 to EN 60204-1:2006.
2014 / 35 / EU	Low Voltage Directive
EN 61010-1:2010	Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1: General requirements (IEC 61010-1:2010 + Cor.: 2011); German version EN 61010-1:2010.
2014 / 30 / EU	EMC Directive
EN 61000-3-2:2014	Electromagnetic compatibility (EMC) – Part 3-2: Limits – Limits for harmonic current emissions (equipment input current ≤ 16 A per phase) (IEC 61000-3-2:2014); German version EN 61000-3-2:2014.
EN 61000-3-3:2013	Electromagnetic compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current ≤ 16 A per phase and not subject to conditional connection (IEC 61000-3-3:2013; German version EN 61000-3-3:2013.
EN 61326-1:2013	Electrical equipment for measurement, control and laboratory use – EMC requirements; Part 1: General requirements (IEC 61326-1:2012); German version EN 61326-1:2013 Emission according to living area, business and industrial undertakings as well as small enterprises interference stability according to industrial area.

4. Safety advice

Attention! Read the following notes before connecting and putting the unit into service. The operating service and the function of the unit can only be guaranteed for, if the general safety and accident prevention instructions of the legislation as well as the safety instructions of the working instructions are observed.

1. This unit must only be used according to the present working instructions. We cannot be held liable for damages which have been caused by improper treatment resp. faulty manners of operation.
2. Put the unit onto a stable and even basis which is suitable for the operation weight of approx. 20 kg of the unit.
3. In order to avoid the intrusion of water into the unit (e. g. splash water) the unit is to be placed in a dry surrounding.
4. Do not store any flammable substances in the close surrounding of the unit.
5. The indicated voltage on the type plate has to meet the voltage of the power source.
6. The unit must only be connected with a socket with protective conductor. Never touch the plug with wet hands.
7. Authorised operators: The operating company of the machine has to hold the working instructions at the operators' disposal and has to make sure that authorised operators have read and understood the manual. Only in that case the operators may operate the unit.
8. The unit has to be examined before operation in view of proper condition and operating safety. If the unit is not in proper condition, it must not be used and has to be marked accordingly.
9. The pressure of the compressed air line must not exceed 145 psi.
10. Do not intrude any objects into the unit.
11. Do not lock-up the operation elements of the unit.
12. Do not block the downwards-moving piston.
13. Do not touch the downwards-moving piston.
14. The swivelling housing of the radiator reaches high temperatures and must therefore not be touched – not even for short moments.
15. If the unit is not used, please disconnect the mains.

16. Plates and labels have to be maintained in a well legible condition and must not be removed.

17. Before cleaning and maintenance of the unit or the exchange of parts the plug has to be disconnected.

18. The opening of the unit and repairs are only to be carried out by accordingly trained experts.

19. Only those appliances and spare parts may be used which are released by

5. Initiation

Before initiation make sure that the voltage indicated on the type plate meets the voltage of the power source.

Connect the unit to an air compressor. Push the pressure hose (included in the delivery) into the compressed air connection (9) until it is locked in place. The hose does not need any further fastening. To release the pressure hose push the blue ring upwards and draw the hose off by pulling downwards.

The pressure of the compressed air line must not exceed 145 psi.

Put the power cord into the female connector (33) and create a connection with the Euro/USA plug.

the manufacturer. We are not responsible for damages caused by the use of foreign parts.

20. Arbitrary rebuilding or modifications are not permitted due to safety reasons.


21. The required operation and maintenance conditions mentioned in these working instructions are to be complied with imperatively.

Important: Before the first initiation please read the working instructions thoroughly. In case of comprehension problems contact your responsible dealer or Dreve Dentamid GmbH, if necessary.


Scope of delivery

- 1 x Drufomat scan
- 1 x drawer for granules
- 1 x F insert incl. insert bowl and nut plate
- 1 x tension ring with exhaust valve
- 1 x base range of blanks
- 1 x bowl for granules
- 1 x power cord
- 1 x air connection hose incl. hose clamps and fast coupling
- 1 x working instruction


6. Functional components

No.	Description
1	scanner
2	display
3	keypad
4	heating element
5	handle heating
6	F insert
7	drawer for granules
8	push button
9	compressed air supply
10	handle plate reception
11	plate reception
12	exhaust valve
13	tension ring
33	mains supply/main breaker


7. Display

Letter	Description
a	indication of heating time
b	indication of product name and material strength
c	indication of cooling time
d	change of menu items (backwards)
e	change of time indications and pressure (upwards)
f	confirmation key/start button heating
g	indication of pressure inside piston
h	control indication for individual time modifications
i	change of time indications and pressure (downwards)
j	change of menu items (forwards)

8. Operation

Turn the unit on by means of the main switch (33) at the back of the unit. By this, the scanner (1) at the bottom and the display (2) at the front of the unit will be activated. An acoustic signal shows that the unit is ready for operation; the display will show the name of the machine Drufomat scan.

8.1 Scan modus

Every outer package of Dreve thermoforming blanks is labelled with a barcode label.


This barcode contains the encrypted product name, material strength as well as the individual heating and cooling times.

In order to scan the material data with the Drufomat scan, the barcode is to be placed centrally and evenly under the scan field (1).

There will be an acoustic signal for the confirmation of the data capture and the display will show all necessary data.

Info

The light emanating from the scanner is harmless. Even in case of direct eye contact it has no negative impact!


After 15 min. without use the scanner turns off automatically. By pressing any key it can be re-activated again.

During an active thermoforming process the scanner is deactivated. After termination of the process it will be re-activated automatically.

8.2 Pressure adjustment

The pressure of the Drufomat scan is controlled digitally and can individually be adapted to the thermoforming blank in question between 29 and 87 psi infinitely variable.

It is possible to control and to modify the adjusted pressure and the actual pressure during the thermoforming process via the bar chart (g) of the display.

8.3 Start of programme

In order to start the selected programme, put a blank that corresponds to the entered barcode resp. programming onto the plate reception (11) and fix it by means of the tension ring (13). The model can either be put onto the nub plate of the F insert (6) or fixed by taking out the nub plate and embedding it into granules directly in the F insert.

Info

For an optimum adaptation of the blank towards the model always place the models in a way that the anterior teeth are directed towards the centre of the blank.

When using hard blanks it is recommended to embed the models into granules. Models that are placed onto the nub plate should not exceed a max. height of 23 mm.

The arrow to the left of the bar chart shows the currently adjusted pressure. The arrow can be moved along the bar upwards and downwards with the +/- keys (e/i). The pressure can also be modified during an ongoing thermoforming process!

Info

The minimum operating pressure is 29 psi (2 bars). If a thermoforming process is interrupted before having reached the 29 psi (2 bar) mark (indication of bar chart (g) in the display), the piston will be ventilated and will return to its original position.


Info

For safety reasons the quartz radiator is equipped with a final switch and only works in swivelled in position after having started a programme.

If the quartz radiator is swivelled outwards after having started a programme, the heating process will be interrupted automatically. The indication in the display (a) will stop and show the remaining time. By swivelling the quartz radiator back in the programme will be re-activated.

Attention! The swivelling housing of the heating element reaches high temperatures.

Therefore it must not be touched – even not for short moments. After the adjusted time has elapsed the indicator (a) in the display will stop and count the duration after termination of the heating time. This process may be stopped by abort (long pressing of button 8) or by successive thermoforming.

Info

For safety reasons the process of an unwanted downwards movement of the piston or an interlock of the escapement can be stopped by the electronics:

- a) The process of moving the piston downwards is carried out too slowly (max. 3 seconds between action right hand and action left hand) or
- b) The two-hand control is aborted before reaching the minimum pressure.

In both cases there will be 5 short acoustic warnings, the piston will be ventilated and remove to its original position. The display

will show the last stored values. In order to start the thermoforming process newly the process can be repeated immediately, as described in this chapter.

The adjusted pressure will be verified by the electronics and regulated, if necessary. The display shows the adjusted cooling time (c) which then counts down to 00:00.

8.4 Thermoforming

After the heating time has elapsed there will be an acoustic signal; the heating time in the display will be 00:00. Heating and LED illumination will be switched off automatically. Pull the heating (4) with the right hand at the knob rod (5) as far as it will go to the front and hold it tight.

Then push the push button (8) with the left hand at the left side of the housing. The piston will move downwards and close the pressure chamber; an acoustic signal will be heard.

Info


The closing pressure of the downwards moving piston is always higher than the forming pressure!

The downwards moving piston is equipped with a dual electronic safety feature and can only be activated, if the unit is operated with two hands at the same time (two-hand control). Always make sure that no objects are in the contact area between piston and F insert while the piston is closing!

Hold both buttons until the pressure bar in the display (g) shows the adjusted minimum of 29 psi (2 bars). Then release the push button (8) and finally the knob rod of the heating (5). After having reached the 29 psi (2 bar) mark the pressure will automatically rise up to the pre-adjustment.

8.5 Opening of the pressure chamber

In order to attain a complete cool down of the thermoforming blank the valve at the tension ring (12) is to be opened slightly after 2/3 of the adjusted cooling time. Through this the warm air inside the piston will be released and replaced by cool air due to the pressure regulation of the electronics. After the cooling time has elapsed there will be an acoustic signal (indication in the display (c) will blink).

Directly before releasing the pressure inside the piston the ventilation valve at the F insert (6) should be opened. Push the push button (8) for approx. two seconds; the pressure inside the piston will be released. After complete ventilation the piston will remove to its original position. Swivel out the plate reception (11) at the corresponding handle (10) and remove tension ring and thermoformed result.

The display (2) will show the last used programme; the scanning device (1) is ready for a new thermoforming process.

9. Deviation to the rule


9.1 Modification of given times

Scanned times may be prolonged or shortened individually, if required.


Info

The performance of the heating might decrease slightly after 1000 operating hours. By prolonging the heating times by 5–10 seconds optimum thermoforming results can be attained again. The times given by the barcodes are average times and designed for general use. In case of extreme undercut areas or very high models a prolongation of the heating time can optimise the results.

After having scanned the pre-adjusted material parameters (see paragraph 8.1) the heating time (1x) or the cooling time (2x) can be selected by pressing the menu key (j). There will be an arrow in the display in front of the chosen time. The heating time can be modified positively and negatively in steps of 5 seconds. For this press the menu keys (e/i) according to the key description. The max. deviation from the scanned time is 35 seconds.


With the same working steps the cooling time can be modified in steps of 10 seconds by max. 70 seconds.

For identification of the modified parameters the display will show a + resp. a - sign. The degree of deviation will be shown by a scale at the bottom of the display. Modified values have to be confirmed by pressing the OK key (f).

All further working steps correspond to 8.2.

9.2 Creation and storage of manual time programmes

Up to three individual time programmes may be created and stored for thermoforming foils that are not equipped with a corresponding barcode.

Turn on the unit – the initial display will be visible. When pushing the menu key (j) one time the display will show ----- at position b.

As described in point 9.1, push menu key (j) and select the desired heating time by the menu keys +/- (e/i). Following push menu key (j) again and enter cooling time.

For single use confirm the two entered values by pushing OK (f). Then go on as described in point 8.2.


If the created time programme is to be stored for repeated use, 3 memory capacities (F1–F3) are available. Do not confirm the created time programme by pushing OK, but press menu key (j) again. The display will show F1 in position b, the entered times will be taken over. The memory capacities F2 and F3 can be selected by pushing the menu key (j) 2 resp. 3 times. Store the times in the desired memory capacity by pushing the OK key (f). Then go on as described in point 8.2.

9.3 Calling of individual time programmes

Turn on the unit – the initial display will be visible. By pushing the menu key (j) several times you may select one after another

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By pressing the menu keys (d/j) you may navigate arbitrarily between the a. m. programme points. Confirm the selected programme point with the OK key (f).

Then go on as described in point 8.2.

9.4 Plasticizing of blanks with unknown heating times

Create an individual time programme with high heating and cooling times (see point 9.2). Start the programme (see point 8.2). You may test the plasticization degree of the blank by means of an edgeless metal instrument (e. g. wax knife).

Info

The complete deformability of the blank can only be tested at the edges, as that is where warmth

will leak through the metal contacts and the necessary flexibility will be attained at last!

After the blank has attained the necessary flexibility initiate the thermoforming process (see point 8.4). The remaining heating time will be erased from the display (a) and the adjusted cooling time will be counted down (c). The cooling time can be terminated at any time by pushing the push button (8).

9.5 Abort/termination of programme

It is possible to abort / terminate all actions of the Drufomat scan by pushing the push button (8) at the left side

of the unit (keep pushed for aprox. 2 seconds!).

9.6 Exceeding of entered heating/cooling times


After elapse of times adjusted regularly (display position a/c) there will be an acoustic signal. The exceeded time

will be blinking in the display as control mechanism.

9.7 Entering a wrong barcode

If a barcode is scanned that does not correspond to the programming of the Drufomat scan, four acoustic warning signals will be heard and the display will show ERROR at position b.

By scanning the correct barcode you can proceed further, as described in point 8.0.

10 Optional additional functions

10.1 Pressure polymerisation

In combination with the optionally available D insert the Drufomat scan can be used for pressure polymerisation of resin repairs or as compressor of plasters, embedding masses and silicones.

Turn on the unit, swivel out plate reception (11) and heating element completely (4) by means of the operating knobs (10/5) and take the F insert (6) out of the bottom plate. The material to be compressed (e. g. cast flask, duplication etc.) is to be placed inside the D insert.

For the polymerisation of resins the D insert is to be filled before with hand warm water (corresponding to the indications of the resin manufacturer).

Fix the D insert in the bottom plate. By moving the piston downwards the D insert will be closed and will build up the pre-adjusted pressure inside.


For this purpose go on as described in point 8.4.

10.2 Flask press

In combination with the optionally available P insert the Drufomat scan can also be used as flask press. One or two flasks may be pressed at the same time.

Turn on the unit, swivel out plate reception (11) and heating element completely (4) by means of the operating knobs (10/5) and take the F insert (6) out of the bottom plate. Insert the P insert into the pressure chamber from underneath.

For one flask: Put back the F insert (6) into the bottom plate.

For two flasks: Remove nub plate and insert bowl from F insert and insert them into the bottom plate.

Place flasks onto the nub plate. By moving the piston with the P insert downwards the necessary pressure will be exerted upon the flasks.

For this purpose go on as described in point 8.4.

Info

In order to attain a homogeneous distribution of the material to be compressed in the flask we recommend an initial pressure of 29 psi (2 bars) which can be increased step by step up to 87 psi (6 bars) during the pressing procedure (see point 8.2).

10.3 Fabrication of cast form parts

With the K-7 insert it is possible to thermoform up to 7 resin bases for the crown and bridge technique at the same time.

Fix the single dies by means of a flexible blocking out material (e. g. Fillin) in the holes of the K-7 insert.

Remove the nub plate from the F insert (6) and replace it by the K-7 insert.

All further working steps as described from point 8.1 on.

11. Care and maintenance

Generally disconnect the unit from mains supply before executing care and maintenance.

Ideally cleaning is to be effected dryly with a soft piece of cloth, if necessary with a slightly humid sponge and a mild cleaning liquid. Water and cleaning liquids must not get into the unit.

The O-rings should be greased regularly with silicone grease.

Removal and cleaning of the filter sieve

Disconnect the power cable from the socket and shut down the compressed air supply. Push the blue ring of the compressed air connection upwards and draw the hose off the unit by pulling it downwards. Now unscrew the connection using a suitable tool. Then carefully remove the filter sieve with a slotted screwdriver by turning anti-clockwise. Clean with stream of compressed air and reassemble in reverse order.


12. Exchange of fuses

Separate power cord from the mains. The fuse housing is to be found on the left back side of the unit.

In order to exchange the fuse press the locking cap slightly inside and turn it to the left in accordance with the arrow.

Take out fuse carrier. Pull the defect fuse carefully out of the clamp and replace it by a new one.

Attention: Only use the fuses indicated on the back of the unit.

Insert fuse carrier again and push it slightly to the right.

13. Accessories

F forming table insert REF D6505590

- F insert with cone plate
- insert bowl
- nub plate
- drain bowl

(In order to empty the F insert filled with granules the insert is to be placed onto the drain bowl; the granules will automatically drain through the lower opening of the F insert.)

Not included in delivery:

D insert REF D3227
Pressure polymerisation insert

P insert REF D3211
Pressure insert for pressing flasks

K-7 insert REF D3228
Die reception

14. General advice

Repairs and overhauling on the Drufomat scan are only to be executed by authorised experts. Disconnect the unit from mains supply at any rate. Dreve company is only responsible for the safety, reliability and performance of the unit, if:

- the unit is exclusively used for the corresponding application.
- operation is only carried out by accordingly trained personnel.

15. Guarantee

The guarantee time amounts to two years after the date of delivery. Guarantee is only accepted in case of material or fabrication defects. The guarantee obligation does not refer to natural wear and tear nor to damage occurring after the transfer of risk due to improper or unsuitable handling, excessive stressing, unsuitable operating material, and chemical, electrochemical or electrical influences of a nature not provided for according to the contract. On improper alterations, repairs or maintenance or breach of seals carried out by the purchaser or third parties our warranty obligations and liability for consequences caused thereby are rendered void.

- the unit is used in accordance with these working instructions.
- extensions, new adjustments, changes or repairs are executed by persons authorised by Dreve company.

As we permanently develop our products further, we reserve the right for technical changes.

Liability

The right of the purchaser to bring actions arising out of defects shall in all cases be in lapse one year after the date of delivery. The warranty obligation does not refer to natural wear and tear nor to damage occurring after the transfer of risk due to improper or unsuitable handling, excessive stressing, unsuitable operating material, and chemical, electrochemical or electrical influences of a nature not provided for according to the contract. On improper alterations, repairs or maintenance or breach of seals carried out by the purchaser or third parties our warranty obligations and liability for consequences caused thereby are rendered void.

16. Trouble shooting



Problem	Cause	Aid	Problem	Cause	Aid
Unit does not work	Power line is not correct	Control power supply	Electronic can be set, but radiator does not heat	Heating element is not swivelled-in completely	Swivel in the heating element
	Main breaker is off	Switch on main breaker		Micro-switch defect	Micro-switch has to be exchanged by an expert
	Miniature fuse is defect	Check miniature fuse and exchange it, if necessary	Heating time elapses, but radiator does not heat	Quartz radiator defect	Exchange quartz radiator
Main breaker on, electronic shows no function	Electronic is out of order	Electronic has to be exchanged by an expert		Temperature monitor was activated	Exchange temperature monitor
Piston cannot move downwards	Two-hand control has not been operated	Operate two-hand control correctly, see working instructions	Heating time elapsed, blank is not plasticized sufficiently	Lifetime of quartz radiator exceeded	Exchange quartz radiator
	Unit is not connected to compressed air	Control compressed air connection – connect to compressed air, if necessary		Heating time was too short	Prolong heating time by 5–10 sec.
	Final switch is out of order	Final switch has to be exchanged by an expert	Scanner not active	Scanner in sleeping position	Press any key
	Magnet valve is out of order	Magnet valve has to be exchanged by an expert		Scanner or electronic defect	Electronic has to be checked by an expert or exchanged, if necessary
	Push button is defect	Push button has to be exchanged by an expert	Piston moves downwards, unit aborts process	Push button was pushed too long	Operate two-hand control correctly, see working instructions
Piston has moved downwards, unit does not take over self-maintenance	Minimum pressure of 29 psi has not been reached	Hold escapement for longer, see working instructions	Unit does not reach the pre-adjusted pressure	Line pressure too low	Increase line pressure
	Electronic is out of order	Electronic has to be checked by an expert	Unit makes fizzy sounds during thermoforming	Ventilation valves not closed	Close ventilation valves
Piston does not move upwards, unit cannot be opened	Wrong operation of push button	Hold push button for aprox. 2 sec.		Tension ring has been placed wrongly into plate reception	Put tension ring in correct position
	Push button is defect	Push button has to be exchanged by an expert		Connection hose leaky between valve and cylinder	Unit has to be controlled by an expert, connection hose has to be exchanged, if necessary
	Two-hand control or magnet valve are defect	Unit has to be controlled by an expert		Quad rings are dry	Grease rings with silicone grease, have them exchanged by an expert, if necessary
Piston does not come down or only very slowly	Filter sieve on connection for compressed air is dirty	Clean filter sieve			

17. Spare part list

No.	Description	REF	Description	REF
5	Handle heating	D6505519	Power cord Europe	51021
6	F insert	D6505590	Power cord UK	51284
7	Drawer for granules	5565	Power cord USA/JP	51283
8	Push button	D6505528	Electronic	
9	Compressed air connection, complete	D6505536	+ display + keyboard	D6505540X1
10	Handle plate reception	D6505519	Piston	51316
11	Plate reception	D6505512	Control card	D6505572X1
12	Ventilation valve	50761	LED lamps, red (4x)	5543X1
13	Tension ring	D6505518	Fuse F3, 15 AT	51321
14	Cover glas scanner	5594	Piston guidance ring	51315
15	Scanner	D6505541X1	Valve set piston	D6505592
16	Magnetic valve unit	D6505525	Quad-ring, large	51098
17	Housing red	D655510R	Quad-ring small	51097
	silver	D655510S	O-ring 112 x 3 mm	51077
18	Control elements heating	D6505521	O-ring 164 x 3 mm	51099
19	Microswitch	5521	Cable harness	5520
20	Heating housing	D6505513	Limit heating element	D6505561
21	Quartz-infrared radiator	D6505542	Screw joint for compressed	
22	Temperature fuse	5542	Air connection	50340
23	Pillar heating element	D6505515	Screw-in filter sieve	57030
24	Power supply	D6505571		
25	Magnetic valve	D6505526		
26	Cylinder	51314		
27	Display foil, adhesive	5582TEK		
28	Limit plate reception	D6555344		
29	Holder plate reception	D6505514		
30	Pillar plate reception	D6505591		
31	Bottom plate	D6505511		
32	Fuse carrier	5530		
33	Mains supply/ main breaker	5531		

When ordering spare parts, please always mention the serial number of the unit.

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