Bedienungsanleitung · Working instructions
Mode d'emploi · Instrucciones de trabajo
Istruzioni d'uso · Gebruiksaanwijzing

Drufosmart scan

Universelles Druckformgerät 4 bar

Universal pressure thermoforming unit 58 psi

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Contents (EN)

1.	Description of unit	26
2.	Technical data	26
3.	Declaration of conformity	27
4.	Safety advice	28
5.	Initiation – scope of delivery	29
6.	Functional components	30
7.	Display	31
8.	Operation	32
8.1	Scan modus	32
8.2	Pressure adjustment	33
8.3	Start of programme	33
8.4	Thermoforming	34
8.5	Opening of the pressure chamber	35
9.	Deviation to the rule	36
9.1	Modification of given times	36
9.2	Creation and storage of manual time programmes	37
9.3	Calling of individual time programmes	38
9.4	Plasticizing of blanks with unknown heating times	38
9.5	Abort / termination of programme	38
9.6	Exceeding of entered heating/cooling times	39
9.7	Entering a wrong barcode	39
10.	Additional optional functions	40
10.1	Fabrication of cast form parts	40
11.	Care and maintenance	40
12.	Exchange of fuses	41
13.	General advice	41
14.	Trouble shooting	42
15.	Spare part list	44
16.	Explosion sketch	125
17.	Operating diagram	126
18.	Pneumatic plan	127
19.	Barcodes	128

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1. Description of unit

The Drufosmart scan is a universally usable pressure thermoforming unit for daily use in the dental laboratory. An integrated scanner allows a direct taking over of product names, blank thicknesses as well as individual heating and cooling times by means of a barcode that can be found on each blank package.

A defined adjusted forming pressure of 4 bar in connection with a vertically forming blank guarantee for very precise results.

2. Technical data

Dimensions (H x W x D):	420 x 320 x 290 mm
Weight:	15 kg
Supply voltage:	100–240 V / 50–60 Hz
Max. power consumption:	115 V–270 W, 230 V–315 W
Fuse of unit:	T 3.15 A
Operating pressure:	4 bar
Min. operating pressure:	2 bar
Max. pressure in pipe:	10 bar

3. Declaration of Conformity

According to Machinery Directive (2006 / 42 / EC) attachment II 1.A

Herewith we declare that the machine described below is in its conception and design and in the shape delivered by us in accordance with the fundamental requirement for safety and health as prescribed in the applicable EC-guidelines. In the case of any change or modification of the machine not authorised by us this declaration becomes invalid.

Denomination of the machine: Drufomat scan
D3400 / D3400A / D34001 / D34001A

Type of machine: Pressure thermoforming unit

The following additional ECdirectives were applied:

2006 / 35 / EC Low Voltage Directive
2014 / 35 / EU EMC Directive

The protection objectives of EC Directive 2006 / 42 / EC are complied with.

Applied harmonised standards:

2006 / 42 / EC

EN ISO 12100:2010 + Corrigendum to DIN EN ISO 12100:2011-03 Safety of machinery – General principles for design – Risk assesment and risk reduction (ISO 12100:2010); German version EN ISO 12100:2010 + Corrigendum to DIN EN ISO 12100:2011-03.

EN ISO 13732-1:2008 Ergonomics of the thermal enviroment – Methods for the assessment of human responses to contact with surfaces – Part 1: Hot surfaces (ISO 13732-1:2006); German version EN ISO 13732-1:2008-12.

EN 60204-1:2006 / A1: 2009 + Corrigendum to 60204-1:2006 Safety of machinery – Electrical equipment of machines – Part 1: General requirements (IEC 60204-1:2005 / A1 2008); German version EN 60204-1:2006 / A1:2009 + Corrigendum to EN 60204-1:2006; German version CENELEC-Cor.: 2010 to EN 60204-1:2006.

2014 / 35 / EU

EN 61010-1:2010 Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1: General requirements (IEC 61010-1:2010 + Cor.: 2011); German version EN 61010-1:2010.

2014 / 30 / EU

EN 61000-3-2:2014 Electromagnetic compatibility (EMC) – Part 3-2: Limits – Limits for harmonic current emissions (equipment input current <= 16 A per phase) (IEC 61000-3-2:2014); German version EN 61000-3-2:2014.

EN 61000-3-3:2013 Electromagnetic compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current <= 16 A per phase and not subject to conditional connection (IEC 61000-3-3:2013; German version EN 61000-3-3:2013..

EN 61326-1:2013 Electrical equipment for measurement, control and laboratory use – EMC requirements; Part 1: General requirements (IEC 61326-1:2012); German version EN 61326-1:2013 Emission according to living area, business and industrial undertakings as well as small enterprises interference stability according to industrial area.

4. Safety advice

Attention! Read the following notes before connecting and putting the unit into service. The operating service and the function of the unit can only be guaranteed for, if the general safety and accident prevention instructions of the legislation as well as the safety instructions of the working instructions are observed.

1. This unit must only be used according to the present working instructions. We cannot be held liable for damages which have been caused by improper treatment resp. faulty manners of operation.
2. Put the unit on a stable and even basis which is suitable for the operation weight of approx. 15 kg of the unit.
3. In order to avoid the intrusion of water into the unit (e. g. splash water) the unit is to be placed into a dry surrounding.
4. Do not store any flammable substances in the close surrounding of the unit.
5. The indicated voltage on the type plate has to meet the voltage of the power source.
6. The unit must only be connected with a socket with protective conductor. Never touch the plug with wet hands.
7. Authorised operators: The operating company of the machine has to hold the working instructions at the operators' disposal and has to make sure that authorised operators have read and understood the manual. Only in that case the operators may operate the unit.
8. The unit has to be examined before operation in view of proper condition and operating safety. If the unit is not in proper condition, it must not be used and has to be marked accordingly.
9. The pressure of the compressed air line must not exceed 10 bars.
10. Do not intrude any objects into the unit.
11. Do not lock-up the operation elements of the unit.
12. Do not block the downwards-moving piston.
13. Do not touch the downwards-moving piston.
14. The swivelling housing of the radiator reaches high temperatures and must therefore not be touched – not even for short moments.
15. If the unit is not used, please disconnect the mains.

16. Plates and labels have to be maintained in a well legible condition and must not be removed.

17. Before cleaning and maintenance of the unit or the exchange of parts the plug has to be disconnected.

18. The opening of the unit and repairs are only to be carried out by accordingly trained experts.

19. Only those appliances and spare parts may be used which are released by the manufacturer. We are not responsible for damages caused by the use of foreign parts.

5. Initiation

Before initiation make sure that the voltage indicated on the type plate meets the voltage of the power source.

Connect the unit to an air compressor. Push the pressure hose (included in the delivery) into the compressed air connection (9) until it is locked in place. The hose does not need any further fastening. To release the pressure hose push the blue ring upwards and draw the hose off by pulling downwards.

The pressure of the compressed air line must not exceed 10 bars.

Put the power cord into the female connector (33) and create a connection with the Euro/USA plug.

20. Arbitrary rebuilding or modifications are not permitted due to safety reasons.

21. The required operation and maintenance conditions mentioned in these working instructions are to be complied with imperatively. When working with the Drufosmart scan the general accident prevention regulations have to be observed.

Important: Before the first initiation please read the working instructions thoroughly. In case of comprehension problems contact your responsible dealer or Dreve Dentamid GmbH, if necessary.


Scope of delivery

- 1 Drufosmart scan
- 1 drawer for granules
- 1 insert bowl and nub plate
- 1 tension ring with ventilation valve
- 1 base range of blanks
- 1 power cord
- 1 air connection hose incl. hose clamps and fast coupling
- 1 working instructions in 6 languages


6. Functional components

- | No. | Description |
|-----|---------------------------|
| 1 | Scanner |
| 2 | Display |
| 3 | Key pad |
| 4 | Heating element |
| 5 | Heating handle |
| 6 | Nub plate and insert bowl |
| 7 | Drawer for granules |
| 9 | Compressed air supply |
| 10 | Handle plate reception |
| 11 | Plate reception |
| 12 | Ventilation valve |
| 13 | Tension ring |
| 33 | Mains supply/fuse carrier |
| 34 | Activation lever |
| 43 | Main switch |


7. Display

- | No. | Description |
|-----|------------------------------------------------------|
| a | indication of heating time |
| b | indication of product name and material strength |
| c | indication of cooling time |
| d | change of menu items (backwards) |
| e | change of time indications (upwards) |
| f | confirmation key / start button heating |
| g | indication of pressure inside piston |
| h | control indication for individual time modifications |

8. Operation

Turn the unit on by means of the main switch (43) at the back of the unit. By this the scanner (1) and the display (2) at the front of the unit will be activated. An acoustic signal shows that the unit is ready for operation; the display will show the name of the machine Drufosmart scan.

8.1 Scan modus

Every outer package of Dreve thermoforming blanks is labelled with a barcode label.


This barcode contains the encrypted product name, material strength as well as the individual heating and cooling times.

In order to scan the material data with the Drufosmart scan the barcode is to be placed as centrally and evenly as possible onto the scan field (1).

There will be an acoustic signal for the confirmation of the data capture and the display will show all necessary data.

Info
The light emanating from the scanner is harmless. Even in case of direct eye contact it has no negative impact!

After 15 min. without use the scanner turns off automatically. By pressing the OK key it can be re-activated again. During an active thermoforming process the scanner is deactivated. After termination of the process it will be re-activated automatically.

8.2 Pressure adjustment

The incoming pressure from the compressor must not exceed 10 bars! If necessary, assemble a pressure reducing valve into the supply resp. decrease the total pressure of the compressor.

The working pressure of the Drufosmart scan is 4 bars (58 psi). The pressure is controlled digitally and will be controlled permanently.

The effective pressure in the pressure chamber during the thermoforming process can be controlled via the bar chart (g) of the display.

8.3 Start of programme

In order to start the selected programme, put a blank that corresponds to the entered barcode resp. programming onto the plate reception (11) and fix it by means of the tension ring (13).

The model can either be put onto the nub plate (6) or fixed by taking out the nub plate and embedding it into granules in the insert bowl (6a).

Info
For an optimum adaptation of the blank towards the model always place the models in a way that the anterior teeth are directed

towards the centre of the blank. When using hard blanks it is recommended to embed the models into granules. Models that are placed onto the nub plate should not exceed a max. height of 23 mm.

Swivel the plate reception inwards at the handle (10) just up to the mechanical interlock in the centre of the unit. Also swivel in the heating element (4) at the corresponding handle (5) just up to the mechanical interlock in the centre of the unit. Close the ventilation valve of the tension ring (12).

Press the OK key (f). There will be an acoustic signal as confirmation for the activation of the heating; the scanner will turn off and the adjusted heating time will elapse in the display.

Info

For safety reasons the quartz radiator is equipped with a final switch and only works in swivelled in position after having started a programme.

If the quartz radiator is swivelled outwards after having started a programme, the heating process will be interrupted automatically. The indication in the display (a) will stop and show the remaining time. By swivelling the quartz radiator back in the programme will be re-activated.

8.4 Thermoforming

After the heating time has elapsed there will be an acoustic signal; the heating time in the display will be 00:00. The heating will be switched off automatically. Pull the heating (4) at the handle (5) away from the centre of the unit as far as it will go with the right hand.

Then push the activation lever (34) with the left hand downwards as far as it will go. By pulling the heating handle (5) again as far as it will go the pressure chamber closes completely and pressurising of the piston will be carried out automatically. This is confirmed by an acoustic signal.


Attention! The swivelling housing of the heating element reaches high temperatures. Therefore it must not be touched – even not for short moments.

After the adjusted time has elapsed the indicator (a) in the display will stop and count the duration after termination of the heating time. This process may be stopped by successive thermo-forming.


Info

For safety reasons the pressure chamber can only be pressurised, if it has been moved downwards completely manually. In order to protect the operating personnel securely against bruises of the fingers, pressurising can only be effected, if the unit is operated with two hands at the same time (two-hand control). Always make sure that no objects are in the contact area between piston and insert bowl while the piston is closing!

There will be an automatic pressurising up to 4 bars. The pressure inside the pressure chamber will be verified permanently by the electronics and regulated, if necessary. The display shows the given cooling time (c) which then counts down to 00:00.

**8.5 Opening of the pressure chamber**


In order to attain a complete cool down of the thermoforming blank the valve at the tension ring (12) is to be opened slightly after 2/3 of the adjusted cooling time. Through this the warm air inside the piston will be released and replaced by cool air due to the pressure regulation of the electronics. After the cooling time has elapsed there will be an acoustic signal (indication in the display (c) will blink).

Pull the heating handle (5) again as far as it will go and maintain it in this position (approx. 2–3 seconds) until the pressure is relieved.

Then swivel the pressure chamber into its original position by leading back the activation lever (34). Swivel out the plate reception (11) at the corresponding handle (10) and remove tension ring and thermoformed result.

The display (2) will show the last used programme; the scanning device (1) is ready for a new thermoforming process.

9. Deviation to the rule


9.1 Modification of given times

Scanned times may be prolonged or shortened individually, if required.

Info


The performance of the heating might decrease slightly after 1000 operating hours. By prolonging the heating times by 5–10 seconds optimum thermoforming results can be attained again. The times given by the barcodes are average times and designed for general use. In case of extreme undercut areas or very high models a prolongation of the heating time can optimise the results.

After having scanned the pre-adjusted material parameters (see paragraph 8.1) the heating time (1 x) or the cooling time (2 x) can be selected by pressing the menu key (j). There will be an arrow in the display in front of the chosen time.


The heating time can be modified positively and negatively in steps of 5 seconds. For this press the menu keys (e/i) according to the key description. The max. deviation from the scanned time is 35 seconds.

With the same working steps the cooling time can be modified in steps of 10 seconds by max. 70 seconds.

For identification of the modified parameters the display will show a + resp. a - sign. The degree of deviation will be shown by a scale at the bottom of the display.

Modified values have to be confirmed by pressing the OK key (f).


All further working steps correspond to 8.2.

9.2 Creation and storage of manual time programmes

Up to three individual time programmes may be created and stored for thermoforming foils that are not equipped with a corresponding barcode.

Turn on the unit – the initial display will be visible. When pushing the menu key (j) one time the display will show ----- at position b.


As described in point 9.1, push menu key (j) and select the desired heating time by the menu keys +/- (e/i). Following push menu key (j) again and enter cooling time.

For single use confirm the two entered values by pushing OK (f). Then go on as described in point 8.2.

If the created time programme is to be stored for repeated use, 3 memory capacities (F1–F3) are available.

Do not confirm the created time programme by pushing OK, but press menu key (j) again. The display will show F1 in position b, the entered times will be taken over. The memory capacities F2 and F3 can be selected by pushing the menu key (j) 2 resp. 3 times. Store the times in the desired memory capacity by pushing the OK key (f).


Then go on as described in point 8.2.

9.3 Calling of individual time programmes

Turn on the unit – the initial display will be visible. By pushing the menu key (j) several times you may select one after another:

- a) enter modus
- b) F1
- c) F2
- d) F3

By pressing the menu keys (d/j) you may navigate arbitrarily between the a. m. programme points. Confirm the selected programme point with the OK key (f).

Then go on as described in point 8.2.

9.4 Plasticizing of blanks with unknown heating times

Create an individual time programme with high heating and cooling times (see point 9.2). Start the programme (see point 8.2). You may test the plasticization degree of the blank by means of an edgeless metal instrument (e. g. wax knife).

Info

The complete deformability of the blank can only be tested at the edges, as that is where warmth will leak through the metal con-

tacts and the necessary flexibility will be attained at last!

After the blank has attained the necessary flexibility initiate the thermoforming process (see point 8.4). The remaining heating time will be erased from the display (a) and the adjusted cooling time will be counted down (c).

The cooling time can be terminated at any time by pulling the heating handle (5).

9.5 Abort/termination of programme

During an active time programme the key pad is deactivated. It is possible to abort/terminate all actions of the

Drufosmart scan by pulling the heating handle (5) as far as it will go (for approx. 3 seconds).

9.6 Exceeding of entered heating/cooling times


After elapse of times adjusted regularly (display position a/c) there will be an acoustic signal. The exceeded

time will be blinking in the display as control mechanism.

9.7 Entering a wrong barcode

If a barcode is scanned that does not correspond to the programming of the Drufosmart scan, four acoustic warning signals will be heard and the display will show ERROR at position b.


By scanning the correct barcode you can proceed further, as described in point 8.0.

10. Optional additional functions

10.1 Fabrication of cast form parts

With the K-7 insert it is possible to thermoform up to 7 resin bases for the crown and bridge technique at the same time. Fix the single dies by means of a flexible blocking out material (e.g. Fillin) in the holes of the K-7 insert. Remove the nub plate from the insert bowl and replace it by the K-7 insert.

All further working steps as described from point 7.1 on.

11. Care and maintenance

Generally disconnect the unit from mains supply before executing care and maintenance.

Ideally cleaning is to be effected dryly with a soft piece of cloth, if necessary


with a slightly humid sponge and a mild cleaning liquid. Water and cleaning liquids must not get into the unit.

The O-rings should be greased regularly with silicone grease.

12. Exchange of fuses

Disconnect the unit from mains supply. The fuse carrier is situated at the rear of the unit, integrated in the mains supply.

In order to control or exchange the fuses push the push catch of the carrier, keep it pressed and remove the carrier from the unit.

Pull the defective fuse out of the clamp carrier applying slight traction and replace it by a new, equivalent one.



Important: Only use the fuse strength indicated on the rear.

13. General advice

Repairs on the Drufosmart scan units may only be executed by experts. Disconnect the unit from mains supply at any rate. Dreve company is only responsible for the safety, reliability and performance of the unit, if:

- the unit is exclusively used for the corresponding field of application.
- operation is only carried out by accordingly trained personnel.

After exchanging the fuse insert the fuse carrier into the housing applying slight pressure until the fixation latches noticeably.

- the unit is operated in accordance with these working instructions.
- extensions, new adjustments, changes or repairs are executed by persons authorised by the company.

As we are currently developing our products, we reserve the right for technical changes.

14. Trouble shooting

Problem	Cause	Aid
Unit does not work	Power line is not correct	Control power supply
	Main breaker is off	Switch on main breaker
	Miniature fuse is defect	Check miniature fuse and exchange it, if necessary
Main breaker on, electronics show no function	Electronics are out of order	Electronics have to be exchanged by an expert
Piston cannot move downwards	Heating is not swivelled outwards	Swivel heating outwards
Piston has moved downwards, unit does not take over self-maintenance	Electronics are out of order	Electronics have to be checked by an expert
Piston does not move upwards	Wrong operation of push button	Pull push button for 3 sec.
	Contact pressure has not been relieved	Relieve contact pressure, see working instructions
	Final switch is defect	Final switch has to be exchanged by an expert
	Magnet valve or electronics are defect	Unit has to be controlled by an expert
Electronics can be set, but radiator does not heat	Heating element is not swivelled-in completely	Swivel the heating element in correctly
	Micro-switch defect	Micro-switch has to be exchanged by an expert

Problem	Cause	Aid
Heating time elapses, but radiator does not heat	Quartz radiator defect	Exchange quartz radiator
	Temperature monitor was activated	Exchange temperature monitor
Heating time elapsed, blank is not plasticised sufficiently	Lifetime of quartz radiator exceeded	Exchange quartz radiator
	Heating time was too short	Prolong heating time by 5–10 sec.
Scanner not active	Wrong heating time was chosen	Scan correct barcode
	Scanner in sleeping position	Press OK key briefly
Unit does not reach the pre-adjusted pressure	Scanner or electronics defect	Electronics have to be checked by an expert or exchanged, if necessary
	Line pressure too low	Increase line pressure
Unit makes fizzy sounds during thermoforming	Ventilation valve not closed	Close ventilation valve
	Tension ring has been placed wrongly into plate reception	Put tension ring into correct position
	Connection hose leaky between valve and cylinder	Unit has to be controlled by an expert, connection hose has to be exchanged, if necessary

Drufosmart scan

15. Spare part list

No.	Description	REF
5	Heating handle	D6505519
6	Insert bowl with nub plate	50591/50592
7	Drawer for granules	5565
9	Compressed air connection, complete	D6505536SM
10	Handle plate reception	D6505519
11	Plate reception	D6505512
12	Ventilation valve	50761
13	Tension ring	D50038
14	Cover glass scanner	1429
15	Scanner	D6505541X1
17	Housing red	D65555RSCAN
	Housing silver	D65555SSCAN
19	Micro-switch	55790
20	Heating housing	D6505513
21	Quartz-infrared-radiator	D6505542
22	Temperature monitor	5542
23	Pillar heating element	55773
24	Power supply	D6505571
25	Magnetic valve	5202124
26	Cylinder	52020
27	Display foil, adhesive	55810
28	Limit plate reception	D6555344
29	Support plate reception	D6505514
30	Pillar plate reception	55772
31	Bottom plate	D6505511
33	Mains supply/fuse carrier	51208
34	Activation lever	D65885
35	Connector lever/axis	D65017
36	Bushings	55708
37	Final switch activation lever	5521
38	Final switch heating	55790
39	Recuperating spring	55861/55862
40	Lever axis	D6555847
43	Main switc	50841
	without drawing no.	
	Power cord Europe	51021
	Power cord England	51284
	Power cord America	51283
	Electronics incl. display, key pad, connection cable	D6505530
	Fuse F3.15 AT	51321
	O-ring 112 x 3 mm	51077
	Cable harness complete	55520
	Limit heating element	D65016
	Control Card	D6505572X1

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