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
1. Device description

The EyeVolution® is an LED-based light polymerisation device for all light-curing materials with the exception of tooth-coloured composites.

The integrated LEDs cover both the requisite UVA and blue light ranges and offer a very high penetration depth of objects to be polymerised at an extremely low thermal load.

A reflective rotary plate guarantees shadow-free polymerisation even of larger objects at short polymerisation times.

The 12-V device is particularly energy-saving and can be powered by any electrical socket.

2. Technical data

Dimensions (H x W x D):

205 x 205 x 255 mm

Weight:

1.5 kg

Mains voltage

Input:

100–240 V / 50–60 Hz, 0.7A

Output:

12 V / 1.5 A

Max. power consumption:

max. 18 W

Polymerisation range:

385–390 + 465–470 nm

Plate rotation speed:

4.5 rpm

Time programmes:

1/2/3/5/10 min

3. Declaration of Conformity

in accordance with the EMC Directive (2004/108/EC)

We hereby declare that the device outlined below complies with the basic safety and health requirements of the EU Directives in terms of its design and construction as well as the version placed on the market by us.

This declaration shall not apply in the event of modifications made to the machine without our approval.

Machine designation: EyeVolution®
D-38100

Type of unit: Light polymerisation device

Applicable EC Directives:
2004/108/EG EMC Directive

Applied harmonised standards:
2004/108/EG EMC Directive

EN 61000-3-2:2006 +
A1:2009 + A2:2009

Electromagnetic Compatibility (EMC) – Part 3-2: Limits – Limits for harmonic current emissions (equipment input current ≤ 16 A per phase) (IEC 61000-3-2:2005 + A1:2008 + A2:2009); German version EN 61000-3-2:2006 + A1:2009 + A2:2009.

EN 61000-3-3:2008

Electromagnetic Compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems for equipment with rated current ≤ 16 A per phase and not subject to conditional connection (IEC 61000-1:2008); German version EN 61000-3-3:2008.

EN 61326-1:2006

Electrical equipment for measurement, control and laboratory use – EMC requirements. Part 1: General requirements (IEC 61326-1:2005); German version EN 61326-1:2006 Emissions by residential area, commercial area and small businesses; interference immunity by industrial area.

Dreve Dentamid GmbH · Max-Planck-Straße 31 · 59423 Unna/Germany

Tel.: +49 2303 8807-0 · Fax: +49 2303 8807-55 · dentamid@dreve.de · www.dreve.com/dentamid

4. Safety instructions

NB: Please read these instructions carefully before connecting and commissioning the device. The operational safety and function of the device can only be guaranteed if both the general statutory safety and accident prevention guidelines and the safety information in the operating instructions are observed.

- 1.** The device may only be used in accordance with these operating instructions. We do not accept any liability for damage incurred as a result of non-designated use or incorrect operation.
- 2.** The device must be set up on a stable and level surface.
- 3.** In order to prevent water (e. g. spray water) ingress into the device, the device should be set up in a dry environment.
- 4.** Do not store any highly-inflammable materials in the direct vicinity of the device.
- 5.** Never touch the mains plug with wet hands.
- 6.** Only use the plug-in power supply provided; other power supply units can cause damage to and failure of the device.
- 7.** Authorised operators: The owner of the machine must provide the operator with the operating instructions and ensure that he has read and understood them. Only then may the operator commission the device.
- 8.** Prior to each operating cycle, the device must be examined to ensure its impeccable condition and operational safety. If it does not display an impeccable condition, the device may not be used and must be marked accordingly.
- 9.** Do not insert any objects into the housing interior.
- 10.** Do not lock the device operating elements in place.
- 11.** Never look directly at the LED lights.
- 12.** The device should be disconnected from the mains supply for longer periods of non-use.
- 13.** Signs and labels must be kept clearly legible and may not be removed.
- 14.** The mains plug must be removed before cleaning and maintaining the device or replacing parts.
- 15.** The device may only be opened and repaired by authorised specialists.
- 16.** Only accessories and spare parts approved by the manufacturer may be

used. We do not accept any liability for damage incurred by using parts which have not been approved.

17. Independent conversions and modifications are not permissible for safety reasons.

18. It is imperative that the above operating and maintenance conditions in these operating instructions are observed. The general accident prevention regulations must be observed when working with the EyeEvolution®.

5. Commissioning

Insert the mains cable in the device socket (1) to establish a mains connection. The power supply unit is supplied with the EU connector plug. The plug adapters supplied for the UK, Japan and the US can be replaced at the mains device at any time.

Do this by pressing the push-button on the inside of the power supply unit and pulling the EU connector unit out upwards. Insert the corresponding plug adapter into the guide groove. Slide it forwards until it audibly clicks into place.


Only use the plug-in power supply provided for connecting to the mains supply; other power supply units can cause damage to and failure of the device.

Important: Read the operating instructions thoroughly before commissioning the device for the first time. If there are any interpretative problems, please contact the responsible depot or Dreve Dentamid GmbH directly.

Scope of supply:


- 1x EyeEvolution® light-curing device
- 1x 12 V plug-in power supply unit
incl. EU connector
- 1x plug adapter for UK
- 1x plug adapter for US/JAPAN
- 1x operating instructions in
5 languages

6. Functional elements

No. Description

- 1 Power supply connector socket
- 2 Main switch
- 3 Display foil
- 4 Time programme push-button
- 5 Rotary plate function push-button
- 6 Function display
- 7 Rotary plate
- 8 Reflective foil
- 9 LED circuit board
- 10 Glare protection cap
- 11 Cover handle
- 12 Device cover
- 13 Inspection window

7. Operation

Activate the main switch (2) on the back of the device; the functional display LED (6) lights up in red. The device is now ready for operation.

Open the polymerisation chamber by raising the device cover (12) at the handle (11). The opening angle is approx. 55° and is fixed in its final position by a magnet.

Position the item(s) to be polymerised on the reflective area (8) of the rotary plate (7).

Release the cover from its magnetic anchoring by pulling slightly on the handle and close the polymerisation chamber.

7.1 Time programmes

The EyeVolution® light-curing device features several pre-programmed time relays. Use the push-button (4) to activate time programmes with 1, 2, 3, 5 and 10 minutes.

Start the selected time programme by pressing the corresponding push-button (4). The polymerisation process starts automatically and is indicated as an active programme by a green LED.

7.2 Rotary plate

The platform for supporting the model is designed as a reflective rotary plate (7) for full and shadow-free polymerisation wherever possible. The rotary function is coupled with the time programmes and is automatically activated when the polymerisation programme starts.

The LED (6) indicates the operating status of the rotary plate:

Green LED = rotary plate in action

Red LED = rotary function deactivated

Use the push-button (5) to activate or deactivate the rotary plate function at any time, even during on-going polymerisation programmes.

7.3 Programme sequence

Important: In order to protect your eyes, polymerisation is only possible when the device cover is closed.

The light used is not suitable for direct eye contact and can cause retinal damage over the long term when unfiltered. The cover features a light-filtering inspection window (13) to facilitate visual inspection of the polymerisation process.

An on-going time programme is always interrupted by opening the cover (12);

the LED assigned to the respective programme flashes. The time remaining is saved in the electronics module and continues once the polymerisation chamber is closed again.

Once a time programme is complete, 2 brief signal tones sound and the polymerisation LEDs and rotary motor switch off automatically. The green LED indicating the programme on the operating foil (3) is also deactivated. A ring of 6 white LEDs stays active and is only dimmed down after the device cover is opened.

7.4 Premature programme abort

An on-going time programme can be ended and deleted at any time by press-

ing the activated programme key.

8. Cleaning and maintenance

Always remove the mains plug from the mains power supply before carrying out any cleaning or maintenance work. Ideally, the device is cleaned using a soft cloth, possibly a slightly dampened sponge and a mild cleaning agent.

Water and cleaning agent must not penetrate the device.

The LED circuit board (9) must not come into contact with damp materials. Minor soiling can be removed using a soft brush.

9. General instructions

Repairs and maintenance work on the EyeVolution® polymerisation device may only be carried out by qualified personnel. In any case, the power supply to the device must be interrupted. The Dreve Group is only responsible for the safety, reliability and performance of the device when

- the device is only used for the corresponding designated application
- operation is only by qualified personnel.

- the device is used in accordance with these operating instructions.
- extensions, new settings, modifications or repairs are only carried out by persons authorised by Dreve.

As we constantly strive to develop our products, we reserve the right to make technical modifications.

10. Warranty

The warranty period is 2 years as of the delivery date. The warranty only applies for material or manufacturing errors. Warranty claims do not apply for natural wear or damage incurred after the transfer of risk as a result of incorrect or non-designated handling, excessive strain, unsuitable operating equipment and such chemical, electro-chemical or electrical influences which are not intended by the contract. Warranty claims are excluded by any non-designated modifications made by the customer, third parties or maintenance and repair work as well as in the event of violation of lead seals.

Responsibility for defects

Claims for defects shall expire one year after delivery. Claims for defects do not apply for natural wear or damage incurred after the transfer of risk as a result of incorrect or non-designated handling, excessive strain, unsuitable operating equipment and such chemical, electro-chemical or electrical influences which are not intended by the contract. Claims for defects are excluded by any non-designated modifications made by the customer, third parties or maintenance and repair work as well as in the event of violation of lead seals.

11. Troubleshooting

Fault	Cause	Troubleshooting
Device shows no function	Incorrect mains connection	Check mains connection
	Main switch off	Switch on main switch
	Internal cabling is defective	Internal cabling must be checked by a specialist.
LED lights up, electronics without function	Plug-in power supply unit is defective	Replace the plug-in power supply unit
	Electronics are defective	Electronics must be checked by a specialist
	Function deactivated	Activate with pressure switch (5)
Rotary plate without function	Device cover is open	Close cover (12)
	Cable connection to rotary motor (16) is defective	Examine the cabling
	Rotary motor (16) is defective	The motor must be examined by a specialist, and replaced if necessary.
	Electronic module is defective	The electronics must be examined by a specialist, and replaced if necessary.
	Motor adapter (25) rotates continuously	Tighten the motor adapter threaded pin
Time programmes can not be activated	Main switch off	Switch on the device
	Device cover is open	Close cover (12)
	Internal cabling is defective	Internal cabling must be checked by a specialist


Fault	Cause	Troubleshooting
	LED circuit board (9) is defective	The LED circuit board must be examined by a specialist, and replaced if necessary.
	Electronic module is defective	The electronics must be examined by a specialist, and replaced if necessary.
	Micro-switch (15) is defective	The micro-switch must be checked by a specialist and replaced if necessary.
Polymerisation time does not end automatically	Electronic module is defective	The electronics must be examined by a specialist, and replaced if necessary.
Polymerisation process is not interrupted when the cover is opened	Both micro-switches (15) are defective	Deactivate the device immediately. The micro-switch must be checked by a specialist and replaced if necessary.
Unsatisfactory polymerisation results	Polymerisation time is too short	Extend the polymerisation time
	Material reacts to fluctuating wavelengths	Replace the material
	LED circuit board (9) is soiled	Clean the LED circuit board carefully
	LED circuit board is damaged	Examine the LED circuit board, replacing if necessary

12. List of spare parts

No.	Designation	Art. no.
1	Power supply connector socket	58516
2	Main switch	50890
3	Adhesive display foil	59086
7	Rotary plate with reflective foil	D-6559063
8	Reflective foil	D-65590611
9	LED circuit board incl. fixing material	D-6559051
10	Glare protection cap incl. fixing material	D-6559064
11	Handle incl. fixing material	D-6559065
13	Inspection window	59085
14	Power supply unit / Device connector international	59053
15	Micro-switch with cable and fixing material	D-6559055
16	Motorised rotary plate	D-6559054
17	Electronics incl. fixing material	D-6559052
18	Housing (cover, handle and housing base)	D-6559061
19	Retention disc magnet incl. fixing material	D-655573
20	Metal base plate incl. fixing material	D-6550809
21	Device feet (5x)	D-6559060
22	Cover hinge (micro-switch, cable, cover plate, cover axle, disc magnet) incl. fixing material	D-6559083
23	Cooler plate	50815
24	Motor flange	50810
25	Motor adapter incl. threaded pin	D-6559081
26	Cover plate incl. fixing material	D-6559059
27	Cover axle	59083

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