20240320-3999

5. **Notification**
The specifications have been validated using the printer software (Autodesk Netfabb®, Composer, Dreve Elements) and the previously mentioned printers and process parameters indicated.

- Caution:**
- To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and/or color variations of the FotoDent® biobased model material.
 - Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® biobased model. It is recommended to regularly take out the container of the production unit to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hours with FotoDent® biobased model.

Processing

- Preparation:**
- During the application of FotoDent® biobased model we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
 - FotoDent® biobased model should be well shaken approx. 1 hour before use.
 - Carefully pour FotoDent® biobased model into the pre-determined container of the production unit.
 - Please see further steps within printer manual.
 - Remove all bubbles with a cleaned object.

Processing temperature	
MAX UV, PRO 4K	35 ± 3 °C
Rapid Shape DII-series, D20+, D30+	23 ± 3 °C
Dreve Phrozen Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)	23 ± 3 °C

- Select the corresponding profile for FotoDent® biobased model in the printer software. Make sure that all software is up to date.
Asiga:
for 50 µm:
FotoDent® biobased model 385nm_D35450_v2.2_50.ini
for 100 µm:
FotoDent® biobased model 385nm_D35450_v2.2_100.ini
Rapid Shape:
for 100 µm:
FotoDent® biobased model #D35450

Dreve Phrozen:
for 100 µm:
Dreve FotoDent biobased model beige v4_1
Dreve FotoDent biobased model beige v4_2
Dreve FotoDent biobased model beige v4_3
for 50 µm:
Dreve FotoDent biobased model beige v4_4_low
Dreve FotoDent biobased model beige v4_4_med
Dreve FotoDent biobased model beige v4_4_xtr
Further information at https://phrozen.dreve.de/profiles_en.

- Further processing:**
- After completion of the building process, direct further processing (see following steps) is recommended.
 - After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
 - Resin coated parts can be cleaned with FotoClean in a heated ultrasonic bath. Please follow the FotoClean instructions for use. Make sure the cleaner is warmed up to 30–35 °C. The cleaning time is 6–8 minutes and the components should then be removed directly from the ultrasonic bath. After cleaning in the ultrasonic bath, the components should be rinsed in a container with at least 2 liters of 23 ± 3 °C warm water for 1 minute by swirling. The components can then be airdried or dried with compressed air.
 - Alternatively, resin coated parts can be cleaned with Isopropanol (97 %) within a time of 2 x 6 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended.
 - For post-curing of the parts please select the parameter setting for your process. You can find more qualified printers on our website www.dentamid.dreve.de

Recommended curing light equipment (Post curing units)	
Manufacturer, Model	PCU LED N ₂ (Dreve)
Supply voltage	100–240 V
Lamp Power	N/A
Light intensity	Adjustable
Wavelength	410 nm

Post curing	
PCU LED N ₂	
Asiga MAX UV, PRO 4K	12 minutes, 80 % light power using nitrogen
Rapid Shape DII-series, D20+, D30+	12 minutes, 80 % light output using nitrogen
Dreve Phrozen Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)	12 minutes, 80 % light output using nitrogen

Disposal

Disposal of contents/containers in accordance with local/regional/national/international regulations and according to the safety data sheet.

Hazard identification

Hazard and safety information can be found on the product label and/or in the corresponding safety data sheet.

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® biobased model packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

Glossary


	Manufacturer		Protect from sunlight
	Batch Code		Temperature limitation
	Item number		Expiry date
	Use by only professionals		Please note instructions for use
	Nanoparticles contained		Made in Germany

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Product description EN

Light-curing material for the fabrication of dental working models.

Intended users: The dental technicians or specialist staff in the dental laboratory is responsible for the production of dental models in 3D printing. The dental technicians or specialist staff use the finished dental models only in the laboratory.

Requirements:

- Printer:
 - Rapid Shape: DII-Serie, D20+, D30 (385 nm)
 - Asiga: MAX UV, PRO 4K (385 nm)
 - Dreve: Phrozen Sonic XL 4K (qualified by Dreve) / Phrozen Sonic XL 4K 2022 (qualified by Dreve) (405 nm)
- Software: Autodesk Netfabb®, Composer
- Post curing unit: PCU LED N₂

Basic Material Properties: Please refer to Technical Data Sheet (TDS). Available on request.

Specific Manufacturing Considerations:

- Digital dental model file**
 - File format: STL file
 - Digital design: dental models
- Asiga (MAX UV, PRO 4K) and software**
 - Hardware: see directions for use in the printer manual
 - Software: see directions for use in the printer manual
 - Printingparameters:
 - Slice thickness: 50 µm, 100 µm
 - Optimal orientation: horizontal or vertical
 - Support point size: no supports are needed
 - Environmental Conditions
 - Print room temperature: 35 °C ± 3 °C
- Rapid Shape (DII-series, D20+, D30+) and Software**
 - Hardware: see directions for use in the printer manual.
 - Software: see directions for use in the printer manual.
 - Printig parameters:
 - Slice thickness: 100 µm
 - Optimal Orientation: horizontal or vertical, directly on the platform
 - Support point size: no supports are needed
 - Environmental Conditions
 - Print temperature: 23 ± 3 °C
- Dreve Phrozen Sonic XL 4K (qualified by Dreve) / Phrozen Sonic XL 4K 2022 (qualified by Dreve) und Software**
 - Hardware: Firmware version XL22_Dreve_V1.62 or higher
 - Software: Dreve ElementS version 1.0.9.7 or higher
 - Printig parameters:
 - Slice thickness: 50 µm, 100 µm
 - Optimal Orientation: horizontal or vertical, directly on the platform
 - Support point size: no supports are need
 - Environmental Conditions:
 - Print temperature: 23 ± 3 °C