Dreve Design Rev.070319-3933

• For post-curing of the parts, please select the parameter settings for your process. You can find more qualified printers on our website www.dreve.dentamid.de.

• Post-curing
PCU LED N₂
Rapid Shape DII / Asiga MAX UV 8 minutes, 80 % light power (vacuum)

FotoDent® flash
Rapid Shape DII / Asiga MAX UV 10 minutes + 1 minute gas pre-flow

RS Cure
Rapid Shape DII / Asiga MAX UV 8 min, 75 % light power (vacuum)

• Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® gingiva. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hour with FotoDent® gingiva.

To avoid damage of material quality does not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and/or colour variations of the FotoDent® gingiva material.

Possible dangers

Please note the safety instructions for processing FotoDent® gingiva:

Hazard statements: **H319** Causes serious eye irritation. **H317** May cause an allergic skin reaction. **H411** Toxic to aquatic life with long lasting effects.

Precautionary statements: **P261** Avoid breathing dust/fume/gas/mist/vapours/spray. **P273** Avoid release to the environment. **P280** Wear protective gloves/eye protection/face protection. **P308+P313** IF exposed or concerned: Get medical advice/attention. **P333+P313** If skin irritation or rash occurs: Get medical advice/attention. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: Aliphatic difunctional Urethane Acrylate, Urethane acrylate in HPMA, phenyl bis(2,4,6-trimethylbenzoyl)-phosphine-oxide, Pigments

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® gingiva packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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Product description

FotoDent® gingiva is a light-curing resin for the production of gingival masks using 385 nm LED based stereolithography. For obligatory post curing of FotoDent® gingiva built parts, the high performance light curing units FotoDent® Flash and PCU LED N₂ are recommended, which are optimal balanced to the production process.

Processing

- During the application of FotoDent® gingiva we recommend wearing personal protective equipment f. i. suitable gloves, safety glasses etc.
- FotoDent® gingiva should be well shaken approx. 1 hour before use.
- Please see further steps within printer manual.
- Pour FotoDent® gingiva carefully into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.
- Processing temperature:
Rapid Shape DII printer: 23 ± 3 °C
Asiga MAX UV: 35 ± 3 °C
- Select the product specific parameter settings for FotoDent® gingiva for your printing process.
- When building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 2x6 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended. Alternatively, cleaning with isopropanol can be carried out in a separate beaker under pouring conditions for approx. 2x6 minutes. Recommendation: Beaker and magnetic stirrer, laboratory shaker or ultrasonic cleaner. In case of insufficient cleaning (such as shiny areas on the object after drying) repeat cleaning process.

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Dreve Dentamid GmbH