- For post-curing of the parts, please select the parameter settings for your process. You can find more qualified printers on our website [www.dentamid.dreve.de](http://www.dentamid.dreve.de).

Post curing FotoDent® model2 beige-opaque	
PCU LED N <sub>2</sub> Rapid Shape DII / Asiga MAX UV Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)	8 minutes, 80 % light power (vacuum)
FotoDent® flash Rapid Shape DII / Asiga MAX UV	8 minutes + 1 minute gas pre-flow / 12 minutes + 1 minute gas pre-flow
PCU 90	8 minutes (incl. 4 minutes gas pre-flow)
RS Cure	According to respective parameter set for FotoDent® model2
Post curing FotoDent® model2 black	
PCU LED N <sub>2</sub> Rapid Shape DII / Asiga MAX UV Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)	12 minutes, 80 % light power (vacuum)
RS Cure	According to respective parameter set for FotoDent® model2
Post curing FotoDent® model2 white	
PCU LED N <sub>2</sub> Rapid Shape DII / Asiga MAX UV Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)	6 minutes, 40 % light power (vacuum)
RS Cure	According to respective parameter set for FotoDent® model2

- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® model2. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hour with FotoDent® model2.

**To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and/or colour variations of the FotoDent® model2 material.**

### Hazards identification

Hazard and safety notes can be found on the product label and the corresponding safety data sheet.

### Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® model2 packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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## Product description

Light-curing material for the fabrication of dental working models.

## Processing

- During the application of FotoDent® model2 we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- FotoDent® model2 should be well shaken approx. 1 hour before use.
- Carefully pour FotoDent® model2 into the pre-determined container of the production unit.
- Please see further steps within printer manual.
- Remove all bubbles with a cleaned object.

Processing temperature	
Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)	23 ± 3 °C
Rapid Shape DII printer	23 ± 3 °C
Asiga MAX UV	35 ± 3 °C

- Select the product specific parameter settings for FotoDent® model2 for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 2x6 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended. Alternatively, cleaning with isopropanol can be carried out in a separate beaker under pouring conditions for approx. 2x6 minutes. Recommendation: Beaker and magnetic stirrer, laboratory shaker or ultrasonic cleaner. In case of insufficient cleaning (such as shiny areas on the object after drying) repeat cleaning process.