Dreve Design Rev.070319-3915

• For post-curing of the parts, please select the parameter settings for your process. You can find more qualified printers on our website [www.dreve.dentamid.de](http://www.dreve.dentamid.de).

• PCU LED N<sub>2</sub>  
Rapid Shape D, DII / Asiga MAX UV 8 minutes, 80 % light power (vacuum)

FotoDent® flash  
Rapid Shape D, DII / Asiga MAX UV 12 minutes + 1 minute gas pre-flow

Otoflash G171  
Rapid Shape D, DII / Asiga MAX UV 6000 flashes (incl. nitrogen and UVB blocker cup)

RS Cure  
Rapid Shape D, DII / Asiga MAX UV 8 minutes, 75 % light power (vacuum)

• Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® model. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hour with FotoDent® model.

**To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and/or colour variations of the FotoDent® model material.**

### Possible dangers

Please note the safety instructions for processing FotoDent® model.

**Hazard statements: H317** May cause an allergic skin reaction. **H412** Harmful to aquatic life with long lasting effects.

**Precautionary statements: P261** Avoid breathing dust/fume/gas/mist/vapours/spray. **P273** Avoid release to the environment. **P280** Wear protective gloves. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P333+P313** If skin irritation or rash occurs: Get medical advice/attention. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

**Contains:** 1,4 Butanediol dimethacrylate, Diurethandimethacrylate, hydroxypropylmethacrylate, monoester with propane-1,2-diol, Diphenyl(2,4,6-trimethylbenzoyl)phosphinoxide

### Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® model packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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## Intended use

FotoDent® model is a material based on (meth-)acrylate resin for DLP systems with 385 nm / 405 nm LED for the production of dental models. For the post curing of FotoDent® model built parts, the light curing units PCU LED N<sub>2</sub>, FotoDent® flash, Otofash G171 and RS Cure are recommended, which are optimal balanced to the production process.

## Processing

- During the application of FotoDent® model we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- FotoDent® model should be well shaken approx. 1 hour before use.
- Please see further steps within printer manual.
- Carefully pour FotoDent® model into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.
- Processing temperature:  
Rapid Shape D and DII printer 23 ± 3 °C  
Asiga MAX UV 35 ± 3 °C
- Select the product specific parameter settings for FotoDent® model for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 2 x 6 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended. Alternatively, cleaning with isopropanol can be carried out in a separate beaker under pouring conditions for approx. 2 x 6 minutes. Recommendation: Beaker and magnetic stirrer, laboratory shaker or ultrasonic cleaner. In case of insufficient cleaning (such as shiny areas on the object after drying) repeat cleaning process.

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