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Intended use

Flexible gingiva mask material for the C&B technique

Material description

GumQuick scan is an addition-vulcanising vinylploisiloxane for the production of flexible gingiva masks on the master model.

GumQuick scan reaches a final hardness of approx. 63 Shore A and adapts the colour naturally. The material is characterised by its high tear-resistance and is easy to treat with rotating instruments.

Field of application

Gingiva masks produced in the laboratory for the C&B technique and for high-precision implant works.

Processing

Indirect method

A matrix of hard kneadable silicone material (e. g. Laborsil) is to be applied onto the dies of the unsawn die model which should be provided with the gingiva mask. After curing the kneadable silicone should be cut so that the top of the dies will appear (approx. 1 mm). The model can be sawn and gingiva parts ground as usual. Afterwards the correct position of the matrix has to be checked. In order to guarantee a faultless surface of the gingiva mask DuoSep is to be applied onto the contact surface of the kneadable silicone / GumQuick scan.

GumQuick scan is to be applied directly from the cartridge onto the model and into the matrix. Now the matrix has to be positioned at the model so that the surplus of material can flow out through the holes at the top of the die. Ideally, the gingiva mask should be cured under pressure.

After a period of 5–6 min. the matrix can be removed and the mask can be treated with scissors, scalpel or with rotating grinding instruments.

Direct method

At first DuoSep has to be applied onto the contact surface impression material / GumQuick scan in order to guarantee a faultless surface of the gingiva mask.

Attention

In order to avoid a damage of the gingiva mask surface when applying GumQuick scan in conjunction with polyether impressions DuoSep has to be used! Before the production of the gingiva mask the polyether impression should ventilate for at least 3 hours.

Afterwards GumQuick scan is to be applied directly from the cartridge into the impression. After a period of 5–6 min. GumQuick scan will be cured and the im-

pression can be filled with plaster. The gingiva mask can easily be treated with scissors, scalpels or with rotating grinding instruments.

Important advices

- Due to technical filling reasons an unequal transportation of the components can occur when using a new cartridge. For this reason the first transported quantity – approx. size of a pea – should not be used.
- Never try to squeeze the material through a cured mixing canula. This may lead to leakage in the cartridge system and therefore to transportation of un-mixed and uncured single components.
- Before each application the cartridge has to be checked for consistent transportation of material.
- The vulcanisation can be disturbed by the use of latex gloves, hand creams, cleaning agents etc.
- The product is exclusively designed for the above mentioned field of application and should only be used by qualified and instructed persons.

Storage conditions

At room temperature, i. e. 64.4 °F–82.4 °F (18 °C–28 °C)

Batch number / use-by date

The batch number and the use-by date are located on both the outer packaging and the product. Both components A and B of a packaging must only be used in this combination. When submitting a complaint about a product, please always quote the batch number. Do not use the product after expiration of the use-by date.

This product was developed for dental use and has to be applied according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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