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Addition-curing duplication silicone especially for use in the resin pouring technique

Product description

Siliform is an embedding material on the basis of A-silicone especially for the resin pouring technique. This material is especially characterised by the approved practical mixing ratio of 1:1, its brilliant reproduction of details and the lowest shrinkage guarantee an exact reproduction of the wax modellation. In addition, it has a further advantage by which it is distinguished considerably from previous materials: After vulcanisation is finished the material shows as requested a humid surface. This humidity creates a film between the silicone and the resin during the later polymerisation whereby a smooth, silky-mat surface of the denture body is created – known from the use of the duplication gel (**Gelon**). Anymore, this surface structure improves the adhesion between tooth and tooth gap when putting back the teeth into the silicone counter-part later.

Technical data according to DIN EN ISO 4823

Mixing time (manually): 45 sec.
Processing time (incl. mixing): ≥ 6 min.
Deflasking (approx. 8 shore): ≥ 35 min.

Preparatory measures

With **Siliform** the ready modelled set-up has to be duplicated together with the plaster model. The master model must not be watered. **Siliform** slightly tends to sediment. Stir or shake set material components before use.

Dosage

Siliform is mixed in the ratio 1:1 of the components A (catalyst / white) and B (base / coloured).gemischt.

Processing

Shake material well before processing. Mix both components in a mixing beaker or by means of the many dosing and mixing aids available (e. g. dosing dispenser).

The mixing and processing time refer to a room temperature of 23 °C (73.4 Fahrenheit) and a relative air humidity of 50 %.

Duplication

If the material is mixed manually, the mass should be poured into the flask in a thin stream. The distance should be approx. 20 cm.

Further processing

The best results regarding homogeneity can be achieved with the pressure duplication system. The hardening of the silicone should be effected under the same pressure (e.g. **Polymax**) which will be used for the curing of the resin later on. This leads to a safer reproduction of the negative form without deformations.

General Notes

Siliform is not suitable for the duplication of plaster models or for the production of embedding masses models. Here we recommend our **Dublisil 15**. The vulcanisation can be disturbed by latex gloves, hand cremes, cleaning agents etc. The product is exclusively designed for the a. m. field of application and is only to be used by qualified and instructed persons.

The general information of this product is based on detailed application practice and intensive research work. Due to further development we reserve the right of modification.

Keep away from children! Only for the a. m. field of application!

Storage conditions

At room temperature, i. e. 18–28 °C.

Durability: see date of expiry.

Do not use **Siliform** after expiration.

Delivery forms	
[REF] D38605	Siliform 2 x 850 ml bottles A+B

This product was developed for dental use and has to be applied according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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