Betriebsanleitung
Operating manual
Mode d'emploi
Instrucciones de operación
Istruzioni per l'uso

deutsch

english

français

español

italiano

Vacformat V9 2.0

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Vacuum thermoforming unit for aligner technique

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EN Contents

1. Device description	24
2. Intended use	24
3. Technical data	24
4. Declaration of Conformity	25
5. Safety information	26
6. Commissioning	28
7. Functional elements	29
8. Operation	30-34
9. Deviations from regular operation	34
10. Exchanging the device fuse	35
11. Resetting the temperature monitor	35
12. General information	36
13. Care and maintenance	36
14. Liability for defects	36
15. Decommissioning of the product	37
16. Troubleshooting	38
17. Exploded view	40
18. Spare parts list	41
19. Wiring diagram	96

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1. Device description

The Vacformat V9 2.0 is a thermoforming unit on the basis of vacuum technology, which was especially designed for the production of aligner splints. In just one single operation, up to 9 aligners can be thermoformed at the same time.

Very simple handling as well as continuous digital and acoustic control of all working steps ensure problem-free working, even with little prior knowledge.

The Vacformat V9 2.0 has been developed especially for Biolon foils with a material thickness of 0.5–1.0 mm and, thanks to its vertical forming as well as a very powerful vacuum pump with integrated vacuum reservoir, generates extremely precise final results. The product was developed for use in the dental sector and must be used according to the user information. For damage resulting from failure to observe the instructions for use or from improper use, the manufacturer will not accept any liability. Moreover, the user is obliged to check the product on their own responsibility for its suitability and possible use for the intended purposes before using it, especially where these purposes are not specified in the user information.

2. Intended use

The Vacformat V9 2.0 is a thermoforming unit for thermal vacuum forming and can be used for the simultaneous production of nine aligners in one operation by a trained user. With the Vacformat V9 2.0 you can produce only aligners or retainers made from Biolon foils (PETG) with a thickness of 0.5–1.0 mm. The dental models used for this may have a maximum height of 20mm and must be made of a short-term resistant material up to 150 °C.

3. Technical data

Dimensions (H x W x D)	415 x 487.6 x 600.7 mm
Weight	45 kg
Mains voltage / input	220–240 V, 50 / 60 Hz
Max. power consumption	2500 W
Heating	medium infrared radiation
Device fuse	2 x 12 A
Thermoforming capacity	Up to 9 aligner models at the same time
Forming	Vacuum pump inclusive vacuum reservoir (until 2,9 psi)
Foil intake:	354 x 354 mm up to max. 1.0 mm
Maximum humidity	80 %
Max. working height	2000 mm
Operating temperature	10°C–30°C
Overvoltage category	II

4. EC-Declaration of Conformity

The manufacturer/distributor

Dreve Dentamid GmbH • Max Planck Str. 31 • 59423 Unna

hereby declares that the following product

Product designation:	Vacformat V9 2.0
Model name:	D3282
Type name:	Laboratory device / Vacuum thermoforming
Trade name:	Vacformat V9 2.0

Description:

Thermoforming device for the parallel production of up to nine dental aligner splints using vacuum

fulfills all the relevant provisions of the legislation (in the following) - including the changes which applied at the time of the declaration. This declaration of conformity is issued under the sole responsibility of the manufacturer. This declaration only refers to the machine in the condition in which it was placed on the market; it does not apply to parts installed by the user and/or modifications made at a later time.

The following legislation have been applied:

Machinery Directive	2006/42/EG
EMC Directive	2014/30/EU
RoHS Directive	2011/65/EU

The protection objectives of the following further legislation are observed:

Low-Voltage Directive	2014/35/EU
------------------------------	------------

The following harmonised standards have been applied:

ENIEC 61000-3-2:2019	Electromagnetic compatibility (EMC) – Part 3-2: Limits – Limits for harmonic current emissions (equipment input current = 16 A per phase)
EN 61000-3-3:2013	Electromagnetic Compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage power supply systems for devices with a rated current ≤ 16 A per phase, the are not subject to any special connection conditions (IEC 61000-3-3:2013)
EN 61010-1:2010/ A1:2019/AC:2019-04	Safety regulations for electrical measurement, control, regulation and Laboratory equipment – Part 1: General requirements (IEC 61010-1:2010)
EN 61010-2-010:2014	Safety regulations for electrical measurement, control, regulation and Laboratory equipment – Part 2-010: Particular requirements for laboratory equipment for heating substances (IEC 61010-2-010:2014)
EN 61326-1:2013	Electrical equipment for measurement, control and laboratory use – EMC requirements – Part 1: General requirements (IEC 61326-1:2012)

Dreve Dentamid GmbH
Max-Planck-Straße 31 • 59423 Unna/Germany • www.dentamid.dreve.de

5. Safety instructions

Attention! Please read these notes carefully before the connection and commissioning of the device. The operational safety and functioning of the device can only be ensured if both the general safety and accident prevention regulations as specified by the legislator and the safety instructions in the instructions for use are observed.

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**Do not reach into the interior of the device.
Risk of burns!**

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Unplug the power cord before cleaning and maintenance work. Risk of electric shock.

Achtung heiß!
Attention hot!

**The device heats up during operation.
Danger of burns!**

1. The device must only be used in accordance with the present instructions for use and by technically instructed personnel. We do not assume any liability for damages which result from improper use and/or incorrect operation.
2. The voltage specified on the name-plate must match the voltage of the power source.
3. The power supply of the Vacfomat V9 2.0 should only be established by means of an especially secured power supply line.
4. The device may only be connected to a socket with a protective conductor.
5. The device must be placed on a stable and level surface which is sufficiently stable for a weight of 50 kg.
6. To prevent the ingress of water (e.g. splash water) into the device, the device should be installed in a dry environment.
7. Do not keep any highly inflammable substances within the immediate vicinity of the device.
8. Authorised operators: The operator of the machine must make the instructions for use accessible to the operator and must also make sure that s/he has read and understood them. Only then may the operator put the device into operation.
9. Please check the device to ensure its proper condition and operating safety before each operation. If its condition is not flawless, the device must not be used anymore and must be labelled accordingly.
10. Avoid direct eye contact with the heating tubes. Infrared beams may cause permanent eye damage!
11. Signs and stickers must always be kept in an easy-to-read condition and must not be removed.
12. The three (3) heating elements of the device reach high temperatures during the heating phase and, therefore, must not be touched, not even for a short time!
13. If the device is used continuously, the surfaces of the heating cover as well as other outer parts of the Vacfomat V9 2.0 may heat up to a temperature of approx. 45 °C.
14. Do not lock the operating elements of the device.
15. Do not insert any objects into the housing.
16. Cleaning and maintaining the device or exchanging parts should only be carried out after the device has cooled down. Before doing this, the mains plug must be pulled out.
17. Never touch the mains plug with wet hands.
18. The opening of the device and repairs may only be performed by authorised professionals.
19. Only accessories and spare parts may be used which are approved by the manufacturer. We will not accept any liability for any damage caused by the use of third-party parts.
20. For safety reasons, unauthorised modifications and changes to the Vacfomat V9 2.0 are not permitted and any warranty claims will expire immediately.
21. The prescribed operating and maintenance conditions specified in these instructions for use must be complied with. Please observe the general accident prevention regulations when working with the Vacfomat V9 2.0.
Since we are constantly working on the further development of our products, we cannot exclude the possibility of technical changes.
22. The device may only be carried by at least two people. Appropriate safety equipment (shoes) must be worn.
23. When setting up the device, make sure that the on/off switch is operable.

Caution! There is a residual risk of burns should one reach into the unit and touch the radiant heaters.

6. Commissioning

Before using the device for the first time, you must have read the user information thoroughly. In case of any uncertainties, please contact the responsible dealer or Dreve Dentamid directly.

Carefully remove the Vacfomat V9 2.0 from its transport packaging and place it on the intended worktop. It should have a minimum depth extension of > 600 mm and be sufficiently stable to accommodate a weight of approx. 50 kg. The ambient temperature at the installation site of the Vacfomat V9 2.0 should not permanently exceed a temperature of 23 °C. Before establishing the power supply to the mains, make sure that all 3 transport locking devices have been removed from the device and that the lever mechanism for lifting and lowering the foil cartridge (10) is fully functional. Check the capacity of your power connection. The Vacfomat V9 2.0 is operated with 220–240 V, 50/60 Hz and has a max. power consumption of 2500 W. Insert the supplied power cord into the device socket and establish a mains connection. Make sure that the on/off switch is easily accessible.


Only use one of the supplied power connection cables for the mains connection.

Scope of delivery:

- 1x Vacfomat V9 2.0
- 2x power connection cable (EU and US version)
- 1x Instructions for use in (5) languages
- 1x magnetic tongs

7. Functional elements

No Designation

- | | |
|------------------------|------------------------------|
| 1. Device socket | 7. Vacuum indicator |
| 2. Main switch | 8. Operating lever |
| 3. Rotary/push button | 9. Removable model reception |
| 4. Display | 10. Cassette for blanks |
| 5. Start / Stop button | 11. Cassette holder |
| 6. Function indicators | 12. Fuse 12A |

8. Operation

Turn on the main switch (2) on the left side of the device housing. This activates the control light on the main switch, as well as the display (4) on the front side. If a new device has been delivered, the display will show 0:00. Otherwise, the heating time which was used last will be called up. An acoustic signal shows that the device is ready for operation.

8.1 Setting the heating time

Any time between 0:05 min. and 4:00 min. can be set by turning the rotary/push button of the time setting (3). The setting is made in steps of five (5) sec. When setting a new heating time, the information shown on the display flashes and must be confirmed before using the device by pressing the rotary/push button (3). For the optimum values for the Dreve aligner thermoforming programme, please refer to the table below or the information printed on the foil packaging.

Continuous operation of the device in heating mode is not possible. The control of the device only allows heating and cooling phases of the same length


Material / Times	Heating min/sec	Cooling min/sec
Biolon 0,5 mm	1:35	1:35
Biolon 0,625 mm	1:35	1:35
Biolon 0,75 mm	1:40	1:40
Biolon 1,0 mm	1:45	1:45



To ensure that the specified heating times of the respective foils also provide an adequate plasticisation degree for good forming, the device must be pre-heated every morning or after a longer time out of use for at least 0:30 min. without a foil. If the V9 2.0 is used on power grids other than those specified connected, the heat output of the heating system changes. Depending on the network strength it may be necessary to interpose a mains converter to 230V / 50 Hz or adjust the heating times accordingly.

8.2. Preparatory measures

The cartridge holder (11) can be moved between two positions by means of the lateral operating levers (8).


- Remove the foil cartridge (10) from the device in lever position A by using the handle to pull it forwards out of the cartridge holder (11) and placing it on a free table surface.
- Open the lid of the foil cartridge (8) at the handle and insert the Biolon foil into the bottom part with the insulating layer facing downwards. To safely fix the thermoforming foils, retaining pins are installed in each corner of the foil cartridge (8). The V9 foils have corresponding holes in the same positions. Slide the foil over the retaining pins and close the lid.

- Leave the cassette holder (11) in the heating position A. Remove the model reception (9) from the device using the enclosed magnetic tongs (29). To do this, fit the two ring-shaped fixing elements in the lower part of the pliers on the counter bearing mounted in the model reception (9). The upper part of the pliers slide over the fixing elements. Remove the model reception (9) from the device and place it on the tabletop. Position the models for thermoforming of the model reception

(9). Put back the equipped model reception into the device by using the magnet tongs (29). Remove the tongs from the model reception (9 by pulling back the upper part of the tongs).

Tip:

The aligner models should not exceed a max. height of 19 mm. The best forming results can be achieved if you align the models with the anterior teeth towards the centre of the model support (9).

- Bring the cartridge holder (11) into position A until you notice that the fixing magnets (20) have engaged and push the foil cartridge (10) back into the device until the stop has been reached.
- It is also possible to change the opening direction of the foil cartridge (10) from the left to the right side if required. For this purpose, dismantle the handle (2 cross-head screws) and fix it again at the holes provided for this at the back of the device.

8.3. Start of programme

To start the heating process, press the "Start / Stop" button (5) and, on the display, the preset heating time (s. chapter 8.1) will run down to 0:00 in the countdown mode. The function of the heater is indicated by the white function light (6).

At the same time, the vacuum pump (15) is activated, The air tank (27) which is located in the device already vented to a value of -0.8 bar. When reaching the max. venting value it stops the pumping process. When vacuum drops down to -0.7 bar the pump automatically restarts again.

Due to the large foil surface, the thermoforming foil will sag significantly, partially until it is in contact with the model, when heating up. This will have no negative effects on the quality and material thickness of the finished aligner! For safety reasons, the heating is equipped with a limit (24) switch and only works in the heating position A after a programme has been started.

! If the position of the cartridge holder (11) is changed using the operating lever (8) during the active heating phase, the heating process will be interrupted automatically. The information shown on the display (14) will stop and show the remaining time in flashing mode. By bringing the operating lever (8) back to the heating position A, the programme will be reactivated.

If the Vacformat V9 2.0 is used continuously, the outer housing parts (13 & 14) may heat up to temperatures of up to 45 °C during the heating process! To avoid heating up to even higher temperatures, the Vacformat V9 2.0 is equipped with active ventilation. During the entire heating and forming time, three (21) fans will run automatically. The temperatures both in the lid area and in the forming area are permanently checked by temperature limiters and monitors. If the temperature exceeds 45 °C, the fans will remain (21) in the cooling mode until the measured temperatures have fallen below a value of 40 °C. Active intervention of the operator will not be required. If a temperature of 65 °C has been exceeded, the heating will be switched off automatically for safety reasons. The automatic (18) circuit breaker responsible for this can be reset to room temperature after the coolingdown phase has expired. (S. chapter 11.1.)

8.4. Thermoforming

After the heating time has expired, an acoustic signal will be sounded, the heating time on the display will flash at 0:00 and the heating will remain active for another ten 10 sec. Pull the operating lever (8) up from position A (heating) to position B (forming). The allround silicone seal closes the thermo forming area and at the same time opens the valve of the vacuum reservoir (27) for a sudden molding of the blank onto the models. The vacuum pump (15) switches on automatically when the vacuum decreases and keeps the forming force at a constant level.

Continuing the heating process after the set time has elapsed by a further 10 sec.. This optimizes the occlusal detail during the molding. At the end of the active pre-heating phase, (10 sec.) the white function light (6) is switched off and the cooling-down phase will be shown on the display in seconds in the countup mode.

The maximal value that can be displayed is 5:00 min. For the optimum cooling-down time for the Dreve foil used, please refer to the table above (chap. 8.1) or the information of the material packaging.

8.5. End of the thermoforming process

After the end of the recommended cooling-down time, press the start/stop button (5) and hold it for approx. two 2 seconds. The heating time used last will be shown on the display. Move the cartridge holder (11) to the heating position by bringing the operating lever (8) back to position A. Remove the foil cartridge (10) as described in point 8.2.

Continuous operation of the device in heating mode is not possible. The control of the device only allows heating and cooling phases of the same length

9. Deviations from regular operation

9.1. Changing specified times

The specified times can be prolonged or shortened individually if required. For this purpose, the new time to be entered is changed by means of the rotary/push button (3) of the time setting until the new time is shown on the display (4) in flashing mode. Before it is used for the first time, it must be approved by the operator by pressing the rotary/push button (3). The time shown will no longer flash. The new time programme is now ready for operation.

Tip:

The performance of the heating may decrease slightly after 2000 operating hours. By prolonging the heating times by five to ten (5–10) seconds, optimum thermoforming results can be attained again.

The specified times are average values and intended for general use.

In case of extreme undercut areas or very high models, a prolongation of the heating time can optimise the results.

9.2. Abort / termination of programme

It is possible to abort / terminate all actions of the Vacformat V9 2.0 by pressing the start/stop button (5) for at least two 2 seconds.

9.3. Overheating error

If the "E01" error is shown in the display, the device has become too hot and has switched itself off for safety reasons. To reset the temperature fuse see 11.1

10. Exchanging the device fuse

Disconnect the device from the mains supply by pulling out (12) the mains plug of the power cable (1). The fuse carrier is located on the left side of the device, integrated in the mains supply. In order to check or exchange the fuses, push and hold the push catch of the carrier and remove the carrier from the device.

Pull the defective fuse out of the clamp bracket by applying slight traction and replace it by a new, equivalent fuse.

Only use the fuse strength specified on the back.

After the fuse has been exchanged, insert the fuse carrier into the housing by applying slight pressure until you notice that the fixation has engaged.

11. Resetting the temperature monitor / exchanging an infrared heating element

Before accessing the inner parts of the device, the mains plug must be pulled out.

11.1. Resetting the temperature monitor

The temperature monitor (18) checks the ambient temperatures of the heating tubes and the heating cover (13). If the temperature rises on the inside of the heating cover (13) to a value > 65°C during the heating process, the temperature monitor will be activated and disconnect the heating tubes from the mains for safety reasons. The function indicator for the heating (6) will go out and the ongoing thermoforming process cannot be completed. In order to reactivate the device be sure to disconnect from power supply and allow the unit to cool down to room temperature. In the centre on the back of the Vacformat V9 2.0, there is an opening closed with a plastic cover between the two fans (21). Carefully remove the plastic cover, insert an approx. 6–7 cm long rod (e.g. screwdriver) into the opening and press it against the reactivation button of the temperature monitor (18) at the end until it has engaged again. Then connect the device again to the mains.

11.2. Exchanging an infrared heating element / temperature limiter

The exchange of single heating elements (17) and / or the temperature limiters (18 + 19) may only be done by qualified electricians. If necessary, separate instructions for exchanging these spare parts can be requested from Dreve Dentamid.

12. General information

Repair and corrective maintenance work at the Vacformat V9 2.0 may only be carried out by authorised experts.

The power supply to the device must be disconnected in any case.

Dreve will only be responsible for the safety, reliability and performance of the device if:

- The device is used as intended.
- The operation is only carried out by professionally trained personnel.
- The device is used in accordance with these instructions for use.
- Extensions, new settings, changes or repairs are only performed by authorised persons.

Since we are constantly working on the further development of our products, we cannot exclude the possibility of technical changes.

13. Care and maintenance

Before carrying out any care and maintenance work, the mains plug must always be disconnected from the mains.

The thermoforming area of the device should be checked on a daily basis for contamination caused by material residues. These residues, if any, should be removed. Dirt may enter the vacuum system and damage the pump in the long term.

Ideally, the device should be cleaned dry with a soft cloth and, if necessary, with a slightly damp sponge and a mild cleaning liquid.

Water and cleaning liquids must not get into the device.

14. Liability for defects

Claims for defects will lapse after a period of one year from the date of delivery. The claims for defects do not refer to natural wear nor to any damage which result after the transfer of risk due to incorrect or unintended use, excessive loading, unsuitable equipment and any such chemical, electrochemical or electrical influences which are not provided in accordance with the contract.

Claims for defects will be excluded in the case of changes or repair and maintenance work which have been incorrectly completed by the purchaser or third parties as well as in the case of the breaching of the seals.

15. Decommissioning of the product


The DIRECTIVE 2012/19/EU OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 4 July 2012 on waste electrical and electronic equipment (WEEE) lays down measures to protect the environment and human health. Waste electrical and electronic equipment marked as such may not be disposed of as unsorted municipal waste within the European Union. For detailed information regarding the professional disposal of WEEE please contact Dreve or your distributor.

16. Troubleshooting

Problem	Cause	Aid
Device does not show any function	Mains connection not correct	Check power connection
	Main switch off	Switch on main switch
	Miniature fuse defective	Check miniature fuse Exchange it if necessary
Main switch is lit, electronics do not show any function	Electronics defective Wiring defect	The electronics must be checked by an expert
Display is lit, electronics cannot be set	Previous programme not completed	Terminate the programme completely
	Rotary / push button defective	Have it exchanged by an expert
Electronics can be set, but heating does not start	Cartridge holder not correctly in position A	Move cartridge holder up until the stop has been reached
	Microswitch defective	Have microswitch exchanged by an expert
	Start/stop button malfunction	Start / stop button control, if necessary exchange by an expert
Heating time expires, but device does not heat up	Quartz radiator is defective	Exchange quartz radiator
The display shows error code E01	Temperature monitor was activated	Reset temperature monitor and have it exchanged by an expert if necessary

Problem	Cause	Aid
Heating time expired, foil is not plasticised evenly and/or sufficiently	Heating time too short	Prolong heating time by 5–10 sec.
	Device was not preheated	Heat up the device after longer periods out of use for 0:30 min.
	Check heating element / temperature monitor for functioning	Reset temperature monitor and have it exchanged by an expert if necessary
	Service life of quartz radiator exceeded	Exchange quartz radiator Prolong heating time by 5–10 sec.
	Cartridge holder not lowered completely	Bring cartridge holder to position B
	Seal defective	Check seal and exchange it if necessary
	Filter of vacuum pump contaminated	Exchange filter
	Vacuum pump / supply lines defective	Have it checked by an expert and exchanged if necessary

17. Exploded view

18. Spare parts list

Nr.	Bezeichnung	REF
1	Plug combination element	51270
2	Main switch	50820
3	Rotary/push button	51322
5	Start/Stop button	570040
6	Control light heating elements	577947
7	Vacuum indicator	570150
8	Operating lever	D57011
9	Removable model reception	D57020
10	Cassette for blanks	D57012
11	Cassette holder	On request
12	Fuse 12A	570167
13	Heating cover	On request
14	Side plate right left	On request On request
15	Vacuum pump, complete	D57013
16	Operating electronics	5702448
17	Heating element	570140
18	Temperature limiter, 60 °C	570168
19	Temperature limiter 45 °C	570169
20	Holding magnet	540416
21	Fan	50840
22	Switch-mode power supply 24V DC	570211
23	Mains filter	570144
24	Microswitch	50844
25	Rubber feet set, 4 pcs.	D57015
26	Handle adaptor, right left	55045 55045links
27	Vacuum reservoir	D57021
28	Vacuum valve	5729910
29	Magnetic tongs	D57022
30	Cable resistance thermometer TR40	50712

Without drawing no.

Silicone seal thermo forming area	570143
Power cable, Europe	51021
Power cable, England	51284
Power cable, America	51283
Sheet metal parts for housing	On request
Air filter for vacuum pump	570170
Seal for blank cassette	D57018

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