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Intended use

Thermoforming blank for the production of dental splints

Indications for use

Suited for interim splints in case of myoarthropathy (muscular and mandibular discomfort) and accompanying dysfunctions resulting from it, like pressing, bruxism, etc. The splint provides for a very good bite lock and does not cause any tension. Furthermore the production is very easy and there is no need for anchorages, the splint fits very well and secure due to the snap-in effect.

Characteristics

flexible, soft-remaining resin
 Ø 120 mm / 125 mm
 □ 120 mm / 125 mm

Material thickness

1.0 / 1.5 / 2.0 / 3.0 / 4.0 / 5.0 mm

Colours

nature-transparent – 3 mm blanks in various colours available

Processing

All Dreve thermoforming units are suitable for the processing of Drufosoft® blanks (Drufoformat scan, Drufoformat scan, Drufoformat, Drufoformat SQ / TE, Vacfoformat U / P, Vacfoformat 2000). For plastification the blank will be placed into the plate reception and with the help of the radiant heater brought to the necessary temperature. The times vary between 1:30 and 2:45 minutes according to the unit type and the surrounding temperature. After end of the heating time the material hangs down → complete the thermoforming process. Before thermoforming plaster models should be watered and insulated. Depending

on the material thickness the cooling time under pressure / vacuum is 3–8:30 minutes. After the cooling time Drufosoft® can be cut by means of scissors or a scalpel. Drufosoft® is suitable for „additive thermoforming“. Through process of lamination various foils can be connected (**only with pre-forming unit under pressure**)

Cleaning

For cleaning, the use of MaxiFresh™ cleaning tablets is recommended.
 Please be aware of the MaxiFresh™ instruction for use.

Information on health

Drufosoft® is marked with a CE label and therefore fulfils the general requirement of the medical device guideline.

Keep away from children!

Only for dental use!

Storage conditions

At room temperature, i. e. 18–28°C.

Durability: see date of expiry.

Drufosoft® is not to be used when expiry date has passed.

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This product was developed for dental use and has to be used according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence before using it subject to its pur-

pose, especially in case the purpose is not mentioned in the instructions.

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