


Drufomat scan

Professional thermoforming

Very easy to handle high-end thermoforming unit with a maximum in user-friendliness and consistently precise final results. Built for all applications in dental thermoforming technology.

Easy to use keypad, integrated bar code scanner identifies blanks, working pressure of up to 6 bar. Additional safety factors for professional results Two-hand operation and LEDs that light up when heating is in operation.

www.dentamid.dreve.de



Technical data

Dimensions (H x W x D)	460 x 290 x 290 mm
Weight	16 kg
Power input	max. 315 W
Pressure variably adjustable	from 2–6 bar (19–87 psi)
Voltage	115–240 V / 50–60 Hz
Usable thermoforming foils	all from Ø 120 mm, round and square

D3300	Drufomat scan (red)
D3300X1	Drufomat scan (silver)
D3300X2	Drufomat scan (RAL-colour chart)

