



Dynax[®]

A whole system for impression taking.

Dreve Dynax[®]

A whole system for impression taking.

Our new Dynax[®] silicones guarantee both, precision and ease of use. The whole impression process takes less than 3 minutes. High **hydrophilic properties** provide excellent impression precision. Our innovative **snap-effect** enables long processing time but short time in the patient's mouth. We guarantee high safety through outmatched **mechanical properties**. Furthermore: All silicones are insertable for the digital process due to **scannable surface**.





Indication



Colour



Final hardness



Processing time



Time in mouth



Texture





Crown and bridge impressions

Inlay and onlay impressions

Implant impressions

Functional impressions

Fixation impressions

Relining impressions

Pressure marks detection

Temporaries production

mono.



Dynax[®] mono. The Allrounder.

Addition curing, stable, thixotropic. A precise impression for all frequent indications. With only one material!



lavender violet



58 shore A



90 sec.




90 sec.



medium type 2



168



Tear resistance in %



Crown and bridge impressions

Inlay and onlay impressions

putty.



Dynax[®] putty. The Problem Solver.

Non-sticky, easily kneadable. A reliable and proven aide suitable for complex impressions, e. g. reduced dentition with periodontal damage.



sunny yellow



60 shore A



90 sec.




90 sec.



kneadable-type 0



1:30



1:45



1:30



1:30



Processing time mins.




Crown and bridge impressions

Inlay and onlay impressions

Implants impressions

Fixations

heavy.



Dynax[®] heavy body. The Mighty.

Addition-curing „tray-type“-silicone. Very thixotropic, very stable cartridge material. Combine it with Dynax[®] light.



sky blue



70 shore A



90 sec.



90 sec.



high-type 1



3:00

Dynax® heavy body

3:30

Competitor A

3:30

Competitor B

4:30

Competitor C

5:30

Competitor D

Removable after... (incl. processing time)



Crown and bridge impressions


Inlay and onlay impressions

Relining impressions

Detection of pressure marks

light.





Dynax[®] light. The Corrector.

Addition-curing, thin flowing corrective silicone for outstanding detail precision. Hydrophilic surface and perfect thixotropy.



emerald green



43 shore A



90 sec.




90 sec.



low-type 3



2.5



Tear resistance




Crowns and bridges impressions

Inlays and onlays impressions

Implants impressions

clear.





Dynax[®] clear. The clear one.

UNIQUE! Crystal clear results with excellent flowing- and perfect hydrophilic properties. Enables checking details during impression taking.



transparent



60 shore A



90 sec.



120 sec.



medium-type 2



Dynax® clear





Dynax® mono

2 x 50 ml double cartridges **D6106**
incl. 6 mixing tips

8 x 50 ml double cartridges **D6107**

Dynax® heavy body

2 x 50 ml double cartridges **D6112**
incl. 6 mixing tips

8 x 50 ml double cartridges **D6113**



Dynax® light

2 x 50 ml double cartridges **D6103**
incl. 6 mixing tips

8 x 50 ml double cartridges **D6104**

Dynax® clear

2 x 50 ml double cartridges **D6121**
incl. 6 mixing tips

8 x 50 ml double cartridges **D6122**

Dynax® putty

2 x 450 ml tubs incl. dosing spoons D6101



Mixing tip green

Mixing ratio	1:1
Inner diameter	6.5 mm
Article numbers	D4970X2 (40 pieces) D49701X2 (100 pieces)
Products	Dynax® clear Dynax® heavy body Dynax® mono

Mixing tip yellow

Mixing ratio	1:1
Inner diameter	4.2 mm
Article numbers	D4970X5 (40 pieces) D49701X5 (100 pieces)
Product	Dynax® light





Dreve Dentamid GmbH

Max-Planck-Straße 31

59423 Unna/Germany

+49 2303 8807-40

www.dreve.de/dentamid