



# Naturally shaped

Aligner foils by Dreve.



*New  
3 thicknesses*

# Naturally shaped

## Aligner foils by Dreve.

### Aligner

The new technology for cosmetic soft-gentle orthodontics. Especially suitable for adult patients with high aesthetic standards, little time and slightly false positioned teeth.

### Our offer

High-transparency, tough elastic thermoforming foils for custom-made manufacturing of aligners.

### Layer thicknesses for each treatment phase

	soft	0.5 mm	Aligner + Retainer
	medium	0.625 mm	Aligner + Retainer
	hard	0.75 mm	Aligner + Retainer

### Your benefit

Biolon is an absolutely established complete foil program. Therefore each aligner treatment becomes safe and simple.

---

D420001	Biolon, soft 0.5 mm	25 Pieces
D420019	Biolon, medium 0.625 mm	20 Pieces
D420003	Biolon, hard 0.75 mm	20 Pieces