



Dreve

Dosper evo

Runs like a charm

Automatic mixing and dosing unit for thin-flowing, addition-curing duplication silicones. Thanks to its innovative flow-adjustment it is suitable for all duplication silicones with a 1:1 mixing ratio.

No unwanted emptying of the reservoirs or mixing up of components. The unit is constantly ready for use and very economical and profitable also for smaller labs.

www.dentamid.dreve.de

Dreve



Technical data

Dimensions (H x W x D) 320x225x370 mm

Weight 12 kg

Power input 250 W

Usable duplication silicones all with a 1:1 mixing ratio

D5050 Dospet evo (red, 230 V)

D5050A Dospet evo (red, 115 V)

D50501 Dospet evo (silver, 230 V)

D50501A Dospet evo (silver, 115 V)

