



Dreve

Polymax 5

Our power pack

Pressure polymerisation unit for hot and cold polymerisation of dental resins, with automatic safety seal system, sturdy and maintenance-free technology. Can take 1 large articulator, 3 brass flasks including clamps, or 4 Castdon flasks. Temperature adjustable to up to 95 °C, a 90 minutes timer and variably adjustable pressure ensure perfect results. Can also be used as pressure compactor.

www.dentamid.dreve.de



Technical data

Dimensions (H x W x D)	310 x 340 x 400 mm
Weight	15.6 kg
Power input	900 W
Pressure	3–6 bar (variable)
Ø Pressure pot	220 mm
Pressure pot depth	174 mm
Temperature range	0–95 °C
Frequency	50–60 Hz

D3423	Polymax 5 (red, 230 V)
D3423A	Polymax 5 (red, 115 V)
D34231	Polymax 5 (silver, 230 V)
D34231A	Polymax 5 (silver, 115 V)

