

Tips



Technical Data Sheet

Product description: For mixing tips for intra-oral application. Flattened form for bite registration materials.

Manufacturer: Drewe Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	chemical term
1	high density polyethylene

Properties:

Colour: transparent



Ratio: 1:1, 2:1, 4:1, 10:1

Ø outside: 7.6 x 2.8 (mm)

Ø inner: 6.4 x 1.6 (mm)

Drewe Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
dentamid@drewe.de
www.drewe.com