



Technical Data Sheet

Intended Use: Dental splinting material short-term

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical Denomination
1	Ethylen und Venylacetat

Physical characteristics:

Colour: transparent

Thickness: 1.00 mm

Diameter: Ø 120 mm ((special diameters on demand))

Density: 0.935 g/cm³ DIN 53479

Elongation at break: > 750 % DIN 53455

E-Modul: 75 MPa DIN 53457

Melting temperature: 92 °C D.S.C.

Flexural strength: 19 MPa DIN 53457

Hardness: 92 Shore A

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
dentamid@dreve.de
www.dreve.com