

# Drufolen H



## Technical Data Sheet

**Intended use:** Dental splinting material short-term

**Manufacturer:** Dreve Dentamid GmbH  
Max-Planck-Straße 31  
59423 Unna / Germany

**Composition:**

Pos	Chemical Denomination
1	polypropylene

**Characteristics:**

Colour: natural coloured  
Thickness: 0,6mm, 0,8mm  
Diameter: Ø 120mm and 42mm  
Density: 0.9 g/cm<sup>3</sup>  
Melting temperature: approx. 160°C  
Thermal decomposition: >300° C

When used according to specifications no harmful impacts are noted.

Dreve Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 8807-55

[info@dreve.de](mailto:info@dreve.de)

[www.dreve.com](http://www.dreve.com)