

Drufolen W



Technical Data Sheet

Product Description: Soft-elastic thermoforming foil for the production of moulded parts in the crown and bridge technique.

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical Denomination
1	polyethylene

Characteristics:

Colour: natural coloured / blue
Thickness: 0,45mm; 0,65mm, 0,8mm, 1,00mm
Diameter: Ø 120mm und 42mm
Density: 0.92 g/cm³
Melting temperature: approx. 100°C -140°C
Thermal decomposition: > 390° C

When used to specifications no harmful impacts are noted.

Dreve Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 8807-55

info@dreve.de

www.dreve.com