

Drufoplast



Technical Datasheet

Intended use: Dental splinting material short-term

Hersteller: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Zusammensetzung:

Pos	Chemische Bezeichnung
1	Polystyrol

Eigenschaften:

colour: rose/ natural
Thickness: 1,5 mm, 2,0 mm, 3,0 mm
Diameter: Ø 120 mm
Density: 1,18 g/cm³
Melting point: ca. 85°C
Burning point: >350° C

When used according to specifications no harmful impacts are noted.