



Technical Data Sheet

Intended use: Dental splinting material long-term

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Zusammensetzung:

Pos	Chemical term
1	Ethylene and vinylacetate-copolymer

Properties

Colour: various

Thickness: 1.0 mm, 1.5 mm, 2.0 mm, 3.0 mm,
4.0 mm, 5.0 mm

Diameter: Ø 120 mm (special diameter on desire)

Melting temp: approx. 74°C (165°F)

Hardness: 82 Shore A

Density: 0.96 g/cm³

Humidity absorption: < 0,2 %

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
dentamid@dreve.de
www.dreve.com