

# Drufosoft pro



## Technical Data Sheet

**Intended use:** Dental splinting material long-term

**Manufacturer:** Dreve Dentamid GmbH  
Max-Planck-Straße 31  
59423 Unna / Germany

### Composition:

Pos	Chemical term
1	Ethylene
2	Venyl acetate copolymer

### Properties:

Colour: natural clear

Thickness: 3.0 mm

Diameter: Ø 120 mm (special diameter on desire)

Melting point: 92°C according D.S.C.

Tensile strenght: 19 Mpa (according to DIN 53455)

E-module: 75 Mpa (according to DIN 53457)

Hardness: 92 Shore A (according to DIN 53505)

Dreve Dentamid GmbH  
Max-Planck-Str. 31  
59423 Unna / Germany  
Tel.: 02303 / 8807-0  
Fax: 02303 / 8807-55  
dentamid@dreve.de  
www.dreve.com