

DS 50 Injector 4:1



Technical Data Sheet

Purpose : DS-50 injector for the application of 4:1 silicones supplied in 50 ml double cartridges.

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Features:

Colour : grey
Cartridge volume : 50 ml
Volume ratio : 4:1 / 10:1
Conversion ratio : 1 / 9
Plunger feed per stroke: 6.5 mm
Volume per stroke : ca. 3.7 ml
Net weight : 276 g

Disinfection of DS dispensers :

Above mentioned DS dispenser was tested with regard to resistance against disinfectants. 1000 disinfection cycles were carried out according to manufacturers instructions. Following the disinfection endurance tests were made with the units. The disinfection process caused a slight reduction of rigidity as well as minor changes of colour to some materials. The disinfectants caused no mentionable impact on the units function and durability.

There is a vast choice of disinfectants in the market. This choice is constantly subject to modifications. Thus Dreve cannot take any general responsibility with regard to resistance against all commercially available disinfectants.

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55

info@dreve.de
www.dreve.com