

DuoSep



Technical Datasheet

Product description: Brushable separating agent for insulating silicones and polyether.

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	chemical designation
1	distilled water
2	polyvinyl alcohol
3	coating additives

Properties:

Colour: achromatic
Form: liquid
Odour: characteristic
Boiling point: 100° C
Viscosity: 600 – 1000 mPas

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
info@dreve.de
www.dreve.com