

Finishing Liquid



Technical Data Sheet

Product description: Liquid for surface conditioning of Drufosoft® thermoforming-blanks.

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	chemical classification
1	eucalyptus oil
2	acetone

Properties:

Colour: achromatic

Odour: eucalyptus (intensive)

Form: liquid

Boiling point: 56°C

Burning point: -17°C

Inflammation point: 540°C

Density: 0.89 g/cm³

Vapour pressure: 233 hPa (at 20° C)

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-40
Fax: 02303 / 8807-55
dentamid@dreve.de
www.dreve.com