

GipsEx



Technical Data Sheet

Intended use: Liquid for the removal of plaster and gypsum bounded investment on polymer-dentures.

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical term
1	Water
2	Tri-sodiumcitrate-5,5.hydrate

Properties:

Colour: clear

Form: liquid

Boiling Point: 100° C

Density: 1.18 g/cm³

Vapour pressure: 23 hPa (at 20° C)

Dreve Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 82909

dentamid@dreve.de

www.dreve.com