

GumQuick scan



Technical Data Sheet

Intended use: Addition-vulcanising vinylpolysiloxane for the production of flexible gingiva masks on the master model.

Manufacturer: Drewe Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical term
1	Vinyl containing polydimethylsiloxane
2	Polymethylhydrogensiloxane
3	Silicic acid
4	Filling materials
5	Pigments
6	Platinum-catalyst

Properties:

Colour: pink (gingiva colour)

Hardness Shore A: 63 ± 3 (DIN 53505)

Consistency: Type 3: light (ISO 4823)

Elastic recovery: $> 99.6 \%$ (ISO 4823)

Linear dimensional change: $< -0.1 \%$ (ISO 4823)

Processing time: 1:30 min at 23°C

Deflasking time: $\geq 5:00$ min at 23°C

Drewe Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 8807-55

dentamid@drewe.de

www.dentamid.drewe.com