

GumQuick scan



Technical Data Sheet

Intended use: Addition-vulcanising vinylpolysiloxane for the production of flexible gingiva masks on the master model.

Manufacturer: Drewe Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical term
1	Vinyl containing polydimethylsiloxane
2	Polymethylhydrogensiloxane
3	Silicic acid
4	Filling materials
5	Pigments
6	Platinum-catalyst

Properties:

Colour: pink (gingiva colour)

Hardness: 63 ± 3 Shore A (DIN 53505)

Consistency: light-bodied – type 3 (ISO 4823)

Elastic recovery: $> 99.6 \%$ (ISO 4823)

Linear dimensional change: $< 0.5 \%$ (ISO 4823)

Processing time: 1:30 min at 23°C

Demoldability: $\geq 6:30$ min at 23°C

These data result from measurements which are executed from time to time
In view of our quality management system.
This document is valid without any signature.

Drewe Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 8807-55

dentamid@drewe.de

www.dentamid.drewe.com

Rev. 2/ 05.2025

Page 1 of 1

TDS_GumQuick scan_GB