

Laborsil clear



Technical Data Sheet

Product description: Addition-vulcanising silicone for the universal use in the dental lab

Manufacturer: Drewe Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical term
1	Polyvinylmethylsiloxane
2	Vinyl-Q Resin
3	Polymethyl hydrogensiloxane
4	Filling materials
5	Platinum catalyst
6	Pigments

Properties:

Colour: transparent

Hardness: 70 ± 2 [Shore D] (DIN ISO 7619-1)

Consistency: Type 2: medium (DIN EN ISO 4823)

Recovery from deformation: > 99.2 % (DIN EN ISO 4823)

Linear dimensional change: < 0.1 % (DIN EN ISO 4823)

Processing time: 30 sec ± 5 sec at 23°C

Demoulding: ≤ 90 sec at 23°C

Drewe Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
dentamid@drewe.de
www.drewe.com