

Laborsil speed

with colorindicator



Technical Data Sheet

Product description: Addition-vulcanising silicone for the universal use in the dental lab

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Dreve Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 8807-55

dentamid@dreve.de

www.dreve.com

Composition:

Pos	Chemical term
1	Polyvinylmethylsiloxane
2	Vinyl-Q Resin
3	Polymethyl hydrogensiloxane
4	Filling materials
5	Platinum catalyst
6	Pigments

Properties:

Colour: transparent

Hardness: 75 ± 2 [Shore D] (DIN ISO 7619-1)

Consistency: Type 2: medium (DIN EN ISO 4823)

Recovery from deformation: > 99.6 % (DIN EN ISO 4823)

Linear dimensional change: < 0.2 % (DIN EN ISO 4823)

Processing time: 30 sec ± 5 sec at 23°C

Demoulding: ≤ 2:30 min at 23°C

These data result from measurements which are executed from time to time in view of our quality management system.
This document is valid without any signature.

Rev. 0 / 02.2021
Page 1 of 1
TDS_Laborsil speed with
colorindicator_GB