

Lightplast base plates



Technical Data Sheet

Product description: Light polymerizing resin base-plates (upper and lower jaw) for the production of custom trays.

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical term
1	polyfunctional methacrylates
2	glass beads, silane coated
3	silicium dioxide, high dispersed
4	camphorquinone
5	pigments

Properties:

Colour : nature, pink
Strength of plates: $2,2 \pm 0,2$ mm
Penetration: 4,0 – 7,0 mm
Bending strength: ≥ 85 N/mm²

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
dentamid@dreve.de
www.dreve.com