

MegaBase[®] 2K Lack



Technical Data Sheet

Intended use: Soft denture relining material

Extended intended use: Silicone lacquer for coating of dental soft relining material.

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
dentamid@dreve.de
www.dentamid.dreve.com

Composition:

| Pos | Chemical term |
|-----|----------------------------|
| 1 | vinyl polydimethylsiloxane |
| 2 | VQM-resins |
| 3 | platinum catalyst |
| 4 | methylhydrogensiloxane |
| 5 | silicon dioxide |

Properties:

Colour: clear – transparent

Consistency: low – type 3 (DIN EN ISO 4823)

Final hardness: 25 ± 2 Shore A (DIN ISO 7619-1)

Processing time: ≥ 2:15 min. at 23 °C

Total setting time: 20:00 min. at 23 °C

At room temperature: 10:00 min. at 23 °C

Afterwards in the water bath: 10:00 min. at 35 – 40 °C

Linear dimensional change: < 0.2% (DIN EN ISO 4823)

Storage: at room temperature (18 °C – 28 °C; 64 °F – 82 °F)

These data result from measurements which are executed from time to time
In view of our quality management system.
This document is valid without any signature.

Rev. 3 / 08.2022
Page 1 of 1
TDS_MegaBase[®] 2K
Lack_GB