

# Mixing Canulas



## Technical Data Sheet

### Product description:

For the mixing of two-component materials from standard customary double cartridges or double syringes.

### Composition:

Pos	Denomination
1	Canula Housing
2	Mixing Element PP resp. POM

Dreve Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 8807-55

[info@dreve.de](mailto:info@dreve.de)

[www.dreve.com](http://www.dreve.com)

Mixing Canula	Cartridge System	Colour Identification	Inner diameter mm	Length mm	Contents/ Volume ml	Mixing Elements
MA 5,4 – 13D	48 ml	--	5,4	79,0	1,20	13
MA 6,3 – 12D	48 ml	--	6,3	58,8	1,82	12
MBT 3,2 – 12 S	50ml	Blau	3,2	49,8	0,36	12
MBT 4,2 – 12 D	50ml	Gelb	4,2	50,1	0,59	12
MBT 5,4 – 12 D	50ml	Pink	5,4	59,4	1,02	12
MBT 6,5 – 12 D	50ml	Grün	6,5	67,6	1,60	12
MBTX 3,2- 16 S	50ml	Blau 4:1	3,2	57,6	0,46	16
MLT 2,5 – 10 S	5ml	Braun	2,5	34,4	0,13	10
MBT 381- 05 HS	380 ml	Gelb 5:1	6,0	52,5	4,80	5

These data result from measurements which are executed from time to time in view of our quality management system.  
This document is valid without any signature.

Rev. 19.04.2017  
Seite 1 von 1  
41e-3928/CW