

NanoVarnish



Technical Data Sheet

Intended use: Light-curing lacquer for plastic dental applications
Extended purpose: Light-curing lacquer for splints, partial and full dentures and temporary dental crowns in C&B technology

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical term
1	methyl methacrylate
2	multifunctional urethane acrylate
3	SiO ₂ -nanoparticles
4	initiator
5	additives

Properties:

Colour: clear transparent
Form: liquid
Viscosity: 0.005 – 0.015 Pa*s
Density: 1.07 g/cm³

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
dentamid@dreve.de
www.dreve.com