

Plaquit



Technical Data Sheet

Intended use: Light-curing lacquer for plastic dental applications
Extended purpose: Light-curing lacquer for splints, partial and full dentures and temporary dental crowns in C&B technology

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical term
1	Methyl methacrylate
2	Aliphatic urethane acrylate
3	Initiators
4	Additive

Properties:

Colour: clear transparent
Form: liquid
Viscosity: 0.01- 0.02 Pa*s
Density: 0.98 g/cm³

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: + 49 2303 / 8807-40
Fax: +49 2303 / 8807-55
dentamid@dreve.de
www.dreve.com