

Regofix

opaque



Technical Data Sheet

Intended use: Bite registration material based on vinyl polysiloxane

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical term
1	Polyvinylmethylsiloxane
2	Polydimethylsiloxane Oil
3	Filling materials
4	Platinum catalyst
5	Polymethylhydrogensiloxane
6	Pigments

Properties:

Colour: blue

Hardness Shore A: 90 ± 5

Consistency: heavy - body - Type 1 (ISO 4823)

Elastic recovery: > 99.8 % (ISO 4823)

Linear dimensional change: < 0.2 % (ISO 4823)

Processing time: 30 ± 5 sec at 23 °C

Deflasking time: ≥ 90 sec at 37 °C

These data result from measurements which are executed from time to time
In view of our quality management system.
This document is valid without any signature.

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
dentamid@dreve.de
www.dreve.com

Rev. 9 / 08.2023
Page 1 of 1
TDS_Regofix opak_GB