

Regofix

transparent



Technical Data Sheet

Intended use: Bite registration material based on vinyl polysiloxane

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical term
1	Polyvinylmethylsiloxane
2	Polymethylsiloxane-Oil
3	Silicic Acid
4	Platinum catalyst
5	Polymethylhydrogensiloxane
6	Pigments

Properties:

Colour: transparent

Hardness Shore A: 70 ± 2

Consistency: medium - Type 2 (ISO 4823)

Elastic recovery: $> 99.6 \%$ (ISO 4823)

Linear dimensional change: $< 0.2 \%$ (ISO 4823)

Processing time: 30 ± 5 sec at $23 \text{ }^\circ\text{C}$

Deflasking time: ≥ 90 sec at $37 \text{ }^\circ\text{C}$

Dreve Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 8807-55

dentamid@dreve.de

www.dreve.com