

Siliform



Technical Datasheet

Product description: Siliform is an embedding material on the basis of A-silicone especially for the resin pouring technique.

Manufacturer: Drewe Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	chemical designation
1	vinyl containing polydimethylsiloxane
2	silicic acid
3	filling materials
4	pigments

Properties:

Colour : yellowish
Final hardness: 14 Shore A
Mixing time: 45 sec
Processing time : 5-7 min (incl. Mixing time) (at 23°C)
Deflasking time: 30 min (at 23°C)
Dosage: Siliform is mixed in a ratio of 1:1 of the components A+B

General notes: Siliform is not suitable for the duplication of plaster models or for the production of embedding masses models.

Drewe Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
info@drewe.de
www.drewe.com