



## Technical Data Sheet

Product description: Bite registration material based on vinyl polysiloxane

Manufacturer: Drewe Dentamid GmbH  
Max-Planck-Straße 31  
59423 Unna / Germany

Drewe Dentamid GmbH  
Max-Planck-Str. 31  
59423 Unna / Germany  
Tel.: 02303 / 8807-0  
Fax: 02303 / 8807-55  
dentamid@drewe.de  
www.drewe.com

### Composition:

Pos	Chemical term
1	Polyvinylmethylsiloxane
2	Vinyl Q-Resin
3	Silicic Acid
4	Polymethylhydrogensiloxane
5	Filling materials
6	Pigments
7	Platinum catalyst

### Properties:

Colour: grey

Hardness Shore D: 32 ± 2

Consistency: Type 2: medium (ISO 4823)

Elastic recovery: > 99.9 % (ISO 4823)

Linear dimensional change: < 0.2 % (ISO 4823)

Processing time: 30 sec ± 5 sec at 23 °C

Deflasking time: 80 sec ± 15 sec at 37 °C