

# Ultra Trimm



## Technical Datasheet


Product description: Flexible grinding wheels to treat soft materials

Manufacturer: Dreve Dentamid GmbH  
Max-Planck-Straße 31  
59423 Unna / Germany

Composition:

Pos	Chemical term
1	aluminium oxide
2	calcium carbonate
5	Cured adhesive
4	Webbed nylon

Processing: A grinding wheel will be screwed on a mandrel inclusive a support disc .  
Max. rotations: 15000 Upm



Properties: three different final hardness for rough up to fine finishing.  $\varnothing$ :28 mm

Light: Granulation: S super fine ; thickness: 7 mm

Medium: Granulation: A fine ; thickness: 7 mm

Heavy: Granulation: CP, A medium ; thickness:14 mm

Dreve Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 8807-55

[info@dreve.de](mailto:info@dreve.de)

[www.dreve.com](http://www.dreve.com)