

Unisol E



Technical Datasheet

Product description: Decalcifier for maintenance of strongly used wax boilingout units.

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	chemical designation
1	citric acid monohydrate

Properties:

Colour: white
Form: crystalline
Odour: odourless
pH-value: 1,6
Melting temperature: 153° C

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
info@dreve.de
www.dreve.com